

Product datasheet for SC328766

KLHDC4 (NM_001184854) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	KLHDC4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PSI00001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC328766 representing NM_001184854. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTGAACCGTCAGAAATTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGCAAGAAGGGCAAGAAGGAGAAGAAGGGCCGCGCGGAGAAGACGGCCCAAGATGGAGAAG
AAGGTGTCTAAGCGCTCGCGGAAGGAGGAGACTTTTTGTATAACGAGCTCTATGTCTACAATACCAGA
AAGGACACCTGGACCAAGTTGACATCCCCAGTCCACCTCCGAGGCGCTGTGCTCACCAGGCGGTGTA
GTGCCTCAAGGTGGCGGACAGCTGTGGTCTTTGGAGGGGAGTTTGCTCTCCCAACGAGAGCAGTTTC
TACCACTACAAGGATCTCTGGTCTCTGCATTTGGCCACCAAGACCTGGGAACAAGTCAATCAACAGGC
GGTCTTCGGGTGCGAGTGGACATCGGATGGTGGCCTGGAAGAGACAATTGATCCTGTTGGTGGCTTC
CATGAAAGTACACGGGATTACATCTACTACAACGACGTGTATGCCTTTAATCTGGACACCTTCACATGG
AGCAAGCTGTCCCGTACGGGACGGGGCCACACCCAGATCAGGCTGCCAGATGTCCGTCACCTCCCGAG
GGCGGCATCGTCTGTATGGGGCTACTCGAAACAGAGAGTTAAGAAAGACGTGGACAAGGGCACACGG
CACTCAGACATGTTCTGCTGAAGCCAGAGGACGGAAGAGAAGACAAGTGGGTTTGGACTCGGATGAAC
CCTTCGGGGGTCAAGCCACCCACGGTCTGGCTTTCCGTGGCCATGGCCCCGAATCACCAGACACTG
TTCTTCGGGGGTGTCTGTGACGAGGAAGAGGAGAGCCTGTGGGCGAGTTCTTCAACGATCTGTAC
TTCTACGACGCCACCAGGAACCGTTGGTTTGAGGGACAGCTGAAGGGACCAAGTCTGAAAAGAAGAAA
CGCAGGCGGGGCAGAAAAGAGGAGCCCGAAGGTGGTAGCAGGCCGGCGTGTGGGGGAGCTGGCACCAG
GGGCTGTGCAGCTGGTCAAGGAGGTGGTGGCCGAGGATGGACCCGTGGTACCATTAAAGCAGGTGCTC
ACCGCGCCAGGCTCGGCGGGGACGCCCGGTCTGAGGACGAAGACAGCCTTGAGGAGGCCCGGACGCCCC
GCACCTGGGCGGTGCCAGCTCCAACGCCATGCTGGCTGTGAAGCATGGGGTGTCTACGTCTATGGG
GGCATGTTTGAAGCCGGCGACCGCCAGGTACCCTCAGCGACCTGCACTGCCTGGACCTGCACAGGATG
GAGGCGTGGAAGGCCTTGGTGGAGATGGACCCAGAACTCAGGAGTGGCTGGAGGAGACGGACTCGGAA
GAGGACAGTGAGGAGGTTGAGGGCGCCGAGGGTGGGTCGACGACGAAGACAGCGGAGAGAGCGGT
GCGGAGGACTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001184854
Insert Size:	1392 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001184854.1</u>
RefSeq Size:	1770 bp
RefSeq ORF:	1392 bp
Locus ID:	54758
UniProt ID:	<u>Q8TBB5</u>
Cytogenetics:	16q24.2
MW:	51.5 kDa