

Product datasheet for **SC328781**

NOSTRIN (NM_001171632) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NOSTRIN (NM_001171632) Human Untagged Clone
Tag:	Tag Free
Symbol:	NOSTRIN
Synonyms:	DaIP2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC328781 representing NM_001171632.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGGGACCCACTGACAGATTGTCGTATAATAAAGTATACAAGAACCTAAAGGAGTTTTCTCAAAT
GGAGAGAATTTTGCAAACAGGTCACATCTGTTCTTCAGCAAAGTTGTGTAGCAGTGCCTGGCCCTGG
GCCTCAGAGGGAATGAAATCCACAGCGGACCTGCATCAAAAACCTGGCAAAGCAATTGAATTGGAAGCA
ATAAAACCGACTTATCAAGTCCTAAATGTACAAGAGAAGAAGAGAAAATCACTTGACAATGAAGTTGAA
AAGACAGCAAATCTTGTCATTAGCAACTGGAATCAGCAAATTAAGGCCAAGAAGAAATTAATGGTTAGT
ACCAAGAAACATGAAGCACTTTTCCAGCTTGTAGAAAGCTCCAAGCAATCTATGACTGAGAAGGAGAAG
CGGAAGCTCCTCAATAAACTGACAAAATCAACTGAAAAGTTGGAAAAGGAAGATGAAAATTACTACCA
AAAAACATGGCGGGTTATTCTACCAGACTGAAATGGGAAAACACACTAGAGAAGTCTACCAGAGCATT
CTGGAGCTGGAGAAGGAAAGAATTCAACTTTTATGCAATAACTTAAACCAGTACAGCCAACATTTTCT
CTTTTTGGCCAAACCCTGACCACATGCCACACGCAGATTCAGTGTCCATCAGCAAGATTGACATTGAA
AAAGATATCCAGGCTGTAAATGGAAGAACTGCAATTTTATCTACAGAAAACAAATCTGAGTTCCTGTTA
ACGGATTACTTTGAAGAAGATCCTAACAGTGCAATGGATAAAGAGAGACGAAAAGTCTTTACTAAAACCA
AAATTATTGAGACTGCAGAGAGACATTGAAAAGCCTCAAAGACAAGGAAGGCCTGGAACGAATGCTT
AAAACGTAATCCAGCACCTCCTCCTTCTCTGATGCAAAGACCCAGAAAGACACAGCAGCGTTAATGGAT
GAGAACAATTTGAAACTAGACCTTTTGGAAAGCGAACTCTACAACTGTCATCAATGTTAGCAGAACTT
GAGCAAAGACCTCAACCCAGCCATCCTTGATAGTAATCCATCTTCAGGTGGAGGGAAAAGGAGCATACT
CATAGCTATGTGAAAATATCTCGGCCTTTTTAATGAAGAGATTAGAGAATATTGTGAGCAAGGCATCT
TCTGGTGGCAGAGCAATCCAGTTCTTCAACTCCAGCCCTGGTGCAGCCAGCTCAGCAGCAGACTT
TGCAAGGCCTTGTATTCTTTTCAAGCCAGGCAAGATGATGAGTTGAATTTGGAAAAGGGTGACATTGTG
ATTATACAGGAAAAAAGAAGGAGGATGGTGGTTGGATCTTTGAATGGAAAAAAGGCCATTTTCT
GCCGCTTATGTGGAGGATTACCTTCAAATGCTGGCAACACAGCTACAAAGGCATTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001171632

Insert Size: 1437 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001171632.1</u>
RefSeq Size:	1863 bp
RefSeq ORF:	1437 bp
Locus ID:	115677
UniProt ID:	<u>Q8IV19</u>
Cytogenetics:	2q24.3
MW:	54.5 kDa
Gene Summary:	<p>Nitric oxide (NO) is a potent mediator in biologic processes such as neurotransmission, inflammatory response, and vascular homeostasis. NOSTRIN binds the enzyme responsible for NO production, endothelial NO synthase (ENOS; MIM 163729), and triggers the translocation of ENOS from the plasma membrane to vesicle-like subcellular structures, thereby attenuating ENOS-dependent NO production.[supplied by OMIM, Apr 2004]</p> <p>Transcript Variant: This variant (3) lacks an exon in the 5' CDS, as compared to variant 2. The resulting isoform (3) lacks an internal segment, as compared to isoform 2.</p>