

Product datasheet for **SC328789**

ZNF177 (NM_001172651) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF177 (NM_001172651) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF177
Synonyms:	PIGX
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001172651, the custom clone sequence may differ by one or more nucleotides

```

ATGGCTGCAGGGTGGCTGACAACCTGGTCACAGAACTCAGTAACCTTCCAGGAAGTGCCA
GTGGACTTTTCCAGGAGGAGTGGGCATTGCTGGACCCTGCTCAAAAAATCTATACAAA
GATGTGATGCTGGAGAACCTTAGGAACCTGGCCTCAGTAGGGTATCAGCTCTGCAGACAC
AGTCTGATCTCCAAGGTGGATCAAGAACAGCTGAAGACAGATGAAAGAGGAATTTTACAA
GGTACTGTGCAGACTGGGAACTCAACTTAAACCAAAAGATACAATTGCTATGCAGAAC
ATTCTGGGGGAAAAACATCCAATGGCATAAACACGGCAGAAAATCAACCTGGTGAGCAC
TCCCTGGAGTGTAAACATTGTGGGAAATCAGAAAGAACAACCTCGCTTTATTTGTACAAGA
TATTGCAAGGGAGAGAAATGCTATAAATATATAAAGTATAGCAAAGTCTTCAACCATCCC
TCAACTCTTAGGAGTCATGTGAGCATTACATTGGAGAGAAAACCTTGAATTTACTGAT
TGTAGAAAAGCTTTCAATCAAGAGTCATCCCTCAGGAAACACTTAAAGAACTCCCACAGGA
CAGAAGTTTCAGGAGTATGAGCAATGTGATATGTCCTTCAGCCTACACTTTCCTGCTCA
GTACGTGAGCAAAATACCTACTGGAGAGAAAAGGTGATGAATGCAGTGACTATGGCAAAATA
TCTCCCTTAGTGTCCACACAAAACCTGGTAGTGTGGAGGAGGGTTTGAATGTAATGAA
CATGAGAAAACCTTCACTGACCTTTGTCCTTCAGAAGTGTGTCAGAAGTCACTCTGGA
GAGATGCCCTATGAATGCAGTGACTGTGGGAAAGCCTTCATTTTTTCAGTCTTCCCTTAAG
AAACACATGAGATCTCATACTGGAGAGAAGCCTTATGAGTGTGATCACTGTGGAAAATCC
TTTAGCCAGACTCTCATCTGAATGTGCACAAAAGAAGTCAACTGGAGAGAAAACCTAT
GACTGTAAGGAATGTGGGAAGGCTTTCACTGTTCCCTTCATCCCTCAGAAAACATGTGAGA
ACCCACACTGGAGAGAAAACCTATGAATGCAGTGACTGTGGAAAAGCCTTCATCGATCAG
TCATCCCTTAAGAAAACACACAGCTCTCACACTGGAGAGAAGCCTTATGAGTGTAAACAG
TGTGGAAAGTCCCTCAGCACAGGCTTTACCTTATTGTGCACAAGAGAAGTCACTGGT
GAGAAAACCTATGAGTGTAAAGAATGTGGGAAGGCTTTAGGAATTCCTCTTGCCTGAGG
GTACACGTGAGAACTCACACTGGAGAGAAGCCTTATAAATGTATTCAAGTGTGAAAAGCC
TTTAGCACAAGCACTAACCTTATAATGCACAAGCGAATCCACAATGGCCAGAAAACCTCAT
GAATGA

```



[View online >](#)

Restriction Sites:	Please inquire
ACCN:	NM_001172651
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001172651.1</u> , <u>NP_001166122.1</u>
RefSeq Size:	2459 bp
RefSeq ORF:	1446 bp
Locus ID:	7730
UniProt ID:	<u>Q13360</u>
Cytogenetics:	19p13.2
Protein Families:	Transcription Factors
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (a). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.