

Product datasheet for **SC328814**

SSH3BP1 (ABI1) (NM_001178116) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SSH3BP1 (ABI1) (NM_001178116) Human Untagged Clone
Tag:	Tag Free
Symbol:	ABI1
Synonyms:	ABI-1; ABLBP4; E3B1; NAP1BP; SSH3BP; SSH3BP1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC328814 representing NM_001178116.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCAGAGCTGCAGATGTTACTAGAGGAGGAGATCCCCTGTGGCAAGAGGGCGCTGATCGAGAGTTAC
CAGAACCTGACTCGGGTGCAGACTACTGTGAAAACAACACTACATACAGAGACATGGTTTTGCTGTGTTG
CTGTGTCTGCTCAAACCTCCTGGCCTCAGGCTACAGACAAGAGAAAAGCTTTAGAGGAGACCAAAGCC
TATAACAACCAATCTCTAGCTAGTGTGCTTATCAAATAAATGCATTGGCCAACAATGACTCCAGTTG
CTGGATATCCAAGCCTCTCAGCTTCGGAGAATGGAGTCTTCCATCAATCATATCTCACAGACTGTGGAT
ATTCATAAGGAGAAAAGTGGCACGAAGAGAGATTGGTATTTTGACAACAAAATAAGAATACATCAAGAACT
CACAAAATAATAGCACCTGCGAATATGGAGCGCCTGTAAGGTATATTCGGAAACCTATCGATTACACA
GTTCTGGATGATGTGGCCATGGTGTCAAGTGGCTAAAAGCCAAGCATGGAAATAACCAGCCTGCAAGA
ACTGGCACACTGTCGAGAACAAATCCTCCTACTCAGAAACCGCCAAGTCTCCATGTCAGGCCGGGGA
ACACTGGGACGGAATACTCCTTATAAAACCTGGAACCTGTTAAACCCCAACAGTTCCTAATGACTAT
ATGACCAGTCTGCTAGGCTTGGAGTCAAGTCCAGGCAGGACAGCATCTTAAATCAGAGACCA
AGGACACACAGTGGAAAGTAGTGGAGGAAGTGGAAAGTCGAGAAAACAGTGGTAGCAGTAGTATTGGCATT
CCCATTGCTGTGCCTACACCTTCGCCACCCACTATTGGACCAGAAAACATTTCTGTCCCTCCTCCTTCT
GGAGCTCCACCAGCACCCTCTGGCACCCTTCTCCAGTGAGCACTGTGATAGCAGCCCCGGGCTCA
GCTCCTGGTCCCAGTATGGACAATGACCAGGCAGATATCTCGACACAACCTCTACTACTTCTTCGACA
TCTTCTGGTGGATACAGACGAACTCCCTCTGTGACTGCTCAATTTCTGCTCAGCCTCATGTTAATGGA
GGTCCACTTTATCTCAAATTC AATTGCTGATAGTCCAACCCACCACCACCTCCACCAGATGAC
ATTCCCATGTTTGTGACTCTCCACCTCCCCACCACCACCAGTGGATTATGAAGATGAGGAGGCT
GCAGTAGTTCAGTATAATGATCCATATGCAGATGGGGATCCTGCTTGGGCCCAAGAATTATATTGAG
AAAGTTGTTGCAATATATGATTATACAAAAGACAAGGATGATGAGCTGTCATTTATGGAGGGTGCAATC
ATTTATGTTATAAAGAAGAATGATGATGGCTGGTATGAAGGAGTCTGCAATCGAGTGACTGGTCTGTTT
CCTGGAACTATGTTGAATCAATCATGCACTATACTGATTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001178116

Insert Size:	1491 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001178116.1</u>
RefSeq Size:	3689 bp
RefSeq ORF:	1491 bp
Locus ID:	10006
UniProt ID:	<u>Q8IZP0</u>
Cytogenetics:	10p12.1
MW:	53.9 kDa

Gene Summary:

This gene encodes a member of the Abelson-interactor family of adaptor proteins. These proteins facilitate signal transduction as components of several multiprotein complexes, and regulate actin polymerization and cytoskeletal remodeling through interactions with Abelson tyrosine kinases. The encoded protein plays a role in macropinocytosis as a component of the WAVE2 complex, and also forms a complex with EPS8 and SOS1 that mediates signal transduction from Ras to Rac. This gene may play a role in the progression of several malignancies including melanoma, colon cancer and breast cancer, and a t(10;11) chromosomal translocation involving this gene and the MLL gene has been associated with acute myeloid leukemia. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome 14. [provided by RefSeq, Sep 2011]

Transcript Variant: This variant (5) contains one alternate exon and lacks a different in-frame coding exon, compared to variant 1. The resulting protein (isoform e) has the same N- and C-termini but is shorter when it is compared to isoform a. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.