

Product datasheet for **SC328816**

MEF2A (NM_001171894) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MEF2A (NM_001171894) Human Untagged Clone
Tag:	Tag Free
Symbol:	MEF2A
Synonyms:	ADCAD1; mef2; RSRFC4; RSRFC9
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC328816 representing NM_001171894.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGGCGGAAGAAAATACAAATCACACGCATAATGGATGAAAGGAACCGACAGGTCACTTTTACAAG
AGAAAGTTTGGATTAATGAAGAAAAGCCTATGAACCTAGTGTGCTCTGTGACTGTGAAATAGCACTCATC
ATTTTCAACAGCTCTAACAAACTGTTTCAATATGCTAGCACTGATATGGACAAAGTTCTTCTCAAGTAT
ACAGAATAAATGAACCTCATGAAAGCAGAACCAACTCGGATATTGTTGAGACTTTAAGAAAGAAAAGGC
CTTAATGGTTGTGAGAGCCCTGATGCTGACGATTACTTTGAGCACAGTCCACTCTCGGAGGACAGATTC
AGCAAATAAATGAAGATAGTGATTTTATTTTCAAACGAGGCCCTCCTGGTCTGCCACCTCAGAATTT
TCAATGTCTGTACAGTTCAGTGACCAGCCCCAATGCTTTGTCTACACTAACCAGGGAGTTCCTG
GTGTCCCATCTTTGGCAGCCAGCTCACGTTAACAGATTCAGCATGCTCTCTCCACCTCAAACCACA
TTACATAGAAATGTGTCTCCTGGAGCTCCTCAGAGACCACCAAGTACTGGCAATGCAGGTGGGATGTTG
AGCACTACAGACCTCACAGTGCCAAATGGAGCTGGAAGCAGTCCAGTGGGAATGGATTTGTAAGTCA
AGAGCTTCTCAAATTTGATTGGAGCTACTGGTGCAAATAGCTTAGGCAAAGTCATGCCTACAAAGTCT
CCCCCTCCACCAGGTGGTGGTAATCTTGAATGAACAGTAGGAAACCAGATCTTCGAGTTGTATCCCC
CCTTCAAGCAAGGGCATGATGCCTCCACTAAATACCCAGAGGATCAGTAGTTCTCAAGCCACTCAACCT
CTTGCTACCCAGTCGTGTCTGTGACAACCCCAAGCTTGCTCCGCAAGGACTTGTGTACTCAGCAATG
CCGACTGCCTACAACACTGATTATTCAGTACCAGCGTGACCTGTGACGCCCTTCAAGGCTTCAACTCG
CCAGGAATGCTGTGCTGGGACAGGTGTCGGCTGGCAGCAGCACCACCTAGGACAAGCAGCCCTCAGC
TCTCTTGTGCTGGAGGGCAGTTATCTCAGGGTTCCAATTTATCCATTAATACCAACCAAAACATCAGC
ATCAAGTCCGAACCGATTTACCTCCTCGGGATCGTATGACCCCATCGGGCTTCCAGCAGCAGCAGCAG
CAGCAGCAGCAGCAGCAGCCGCCACCACCGCAGCCAGCCAGCCACAACCCCGCAGCCCGCAGCCCGA
CAGGAAATGGGGCGCTCCCTGTGGACAGTCTGAGCAGCTCTAGTAGCTCTATGATGGCAGTGATCGG
GAGGATCCACGGGGCGACTTCCATTCTCCAATTGTGCTTGGCCGACCCCAAACTGAGGACAGAGAA
AGCCCTTCTGTAAGCGAATGAGGATGGACGCGTGGTGACCTAA
AGGCGACCGACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
  
```

Restriction Sites: SgfI-RsrII

Plasmid Map:



ACCN: NM_001171894

Insert Size:	1494 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001171894.2
RefSeq Size:	5838 bp
RefSeq ORF:	1494 bp
Locus ID:	4205
UniProt ID:	Q02078
Cytogenetics:	15q26.3
Protein Families:	Transcription Factors
MW:	53.5 kDa
Gene Summary:	<p>The protein encoded by this gene is a DNA-binding transcription factor that activates many muscle-specific, growth factor-induced, and stress-induced genes. The encoded protein can act as a homodimer or as a heterodimer and is involved in several cellular processes, including muscle development, neuronal differentiation, cell growth control, and apoptosis. Defects in this gene could be a cause of autosomal dominant coronary artery disease 1 with myocardial infarction (ADCAD1). Several transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Jan 2010]</p> <p>Transcript Variant: This variant (5) lacks an in-frame coding exon and contains an additional 5' non-coding exon, compared to transcript variant 6. These differences result in a shorter isoform (2), compared to isoform 5. Variants 2 and 5 both encode isoform 2. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>