

Product datasheet for **SC328835**

BANP (NM_001173539) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	BANP (NM_001173539) Human Untagged Clone
Tag:	Tag Free
Symbol:	BANP
Synonyms:	BEND1; SMAR1; SMARBP1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC328835 representing NM_001173539.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGATGTCCGGAACACGACCTGGCCGATGTGGTTCAGATTGCAGTGAAGACCTGAGCCCTGACCACCCA
GTTGTTTTGGAGAATCATGTAGTGACAGATGAAGACGAACCTGCTTTGAAACGCCAGCGACTAGAAATC
AATTGCCAGGATCCATCTATAAAGTCATTCTGTATTCCATCAACCAGACAATCTGCTTGCGGTTGGAT
AGCATTGAAGCCAAATTGCAAGCCCTGGAGGCTACTTGTAAATCCTTAGAAGAAAAGCTGGATCTGGTC
ACGAACAAGCAGCACAGCCCATCCAGGTCCCCATGGTGGCCGGCTCCCCTCTCGGGGCAACCCAGACG
TGCAACAAAGTGCATGCGTTCGCCCCAGACTACAGTAATACTCAACAATGATCGGCAGAACGCCATT
GTAGCCAAGATGGAAGACCCCTTGAGCAACAGGGCACCAGGATTCCCTGAAAAATGTCATTAGCAACGCT
GTGCTGGGCGTCGCAGAACACCATTGTGGTGAAGGTGCCGGGCAAGAAGACAGCCACCAGGAGAC
GGGAGAGCGGCTCGGAGGCCAGCGACTCTGTGTCCAGCTGTGGCAGGGGGCAGTCAGAGCATCGGG
AGCAACGTCACGCTCATACCCTGAACCTGGAAGAGGACTACCCCAATGGCACCTGGCTGGGCGACGAG
AACAAACCCGAGATGCGGGTACGCTGCGCCATCATCCCCTCCGACATGCTGCACATCAGCACCAACTGC
CGCACGGCCGAGAAGATGGCGCTCACGCTGCTGGACTACCTCTCCACCCGAGGTGCAGGCTGTGTCC
AACCTCTCGGGGCAAGGCAAGCAGCGGAAGAAGCAGCTGGACCCGCTCACCATCTACGGCATCCGGTGT
CACCTTTTCTATAAATTTGGCATCACAGAATCCGACTGGTACCGAATCAAGCAGAGCATCGACTCCAAG
TGCCGCACGGCGTGGCGGCGCAAGCAGCGGGGCCAGAGCCTGGCGGTCAAGAGCTTCTCGCGGAGAACG
CCCAACTCGTCTCTACTGCCCTTCAGAGCCGATGATGAGCACCCACCTCTGCCAGCGAGCTCCCG
CAGCCACAGCCGAGCCGAGGCCCTGCACTACGCGCTGGCCAACGCACAGCAGGTGCAGATCCACCAG
ATCGGAGAAGACGACAGGTGCAAGTAGGACACCTCCACATCGCCAGGTGCCGAGGGGAGCAAGTC
CAGATCACGCAGACAGCGAGGGCAACCTCCAGATCCATCACGTGGGGCAGGACGGTCAGGTGCTGCAG
GGTGACACAGTGTGCGCGTGGCCCTCTCGGACCCCGCGCGGGGGCGTGGATGGGTGCCACTCCAG
GGCAGGCAGATCCAGTTTCACTGACGCTGCGCCAGTGAGTGACCACAGCCCGGGGCACAGACG
GCCGAGCCCTGCAGCCACGCTACAGCCGGAGATGCAGCTCGAGCACGGGGCCATCCAGATTCAGTGA
```

```
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN:

NM_001173539

Insert Size:	1518 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001173539.1</u>
RefSeq Size:	2345 bp
RefSeq ORF:	1518 bp
Locus ID:	54971
UniProt ID:	<u>Q8N9N5</u>
Cytogenetics:	16q24.2
MW:	55.2 kDa
Gene Summary:	<p>This gene encodes a protein that binds to matrix attachment regions. The protein forms a complex with p53 and negatively regulates p53 transcription, and functions as a tumor suppressor and cell cycle regulator. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Mar 2010]</p> <p>Transcript Variant: This variant (3) represents use of an alternate promoter and 5' UTR, uses one alternate in-frame splices site, and lacks one alternate in-frame exon in the coding region, compared to variant 7. The resulting isoform (c) differs at two internal regions, compared to isoform g.</p>