

Product datasheet for **SC328852**

TTC39B (NM_001168342) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TTC39B (NM_001168342) Human Untagged Clone
Tag:	Tag Free
Symbol:	TTC39B
Synonyms:	C9orf52
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC328852 representing NM_001168342.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTACCATGCCTTGGGCTACAGTACCATTGTGGTGTTCAGGCTGTCTGACCTTCGAGCAACAGGAC
ATCCAGAACGGCATTCTGCCATGAAGGACGCTTTACAAACCTGCCAAAAATACAGGAAAAATACACA
GTTGTAGAATCTTCTCAAGTCTTCTTTCTAGAGGATCCCTGGAGCAACTGAGTGAAGAGGAAATGCAT
GCTGAAATCTGTTATGCCGAGTGTCTCCTACAGAAAGCTGCACTCACGTTTGTGCAGGATGAAAATATG
ATCAACTTCATCAAAGGTGGACTCAAATTAGAACAAAGTTACCAAATATATAAAGAATGTCTTCTATC
CTGCATGAAATTCAGAAGAACAACTTCAGCAGGAATTCTTCTATGAGTTTGAAGGGGTGTCAAGCTT
GGCAGTGGGGCTTTAATTTGATGCTGTCCCTTTACCAGCAAGGATTATCAGACTACTAGAATTTATT
GGATTTCTGGGAATAGGGAGTTGGGCCTCTGCAGTTACGTGAAGGTGCTTCAGGAAGGAGCATGAGG
TCTGCACTGTGCTGCTTAACCATCTTAGCTTTTCATACTACATCTCCCTAATACTTGGTACAGGAGAA
GTGAATGTTGCGGAAGCAGAGCGTCTCCTGGCACCTTCTCCAGCAGTTTCAAATGGCTCACTCGTG
TTGTTTTATCATGCCGAATAGAGTTGCTGAAAGGGAATCTTGAAGAGGCACAAGAAGTATTTCAAAA
TGCATTTCAAGTCAAGAAGAATGAAAACAGTTTACCATCTCTGCTACTGGGAGCTAATGTGGATTAAT
GTTTTCCAACAAAAGTGGATGCAGGCATATTACTATTCAGATCTGCTTTGCAAAGAGAGTAAATGGTCC
AAGGCACTTATGTGTTCTTGAAGCAGCAATTTTGTAGTATGCTTCCAGAGGAGGATGTAGTAGCAACG
AATGAGAATGTGGTAACTTTATTCAGACAGGTGGCAGCTTGAAGCAGAGAATTGCCGGGAAATCTATC
CCTACTGAGAAGTTGCTGTGAGGAAGGCTCGGCGTTACTCCGCCTCCTACCTGCACCTGTGAAGCTC
ATCTTACCTGCCCTGAAATGATGTATGTCTGGAATGGTTTTCAATAGTGAGCAAAAAGAAAAGACCTT
TCTGAAAATCTGTTAGTAACTGTTGAAAAGCTGAGGCAGCTTTACAAAGTCAAAATTTTAAACAGCTTC
TCTGTTGATGATGAGTGCTTAGTGAAGTTACTTAAAGGATGTTGCCTCAAGAAGTTACAGCGGCCCTTG
CAAGCTGAACTATGTTACAATCATGTTGTGAAAAGTAAAAGCTACTGAAGTATGACCACTACCTAGTG
CCGTTCACTCTATTTGAATTGGCATCTTTGTATAAAAGCCAGGGGGAAATTGACAAGGCCATAAAGTTC
CTAGAACTGCAAGGAACAACACTACAAAGATTACTCCCTGGAGTCCAGACTACACTTCAGAATTCAGGCA
GCTTTCACCTCTGGAGAAAACCTTCTCAGATTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001168342

Insert Size:	1554 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001168342.1</u>
RefSeq Size:	3151 bp
RefSeq ORF:	1554 bp
Locus ID:	158219
UniProt ID:	<u>Q5VTQ0</u>
Cytogenetics:	9p22.3
MW:	59.5 kDa
Gene Summary:	<p>Regulates high density lipoprotein (HDL) cholesterol metabolism by promoting the ubiquitination and degradation of the oxysterols receptors LXR (NR1H2 and NR1H3). [UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (5) differs in the 5' UTR and has multiple coding region differences, compared to variant 1. These differences cause translation initiation at an in-frame downstream AUG and result an isoform (5) with a shorter N-terminus, compared to isoform 1.</p>