

Product datasheet for **SC328873**

RAI2 (NM_001172739) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RAI2 (NM_001172739) Human Untagged Clone
Tag:	Tag Free
Symbol:	RAI2
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF:	>NCBI ORF sequence for NM_001172739, the custom clone sequence may differ by one or more nucleotides <pre> ATGGACGACCTGCAGTCCCAGAACCTCTCCATGGACATGACTGACTCCCCTCCTGCCTTG GCTAATAACAGACTGGAGAATGGCATGGCTCAGCTGATCACCACCGAGGCTGGAACATC AACTCCACTGACCTGGTAAAGAAGGCCCTGGTGACCGTGCCAGCCCCATCATTCTGAAC CCCCCTGCCGAGTCTCAGAGTGGCATGGCTCTGAAGGTGGCGGCACTGTGTTGCAGCCC CTGTGCCCTCGGGAGAGCCCAAGTGGTATGCCCATTCACATGCAGGTGGAGGGAAGCTCC GCACCAGAGCTCAATCCGAATGGCAATGCCACCTACGTATGACCACCCAGGGCCCCGTG CAACTGCCCGTGGTGCTGGAGCAGCACGTCTTTCAGCACCTCAACTCCCCCTCTGGTCCTG CCGCAGGAGGCCCATGCTCCTCCAGTACCATCCACAACAACCTCTTCCAGGGAGCGGAG GACCCCGAGGCCAGCCCCAGCTCCTGGACCTGAGGATCCCCAGCCAGCCGAGGAGCCC ACTTTGCCATTTGAAGCTGTGCTCCAGAATTTGTTTCCCTCCCAGGGCACTCTCGGGCCC CCACCCTGTCAGCCTCCTCCTGGCTATGCCCCTGTGCCCCACAGCCTTTTAGCTCCCCC TTGTCCCCCTGGTCCCACCAGCCACCCTCTTGGTGCCGTATCCTGTAATCGTCCCCTTG CCTGTGCCAGTCCCTATTCCTATCCCCATCCCGGTGCCTCAGAGTTCTGAATCCAAGTTC AGCTCCAGTTTCCCAAGCCACCATCTTCTTCGGCTGCACCCCTTAAAGGCACCCAG ACCCCTCTGGAAGAGATGAAGTGAAGCCCTTTGACATCTCCAGCCTAAGGAGTACTTC CAGCTCAGCCGCCACACGGTCATTAAGATGGGAAGTGAGAACGAGGCCCTGGATCTCTCC ATGAAGTCAGTGCCCTGGCTCAAGGCTGGTGAAGTCAGTCCCCCAATCTTCCAGGAAGAT GCAGCCCTAGACCTGTCAGTGGCAGCCACCGAAATCCGAGCCTCCCCCTGAGACACTG TATGACAGTGGTGCATCAGTGGCAGCTCAGGTACACAGTGTGAGAGAACTTCCCAGT GGCATGGAAATTTCTTTTGGCCCTGCCACGTCCCATGAGGCCCCAGCCATGATGGATAGT CACATCAGCAGCAGTGTGCTGCTACCGAGATGCTCAGCCAGCCCAACCACCCAGCGGC GAAGTCAAGGCTGAAAATAACATTGAGATGGTGGGCGAGTCCCAGGCGGCCAAGGTCATT GTCTCTGTGCAAGATGCTGTGCCTACCATATTCTGTGGCAAGATCAAAGGCCTCTCAGGG GTGTCCACCAAAACTTCTCCTTCAAAAGAGAAGACTCCGTGCTCAGGGCTATGACATC AACAGCCAAGGGGAAGAGTCCATGGGAAATGCAGAGCCCTTAGGAAACCCATCAAAAAAC CGGAGCATAAAGTTAAAGAAAGTGAAGTCCCAGGAAATACACATGCTCCCAATCAAAAAA CAACGGCTGGCCACCTTTTTTCCAAGAAAGTAA </pre>
Restriction Sites:	Please inquire
ACCN:	NM_001172739
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001172739.1, NP_001166210.1</u>
RefSeq Size:	2423 bp
RefSeq ORF:	1593 bp
Locus ID:	10742
UniProt ID:	<u>Q9Y5P3</u>
Cytogenetics:	Xp22.13
Gene Summary:	Retinoic acid plays a critical role in development, cellular growth, and differentiation. The specific function of this retinoic acid-induced gene has not yet been determined but it may play a role in development. The chromosomal location of this gene designates it to be a candidate for diseases such as Nance-Horan syndrome, sensorineural deafness, non-specific X-linked cognitive disability, oral-facial-digital syndrome, and Fried syndrome. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Feb 2010] Transcript Variant: This variant (3) differs in the 5' UTR compared to variant 1. Variants 1, 2 and 3 encode the same isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.