

Product datasheet for **SC328880**

HPSE2 (NM_001166244) Human Untagged Clone

Product data:

| | |
|---------------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | HPSE2 (NM_001166244) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | HPSE2 |
| Synonyms: | HPA2; HPR2; UFS; UFS1 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |



[View online »](#)

Fully Sequenced ORF: >SC328880 representing NM_001166244.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGGGTGTCTTGTGCCTTCCCTGAAGCCATGCCCTCCAGCAACTCCCCCCCCCGCGTGCCTAGCC
CCGGGGGCTCTACTTGGCTCTGTTGCTCCATCTCCCTTTCTCCAGGCTGGAGACAGGAGACCC
TTGCCTGTAGACAGAGCTGCAGGTTTGAAGGAAAAGACCCTGATTCTACTTGATGTGAGCACCAAGAAC
CCAGTCAGGACAGTCAATGAGAACTTCTCTCTCTGCAGCTGGATCCGTCCATCATTATGATGGCTGG
CTCGATTTCTAAGCTCCAAGCGCTTGGTGACCCTGGCCCGGGGACTTTCGCCCGCTTTCTGCGCTTC
GGGGGCAAAGGACCGACTTCTGCAGTTCAGAACCTGAGGAACCCGGCGAAAAGCCGGGGGGCCCG
GGCCCGGATTACTATCTCAAAAATGAGGATGACATTGTTTGAAGTATGTTGCCTTAGATAAACAG
AAAGGTGCAAGATTGCCAGCACCTGATGTTATGCTGGAGCTCAAAGGGAGAAGGCAGCTCAGATG
CATCTGGTCTTCTAAAGGAGCAATTCTCAATACTTACAGTAATCTCATATTAACAGAGCCAAATAAC
TATCGGACCATGCATGGCCGGGAGTAAATGGCAGCCAGTTGGGAAAGGATTACATCCAGCTGAAGAGC
CTGTTGCAGCCATCCGGATTTATTCCAGAGCCAGCTTATATGGCCCTAATATTGGCGCGCCGAGGAAG
AATGTCATCGCCCTCTAGATGGATTATGAAGGTGGCAGGAAGTACAGTAGATGCAATACCTGGCAA
CATTGCTACATTGATGGCCGGTGGTCAAGGTGATGGACTTCTGAAAATCGCCTGTTAGACACACTC
TCTGACCAGATTAGGAAAATTCAGAAAGTGTTAATACATACACTCCAGGAAAGAAGATTGGCTTGAA
GGTGTGGTGACCACCTCAGCTGGAGGCACAAACAATCTATCCGATTCCTATGCTGCAGGATTCCTATGG
TTGAACACTTTAGGAATGCTGGCAATCAGGGCATTGATGTCGTGATACGGCACTATTTTTGACCAT
GGATAAATCACCTCGTGGACCAGAATTTAACCATTACCAGACTACTGGCTCTCTCTCTACAAG
CGCCTGATCGGCCCAAAGTCTTGGCTGTGCATGTGGCTGGCTCCAGCGGAAGCCACGGCCTGGCCGA
GTGATCCGGGACAAACTAAGATTTATGCTCACTGCACAACCACCAACCAACTACGTTTCGTGGG
TCCATTACACTTTTTATCATCAACTTGCATCGATCAAGAAAAGAAAATCAAGCTGGCTGGACTCTCAGA
GACAAGCTGGTTCACCAGTACCTGCTGCAGCCCTATGGGCAGGAGGGCTAAAGTCCAAGTCAGTGCAA
CTGAATGGCCAGCCCTTAGTGATGGTGGACGACGGGACCCTCCAGAAATTGAAGCCCCGCCCTTCGG
GCCGGCCGACATTGGTCATCCCTCCAGTACCATGGGCTTTTATGTGGTCAAGAATGTCAATGCTTTG
GCCTGCCCTACCGATTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



| | |
|-------------------------------|--|
| ACCN: | NM_001166244 |
| Insert Size: | 1605 bp |
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). |
| OTI Annotation: | This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>NM_001166244.1</u> |
| RefSeq Size: | 4175 bp |
| RefSeq ORF: | 1605 bp |
| Locus ID: | 60495 |
| UniProt ID: | <u>Q8WWQ2</u> |
| Cytogenetics: | 10q24.2 |
| Protein Pathways: | Glycosaminoglycan degradation, Metabolic pathways |
| MW: | 60.1 kDa |
| Gene Summary: | <p>This gene encodes a heparanase enzyme. The encoded protein is a endoglycosidase that degrades heparin sulfate proteoglycans located on the extracellular matrix and cell surface. This protein may be involved in biological processes involving remodeling of the extracellular matrix including angiogenesis and tumor progression. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Oct 2009]</p> <p>Transcript Variant: This variant (2) lacks an in-frame exon in the coding region, compared to variant 1. The encoded isoform (2) is shorter than isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p> |