

Product datasheet for **SC328886**

C2 (NM_001178063) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C2 (NM_001178063) Human Untagged Clone
Tag:	Tag Free
Symbol:	C2
Synonyms:	ARMD14; CO2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC328886 representing NM_001178063.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCGTGCCTGTGCATCAGGGAGACCTGCAGCTCCGAGCTGGGGTTCAGCAGGAACTGGTCCCGCAGG
AAGGCTGGCCACTGCCCAACCCAGGCATTTCACTGGGCGCAGTGCAGCAGGCTTCCGCTTTGGTCAT
GGGACAAAGTCCGCTATCGCTGCTCCTCGAATCTTGTGCTCACGGGGTCTTCGGAGCGGAGTGCCAG
GGCAACGGGGTCTGGAGTGAACGGAGCCCATCTGCCGCCAACCTACTCTTATGACTTCCCTGAGGAC
GTGGCCCTGCCCTGGGCACTTCTTCTCCACATGCTTGGGGCCACCAATCCCACCCAGAAGACAAAG
GATCATGAAAAATGGAAGTGGGACTAACACCTATGCGGCCTTAACAGTGTCTATCTCATGATGAACAAC
CAAATGCGACTCCTCGCATGGAACGATGGCCTGGCAGGAAATCCGACATGCCATCATCCTTCTGACA
GATGGAAGTCCAATATGGGTGGCTCTCCCAAGACAGCTGTTGACCATATCAGAGAGATCCTGAACATC
AACCAGAAGAGGAATGACTATCTGGACATCTATGCCATCGGGTGGGCAAGCTGGATGTGGACTGGAGA
GAACTGAATGAGCTAGGGTCCAAGAAGGATGGTGAGAGGCATGCCTTATTCTGCAGGACACAAAGCT
CTGACCAGGTCTTTGAACATATGCTGGATGTCTCCAAGCTCACAGACACCATCTGCGGGTGGGGAAC
ATGTCAGCAAACGCCTCTGACCAGGAGAGGACACCCTGGCATGTCACTATTAAGCCCAAGAGCCAAGAG
ACCTGCCGGGGGCCCTCATCTCCGACCAATGGGTCTGACAGCAGCTCATTGCTTCCGGATGGCAAC
GACCACTCCCTGTGGAGGGTCAATGTGGGAGACCCCAATCCCAGTGGGGCAAAGAATTCCTTATTGAG
AAGCGGTGATCTCCCCAGGGTTTGTATGCTTTGCCAAAAAGAACCAGGGAACTCTGGAGTTCTATGGT
GATGACATAGCTCTGCTGAAGCTGGCCAGAAAGTAAAGATGTCCACCATGCCAGGCCATCTGCCTT
CCCTGCACGATGGAGCCAATCTGGCTCTGCGGAGACCTCAAGGCAGCACCTGTAGGACCATGAGAAT
GAACTGCTGAACAAACAGAGTGTCTCTGCTCATTGTCGCTTGAATGGGAGCAAACATGAACATTAAC
CTTAAGATGGGAGTGGAGTGGACAAGCTGTGCCGAGTTGTCTCCAAGAAAAAACCATGTTCCCAAC
TTGACAGATGTCAGGGAGGTGGTACAGACCAGTTCTATGCAAGTGGGACCCAGGAGGATGAGAGTCCC
TGCAAGGGAGAATCTGGGGGAGCAGTTTTCTTGAGCGGAGATTAGGTTTTTTCAGGTGGGTCTGGTG
AGCTGGGGTCTTTACAACCCCTGCCTTGGCTCTGCTGACAAAACTCCCGCAAAGGGCCCTCGTAGC
AAGGTCCCGCCGCCACGAGACTTTCACATCAATCTTCCGCATGCAGCCCTGGCTGAGGCAGCACCTG
GGGATGTCTGAATTTTTTACCCTCTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites:

Sgfl-Mlul

Plasmid Map:



ACCN:	NM_001178063
Insert Size:	1617 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001178063.2</u>
RefSeq Size:	2012 bp
RefSeq ORF:	1617 bp
Locus ID:	717
UniProt ID:	<u>P06681</u>
Cytogenetics:	6p21.33
Protein Families:	Druggable Genome, Protease, Secreted Protein
Protein Pathways:	Complement and coagulation cascades, Systemic lupus erythematosus
MW:	60.3 kDa

Gene Summary:

Component C2 is a serum glycoprotein that functions as part of the classical pathway of the complement system. Activated C1 cleaves C2 into C2a and C2b. The serine proteinase C2a then combines with complement factor 4b to create the C3 or C5 convertase. Deficiency of C2 has been reported to associated with certain autoimmune diseases and SNPs in this gene have been associated with altered susceptibility to age-related macular degeneration. This gene localizes within the class III region of the MHC on the short arm of chromosome 6. Alternative splicing results in multiple transcript variants encoding distinct isoforms. Additional transcript variants have been described in publications but their full-length sequence has not been determined.[provided by RefSeq, Mar 2009]

Transcript Variant: This variant (3) has an alternate 5' sequence and lacks an in-frame internal segment, as compared to variant 1. The resulting isoform (3) is shorter; it has a distinct N-terminus and lacks an internal segment, as compared to isoform 1.