

Product datasheet for **SC328896**

PIP5K3 (PIKFYVE) (NM_001178000) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PIP5K3 (PIKFYVE) (NM_001178000) Human Untagged Clone
Tag:	Tag Free
Symbol:	PIKFYVE
Synonyms:	CFD; FAB1; HEL37; PIP5K; PIP5K3; ZFYVE29
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC328896 representing NM_001178000.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCCACAGATGATAAGACGTCCCAACACTGGACTCTGCTAATGATTTGCCTCGATCTCCTACTAGT
CCTTCTCATCTCACACACTTTAAACCTTTGACTCCTGATCAAGATGAGCCCCCTTTTAAATCAGCTTAT
AGTTCTTTTGTAAATCTCTTTTCGTTTTAACAAAGAGAGAGCAGAAGGAGGCCAGGAGAACAGCAGCCT
TTGAGTGGAAAGTTGGACCAGCCCTCAGCTCCCTTCGAGGACACAGTCTGTTAGGTCACCCACACCTTAT
AAAAAGCAGCTTAATGAGGAACTCCAGCGGCGCTCTTCAGCATTAGACACAAGAAGGAAAGCAGAACCT
ACCTTTGGAGGTCATGACCCTCGTACAGCTGTTGACTTCGAAGCCTCAGCACAGTATTAACACGCCTC
AAGGAAATCATGGAGGGGAAAAGCCAGGATAGTGACCTGAAACAATACTGGATGCCAGATAGCCAATGT
AAAGAGTGCTATGACTGTAGTGAATAATCAACCTTTAGGCGCAGACACCATTGCCGACTATGTGGG
CAGATTTTCTGCAGTCGTTGCTGTAATCAAGAAATCCCTGGAAAATTTATGGGCTATACAGGAGACCTC
CGAGCTTGACATATTGTAGAAAAATAGCCTTAAGTTATGCTCATTCCACAGACAGTAATTCTATTGGG
GAAGACTTGAATGCTCTTTCAGATTCTGCTTGCTCTGTGCTGTGCTTGAATCCAAGTGAACCCCGAACA
CCTGTTGGGAGTAGGAAAGCCAGCCGTAACATATTTTTAGAGGATGATTTGGCCTGGCAAAGTTTGATT
CATCCAGATTCCTCAAATACTCCTCTTTCAACAAGACTTGATCTGTGCAAGAGGATGCTGGGAAATCT
CCTGCTCGAAATAGATCAGCCAGCATTACTAACCTGCTACTGGATAGATCTGTTCTCCTATGGTACCT
TCATATGAGACATCTGTCACTCCCAGGCTAACCGAACATATGTTAGGACAGAGACCCTGAGGATGAA
CGCAAAATTTCTTGGACAGTGTGCACTAAAAGACCTGTGGAAAAAATCTGCCATCACAGCAGTGGA
ATGGAGTTTCAGGATCACCGCTACTGGTTGAGAACGCATCCCACTGCATTGTAGGAAAGGAATTAGTC
AACTGGCTAATCCGAAATGGGCATATTGCCACAAGGCCACAAGCTATAGCAATTGGACAAGCAATGGTT
GATGGACGTTGGCTGGATTGTGTTAGTCATCACGACCAGCTTTTCAGAGATGAGTATGCGCTGTATAGA
CCTACTGCAGAGTACAGAATTTCTGAGACGCCTTCTCCGACAGTGACTCAGTGAACCTCCGTGGAAGGA
CACTCTGAGCCATCCTGGTTTAAAGACATAAAGTTTGTGACAGTGACACAGAACAGATAGCTGAAGAA
GGTGACGATAATTTGGCTAATTCTGCCAGTCTAGCAAGCGCACATCAGTCAGCAGTTTCCAGTCCACA
GTGGACAGTGACTCAGCCGTTCTATCAGCCTGAACGTGGAGCTGGACAACGTGAACCTCCATATCAAG
AAGCCCTCCAAGTACCCACATGTGCCCCCTCACCTGCTGACCAAAAAGGTAGGAGGTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN:	NM_001178000
Insert Size:	1647 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001178000.1</u>
RefSeq Size:	2421 bp
RefSeq ORF:	1647 bp
Locus ID:	200576
UniProt ID:	<u>Q9Y2I7</u>
Cytogenetics:	2q34
Protein Families:	Druggable Genome
Protein Pathways:	Endocytosis, Fc gamma R-mediated phagocytosis, Inositol phosphate metabolism, Metabolic pathways, Phosphatidylinositol signaling system, Regulation of actin cytoskeleton
MW:	61.6 kDa

Gene Summary:

Phosphorylated derivatives of phosphatidylinositol (PtdIns) regulate cytoskeletal functions, membrane trafficking, and receptor signaling by recruiting protein complexes to cell- and endosomal-membranes. Humans have multiple PtdIns proteins that differ by the degree and position of phosphorylation of the inositol ring. This gene encodes an enzyme (PIKfyve; also known as phosphatidylinositol-3-phosphate 5-kinase type III or PIPKIII) that phosphorylates the D-5 position in PtdIns and phosphatidylinositol-3-phosphate (PtdIns3P) to make PtdIns5P and PtdIns(3,5)biphosphate. The D-5 position also can be phosphorylated by type I PtdIns4P-5-kinases (PIP5Ks) that are encoded by distinct genes and preferentially phosphorylate D-4 phosphorylated PtdIns. In contrast, PIKfyve preferentially phosphorylates D-3 phosphorylated PtdIns. In addition to being a lipid kinase, PIKfyve also has protein kinase activity. PIKfyve regulates endomembrane homeostasis and plays a role in the biogenesis of endosome carrier vesicles from early endosomes. Mutations in this gene cause corneal fleck dystrophy (CFD); an autosomal dominant disorder characterized by numerous small white flecks present in all layers of the corneal stroma. Histologically, these flecks appear to be keratocytes distended with lipid and mucopolysaccharide filled intracytoplasmic vacuoles. Alternative splicing results in multiple transcript variants encoding distinct isoforms.[provided by RefSeq, May 2010]

Transcript Variant: This variant (4) lacks all of the 3' coding region and UTR, compared to variant 2. The encoded protein (isoform 4) has an identical N-terminus but a short novel C-terminus, compared to isoform 2.