

Product datasheet for **SC328903**

PPP2R1B (NM_001177562) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PPP2R1B (NM_001177562) Human Untagged Clone
Tag:	Tag Free
Symbol:	PPP2R1B
Synonyms:	PP2A-Abeta; PR65B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC328903 representing NM_001177562.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGGGCGCATCAGAGCTCGGACCGGCCAGGAGCAGCGGGTGGAGATGGAGATGATTTCGTATAC
CCGATCGCGGTTTTAATCGACGAGCTCCGCAATGAAGACGTGCAGCTCCGACTCAACAGTATTAAGAAG
TTATCAACAATTGCCCTAGCACTTGGAGTAGAAAAGGACCCGAAGTGAATTGTTGCCATTTCTACAGAT
ACAATTTATGATGAAGATGAGGTAATAGTCTTGGTGGAGCAGCTGGGAAATTTCACTGGCCTAGTG
GGAGGTCCTGACTTTGCCACTGTCTGCTGCCTCTTTGGAAAATCTGGCAACTGTGGAAGAGACTGTT
GTTCTGACAAAGGCTGTGGAGTCCCTGAGACAGATCTCCAGGAGCATACTCCTGTTGCTCTGGAAGCT
TATTTGTACCTCTGGTGAACGCTTAGCAAGTGGGGATTGGTTCACCTCTGCACATCTGCATGTGGT
TTGTTGAGCGTTTGTATCCAGGGCATCAAATGCTGTTAAAGCAGAAATCAGACAGCAATTCGGTTCC
TTGTGCTCAGATGACACACCAATGGTACGACGTGCTGCTTCCAAATGGGTGAATTTGCAAAAGTT
TTGAATTAGACAGTGTAAAAGTAAAATTTGTTCCACTGTTCACTAGTCTAGCTTCAAGTGAACAGGAT
TCAGTGGCCTCCTTGTGTGGAAGCTTGTGTCAGTATTGCCAGTATTGTCTCAGGATGACCTTGAG
ACTTTGGTGTATGCCTACACTTCGACAAGCAGCAGAAGATAAATCTTGGCGCTTCGCTATATGGTGGCT
GACAGATTTTCAGAGCTCCAGAAAGCCATGGGTCTAAAATCACCTAAATGACCTCATCCCCGCCCTTT
CAGAACCTACTTAAAGACTGTGAAGCTGAAGTCCGGGCAGCTGCTGCCACAAAGTAAAAGAAGTGGT
GAGAAGTTGCCATTGAAGATAGAGAGACCATAATTATGAATCAAATCTGCCTTATATAAAGTGTCTCT
GACGTTTCGTTTGAATATCATCTCCAATTTGGATTGTGTAATGAAGTATTGGAATCCGTCAGCTCTCT
CAGTCTCTCCTTCTGCCATAGTGGAGCTGGCAGAAGATGCCAAATGGAGGGTCCGCCTGGCCATCATT
GAGTATATGCCGCTGCTGGCAGCCAGCTGGGTGTGGAATCTTTGATGAAAAGCTGAATCTTTATGT
ATGGCTTGGCTCGTGGACCATGTATACGCCATCCGAGAAGCTGCCACCAACAACCTCATGAAACTAGTT
CAGAAGTTTGGTACAGAGTGGGCCAAAATACTATTGTTCCAAAAGTGTAGTAATGGCAAATGATCCT
AATTACTTGCATAGAATGACCACTTTATTCTGCATTAATGCACTGTCTGAGGCCTGTGGTCAAGAAATA
ACTACTAAGCAAATGCTGCCATCGTATTAATAAATGGCAGGAGACCAAGTAGCAAATGTTTCGCTTCAAT
GTGGCCAAATCTCTACAAAAGATTGGACCAATCTAGATACCAATGCTTTACAGGGAGAAGTGAAGCCA
GTACTACAGAAGTTAGGTCAAGATGAAGACATGGATGTCAAATACTTTGCACAGGAAGCTATAAGTGTT
CTTGCAATTGGCATAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001177562

Insert Size: 1671 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001177562.1](#)

RefSeq Size: 5470 bp

RefSeq ORF: 1671 bp

Locus ID: 5519

UniProt ID: [P30154](#)

Cytogenetics: 11q23.1

Protein Families:	Druggable Genome, Phosphatase, Transcription Factors
Protein Pathways:	Long-term depression, Oocyte meiosis, TGF-beta signaling pathway, Tight junction, Wnt signaling pathway
MW:	61.3 kDa
Gene Summary:	<p>This gene encodes a constant regulatory subunit of protein phosphatase 2. Protein phosphatase 2 is one of the four major Ser/Thr phosphatases, and it is implicated in the negative control of cell growth and division. It consists of a common heteromeric core enzyme, which is composed of a catalytic subunit and a constant regulatory subunit, that associates with a variety of regulatory subunits. The constant regulatory subunit A serves as a scaffolding molecule to coordinate the assembly of the catalytic subunit and a variable regulatory B subunit. This gene encodes a beta isoform of the constant regulatory subunit A. Mutations in this gene have been associated with some lung and colon cancers. Alternatively spliced transcript variants have been described. [provided by RefSeq, Apr 2010]</p> <p>Transcript Variant: This variant (4) lacks an in-frame exon in the coding region compared to variant 1. The resulting isoform (d) is shorter than isoform a. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>