

Product datasheet for **SC328904**

MFSD5 (NM_001170790) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MFSD5 (NM_001170790) Human Untagged Clone
Tag:	Tag Free
Symbol:	MFSD5
Synonyms:	hsMOT2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC328904 representing NM_001170790.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGAGGTCAGTGCAGTGGGAGAGGAGGCAAACCAAGGAGACCTTCTCGGCTTCGGGCAGCGGTGCTGT
TCCTTAGGAAGCTCCAACCTTGGCTCTACCGGACTGCGGGCATCTTTCCGAGCACTGCACGGAAGAA
AAGAAAATTGGCACTTATTTAGGGTGGCCCTCAGGCATCCAAGAGCCAGAGGAGCAGGTGCACGGG
ATGGCAGCGAGCGAGGCCAGGCCAGGAGAAGCAACCGCAGCCTCCTATGGTCGCCCTTCCCGACCCAC
CCCAGGGCTCAGGAGTTCGCGGCCCGGGTCTCCGGGGGCCACCATGCTGGTACTGCCTACCTTGCT
TTTGTAGGCTCCTGGCTCCTGCCTGGGGCTGGAAGTCAAGATGCCGGGCTAAACCCCTGGAAGG
GCCTGCAGCAATCCCTCCTTCTCGGTTTCAACTGGACTTCTATCAGGTCTACTTCTGGCCCTGGCA
GCTGATTGGCTCAGGCCCTACCTCTATAAACTCTACCAGCATTACTTCTGGAAGTCAAATT
GCCATCCTCTATGTCTGTGGCTTGCTCTACAGTCTCTTTGGCTAGTGGCTCCTCCCTTGTTGAT
TGCTGGGTGCGAAGAATTCTTGTGCTCTTCTCCCTGACTTACTCACTATGCTGCTTAACCAAATC
TCTCAAGACTACTTTGTGCTGCTAGTGGGGCAGCACTTGGTGGCTGTCCACAGCCCTGCTCTTCTCA
GCCTTCGAGGCCTGGTATATCCATGAGCAGTGAACGGCATGACTTCCCTGCTGAGTGGATCCCAGCT
ACCTTTGCTCGAGCTGCCTTCTGGAACCATGTGCTGGCTGTAGTGGCAGGTGTGGCAGCTGAGGCTGTA
GCCAGCTGGATAGGGCTGGGCCGTAGCGCCCTTTGTGGCTGCCATCCCTCTCCTGGCTCTGGCAGGG
GCCTTGGCCCTTCGAAACTGGGGGAGAACTATGACCGGCAGCGTGCCTTCTCAAGGACCTGTGCTGGA
GGCCTGCGCTGCCTCCTGTGGACCGCCGCGTGTGCTGTTGGGCACCATACAAGCTCTATTTGAGAGT
GTCATCTTCATCTTTGTCTTCTCTGGACACCTGTGCTGGACCCACACGGGGCCCTCTGGGCATTATC
TTCTCCAGCTTCATGGCAGCCAGCCTGCTTGGCTCTTCCCTGTACCGTATCGCCACCTCCAAGAGGTAC
CACCTTCAGCCCATGCACCTGCTGCCCTTGTGCTCATCGTCGTCTTCTCTCTCTTCATGTTGACT
TTCTCTACCAGCCAGGCCAGGAGAGTCCGGTGGAGTCTTCATAGCCTTCTACTTATTGAGTTGGCT
TGTGGATTATACTTTCCAGCATGAGCTTCTACGGAGAAAGGTGATCCCTGAGACAGAGCAGGCTGGT
GTACTCAACTGGTTCGGGTACCTCTGCACTCACTGGCTTGCTAGGGCTCCTTGTCTCCATGACAGT
GATCGAAAACAGGCACTCGAATATGTTGAGCATTGCTCTGCTGTGATGGTATGGCTCTGCTGGCA
GTGGTGGGACTCTCACCGTGGTAAGGCATGATGCTGAGCTGCGGGTACCTTACCTACTGAGGAGCCC
TATGCCCTGAGCTGTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCCGC
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:


ACCN: NM_001170790

Insert Size: 1674 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001170790.1](#)

RefSeq Size: 2054 bp

RefSeq ORF: 1674 bp

Locus ID: 84975

UniProt ID: [Q6N075](#)

Cytogenetics: 12q13.13

Protein Families: Transmembrane

MW: 61.5 kDa

Gene Summary: Mediates high-affinity intracellular uptake of the rare oligo-element molybdenum.
[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1).