

Product datasheet for **SC328920**

NOSTRIN (NM_001171631) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NOSTRIN (NM_001171631) Human Untagged Clone
Tag:	Tag Free
Symbol:	NOSTRIN
Synonyms:	DaIP2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC328920 representing NM_001171631.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGAGGGACCCACTGACAGATTGTCGTATAATAAAGTATACAAGAACCCTAAAGGAGTTTCTCAAAT
GGAGAGAATTTCTGCAAAACAGGTCACATCTGTTCTTCAGCAAAGGGCAAACCTGGAAATTAGCTATGCC
AAAGGACTTCAGAAACTGGCAAGCAAGCTGAGCAAAGCATTACAGAACACGAGAAAAAGTTGTGTTAGC
AGTGCCTGGGCCTGGGCCTCAGAGGGAATGAAATCCACAGCGGACCTGCATCAAAAACTGGCAAAGCA
ATTGAATTGGAAGCAATAAAACCGACTTATCAAGTCTAAATGTACAAGAGAAGAAGAGAAAACTCACTT
GACAATGAAGTTGAAAAGACAGCAAATCTGTCTTAGCAACTGGAATCAGCAAATTAAGGCCAAGAAG
AAATTAATGGTTAGTACCAAGAAACATGAAGCACTTTTCCAGCTTGTAGAAAAGCTCAAGCAATCTATG
ACTGAGAAGGAGAAGCGGAAGCTCCTCAATAAACTGACAAAATCAACTGAAAAGTTGGAAAAGGAAGAT
GAAAATTACTACCAAAAAACATGGCGGGTTATTCTACCAGACTGAAATGGGAAAACACACTAGAGAAC
TGCTACCAGGTAACCCACAGCATCTGCCTATATGCCTTTTGGGTCAAAGAGCCTGGGAAAAGTGTGTG
AGTGATTTAAGATACCAAGACACATTCTTCTCGAAATCTCCCTCCACTATGGTTTGGATATGACATA
GTCAAGAGGCTAATCATGAGACTTTGTTCACTGTCTTTCAGAGCATTCTGGAGCTGGAGAAGGAAAAGA
ATTCAACTTTTATGCAATAACTTAAACCGTACAGCCAACATATTTCTCTTTTGGCCAAACCCTGACC
ACATGCCACACGCAGATTCACTGTGCCATCAGCAAAGATTGACATTGAAAAAGATATCCAGGCTGTAATG
GAAGAACTGCAATTTTATCTACAGAAAACAACTCTGAGTTCTGTTAACGGATTACTTTGAAGAAGAT
CCTAACAGTGCAATGGATAAAGAGAGACGAAAGCTTTACTAAAACAAAATTATTGAGACTGCAGAGA
GACATTGAAAAGCCTCAAAGACAAGGAAGGCCTGGAACGAATGCTTAAACGTACTCCAGCACCTCC
TCCTTCTCTGATGCAAAGAGCCAGAAAAGACACAGCAGCGTTAATGGATGAGAACAATTTGAACTAGAC
CTTTTGGAAAGCGAACTCTACAACTGTCATCAATGTTAGCAGAACTTGAGCAAAGACCTCAACCCAGC
CATCCTTGATGTAATTCCATCTTCAGGTGGAGGAAAAGGAGCATACTCATAGCTATGTGAAAATATCT
CGGCCTTTTAAATGAAGAGATTAGAGAATATTGTGAGCAAGGCATCTTCTGGTGGGCAGAGCAATCCA
GGTTCTTCAACTCCAGCCCTGGTGCAGCCAGCTCAGCAGCAGACTTTGCAAGGCCTTGATTCTTTT
CAAGCCAGGCAAGATGATGAGTTGAATTTGAAAAGGGTGACATTGTGATTATACACGAGAAAAAGAA
GGAGGATGGTGGTTTGGATCTTTGAATGGGAAAAGGCCATTTCTGCCGCTTATGTGGAGGAGTTA
CCTTCAAATGCTGGCAACACAGCTACAAAGGCA TAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:


ACCN: NM_001171631

Insert Size: 1692 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001171631.1](#)

RefSeq Size: 2819 bp

RefSeq ORF: 1692 bp

Locus ID: 115677

UniProt ID: [Q8IIV9](#)

Cytogenetics: 2q24.3

MW: 64.4 kDa

Gene Summary: Nitric oxide (NO) is a potent mediator in biologic processes such as neurotransmission, inflammatory response, and vascular homeostasis. NOSTRIN binds the enzyme responsible for NO production, endothelial NO synthase (ENOS; MIM 163729), and triggers the translocation of ENOS from the plasma membrane to vesicle-like subcellular structures, thereby attenuating ENOS-dependent NO production.[supplied by OMIM, Apr 2004]
Transcript Variant: This variant (4) has additional 5' UTR exons and an internal coding exon, as compared to variant 2. The resulting isoform (4) is the longest, and has an additional segment in the middle region, as compared to isoform 2. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.