

Product datasheet for **SC328941**

SCMH1 (NM_001172220) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SCMH1 (NM_001172220) Human Untagged Clone
Tag:	Tag Free
Symbol:	SCMH1
Synonyms:	Scml3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC328941 representing NM_001172220.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGAAATTGGAAGCACAGGACCCAGGAACACCACATCCACCTGTATTGCCACAGTAGTTGGACTGACA
GGTGCCCGCCTTCGCCTGCGCCTTGATGGGAGCGACAACAAAAATGACTTCTGGCGGCTGGTTGACTCA
GCTGAAATCCAGCCTATTGGGAACTGTGAAAAGAATGGGGGTATGCTACAGCCACCTCTTGATTTCGG
CTGAATGCGTCTTCTTGCCCATGTTCTTTTGAAGACGCTAAATGGAGCAGAGATGGCTCCCATCAGG
ATTTTCCACAAGGAGCCACCATCGCCTTCCACAACCTTCTTCAAATGGGAATGAAGCTAGAAGCTGTG
GACAGGAAGAACCCTCATTTTCAATTTGCCAGCCACTATTGGGGAGGTTTCGGGGCTCAGAGGTGCTTGTG
ACTTTTGTGATGGGTGGCGAGGGGCTTTGACTACTGGTGCCGCTTCGACTCCCGAGACATCTTCCCTGTG
GGCTGGTGTTCCTTGACTGGAGACAACCTGCAGCCTCCTGGCACCAGTTGTGATTCCAAAGAATCCC
TATCCTGCCTCCGATGTGAATACTGAGAAGCCAGCATCCACAGCAGCACCAAACTGTCTTGAACAT
CAACCAGGGCAGAGGGGGCGTAAACCAGGAAAGAAGCGGGGCCGGACACCCAAGACCCTAATTTCCCAT
CCCATCTCTGCCCATCCAAGACAGCTGAACCTTTGAAATTCCTCAAGAAGAGAGGTCCCAACCTGGC
AGCAAGAGGAAACCTCGGACTTTGCTGAACCCACCACCTGCCTCACCACAACCCAGCACTCCTGAACCG
GATACCAGCACTGTACCCAGGATGCTGCCACCATCCCCAGCTCAGCCATGCAGGCCCAACAGTTTGT
ATCTACTTGAACAAGAATGGCAGCACAGGCCCCCACTTAGATAAGAAGAAGTCCAGCAACTCCCTGAC
CATTTTGGACCAGCCCGTGCCTCTGTGGTGTGTCAGCAGGCTGTCCAGGCCTGTATCGACTGTGCTTAT
CACCAGAAAACCGTCTTCAGCTTCTCAAGCAAGGCCATGGTGGTGAAGTTATCTCAGCCGTGTTTAC
CGGGAACAGCATACCCTCAACCTCCAGCAGTCAACAGCATCACCTACGTCCTCCGCTTCTGGAGAAA
CTCTGCCACAACCTTCGTAGTGACAATCTGTTTGGCAACCAGCCCTTTACACAGACTCACTTGTCACTC
ACTGCCATAGAGTACAGCCACAGCCACGACAGGTACCTACCAGGTGAAACCTTTGCTCCTGGGAATAGT
CTGGCCCGCTCCTTGGAACCACACTCAGACTCAATGGACTCTGCCTCAAATCCCACCAACCTTGTGAGC
ACCTCCCAGGCAACCGCCCTTCTTTTATCCTGTGGCCTCCCACCAAGCACTGCCTCAGCTGTGCGC
AGGCTATGCTCCAGGGGTTCGGACCGATACCTGGAGAGCCGCGATGCCTCTCGACTGAGTGGCCGGGAC
CCCTCCTCATGGACAGTCGAGGATGTGATGCAGTTTGTCCGGGAAGCTGATCCTCAGCTTGGACCCAC
GCTGACCTGTTTCGAAACACGAGATCGATGGCAAGGCCCTGCTGCTGCTGCGCAGTGACATGATGATG
AAGTACATGGGCCTGAAGCTGGGCCTGCACTCAAGCTCTCTACCACATTGACCGGCTGAAGCAGGGC
AAGTTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001172220

Insert Size: 1734 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001172220.1](#)

RefSeq Size: 3084 bp

RefSeq ORF: 1734 bp

Locus ID: 22955

UniProt ID: [Q96GD3](#)

Cytogenetics: 1p34.2

Protein Families: Transcription Factors

MW: 63.9 kDa

Gene Summary: Associates with Polycomb group (PcG) multiprotein complexes; the complex class is required to maintain the transcriptionally repressive state of some genes.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (5) lacks two alternate exons in the 5' UTR, compared to variant 3. Both variants 3 and 5 encode the same isoform (c).