

## Product datasheet for **SC328951**

### **FAM130A2 (CSRNP3) (NM\_001172173) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	FAM130A2 (CSRNP3) (NM_001172173) Human Untagged Clone
Tag:	Tag Free
Symbol:	CSRNP3
Synonyms:	FAM130A2; PPP1R73; TAIP-2; TAIP2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC328951 representing NM\_001172173.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGAGTGAATTTAAAGAGGAAGTTGAAGAAGTTGACGGCTCCTCACCTGCTCCTGTGTGAGGGAA
TCAGATGATGAAGTTTCCAGCAGTGAAAGTGCTGACAGTGGGGACAGTGTCAATCCATCCACTTCTAGT
CATTTTACCCCTTCTCCATTCTCAAAGGGAGAAACGACTGAGGACAAAGAATGTACATTTTAGTTGT
GTCACCGTGTACTACTTACCAGGAGGCAAGGCTTACAAGTGTGCCAGTCAAGGGGAAGCACCTG
GGGATGTCAGCCGCCATAACAGCGTGCACCAGTACACTTTGGCGAGTTTGCAAGGGAGCAGGAGAGG
CTCCACCGGGAGATGTTGAGAGAACACCTTAGGGAGGAAAAGCTGAACTCCTTAAACTAAAGATGACT
AAGAATGGCACAGTAGAATCAGAAGAAGCCAGCACTTTACACTGGATGACATTTCTGATGATGACATT
GACCTGGACAACACAGAGGTAGATGAGTACTTCTTCTACAACCTTTGCCAACAAAAACAAGAGCT
CTGCTGCGTGCCTCTGGAGTGAAAAGATTGACGTGGAAGAAAAGCACGAACTCCGAGCCATCCGCCTC
TCACGAGAGGACTGTGGCTGTGACTGCCGAGTGTCTGTGATCCAGACACGTGCACCTGCAGCCTGGCT
GGCATTAAAGTGCCAGGTGGATCGTATGTCTTCCCATGCCGCTGCACTAAAGAAGGATGTAGTAACACA
GCAGGTAGAATTGAATTAATCCTATCCGTGTTCCGACTCACTTTTTGCACACAATAATGAACTTGAA
CTGGAGAAAAACCGAGAGCAGCAATCCCACGCTGAATGGCTGCCACAGTGAGATAAGTGTCCACAGT
AGTTCTATGGGCCTGTGCTCACTCCGTAGAATATTCATCGCAGACAGTTTTGAGATTGAACTGAG
CCCCAGGCTGCAGTGTGCACCTGCAGTCGGCTGAAGAATTAGATTGCCAAGGAGAGGAGGAGGAAGAA
GAGGAGGATGGGAGCAGCTTTTGCAGCGGAGTACAGATTCTAGCACGCAAAGCTTGGCACCTAGTGAG
TCAGACGAGGAGGAGGAGGAAGAAGAAGAGGAAAGAGGAGGAGGATGACGATGATGACAAAGGAGAT
GGCTTCGTGGAAGTTTTGGGCACCCATGCCGAAGTTGTCCCTTCTCCTTCAAGTTCTTTATTCTGAT
GGCACCGCGTTACGAAAGCCATGCAAAGAATGCTTCTTTTTATGCCAACTCTTCAACTCTGTATTAC
CAAATAGATAGCCACATTCCAGGAATCCAAATCAGATCTCTGAGAATTTCTGAAAGAGACACTGTC
AAAAATGGTACCTTTTCTGCTGGTGCCTTACACCATGACCCCGGAGCAATTCGTTGACTATGCCCGACAA
GCAGAAGAGGCCATGGTGCCTCCCACTACCCAGCTGCCAACCCCTCTGTAATCGTTTGTGCTCCTCT
TCCGAAAATGATAGCGGTGTGCCCTGCAATAGTTTATATCCTGAACACAGGTCCAATCACCTCAAGTG
GAATTTCACTCATACTTGAAGGCCCTCCCAAGAAGGTTTGTCTCTGCATTGAATGGTGACAGTCAC
ATTTAGAGCATCTGCTGAAAATCTTTGAGCCTTGCAAGAAAGAGCATATTGCATGAAGAGTGCATC
AAATCACCGTGGTTGAGACAGTCCCTGTTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_001172173

**Insert Size:** 1758 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001172173.1](#)

**RefSeq Size:** 11788 bp

**RefSeq ORF:** 1758 bp

**Locus ID:** 80034

**UniProt ID:** [Q8WYN3](#)

**Cytogenetics:** 2q24.3

**MW:** 64.9 kDa

**Gene Summary:** Binds to the consensus sequence 5'-AGAGTG-3' and has transcriptional activator activity. Plays a role in apoptosis (By similarity).[UniProtKB/Swiss-Prot Function]  
Transcript Variant: This variant (1) represents the longer transcript. Both variants 1 and 2 encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.