

Product datasheet for **SC328956**

DMT1 (SLC11A2) (NM_001174125) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DMT1 (SLC11A2) (NM_001174125) Human Untagged Clone
Tag:	Tag Free
Symbol:	SLC11A2
Synonyms:	AHMIO1; DCT1; DMT1; NRAMP2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC328956 representing NM_001174125.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGGAAGAAGCAGCTGAAGACGGAGGCAGCTCCACACTGTGAACTAAAATCATATTCTAAGAACTCA
GCCACTCAGGTATCCACCATGGTGTGGTCTGAACAGAAGATGTCAGATGACAGTGTCTTCTGGAGAT
CATGGGGAGTCTGCCAGTCTTGGTAACATCAACCCTGCCTATAGTAATCCCTCTCTTTCACAGTCCCT
GGGACTCAGAGGAGTACTTCGCCACTTACTTTAATGAGAAGATCTCCATTCTGAGGAGGAGTACTCT
TGTTTTAGCTTTTCGTAACCTCTGGGCTTTCACCGACCAGGTTTTCTTATGAGCATTGCCTACCTGGAT
CCAGGAAATATTGAATCCGATTTGCAGTCTGGAGCAGTGGCTGGATTTAAGTTGCTCTGGATCCTCTG
TTGGCCACCCTTGTGGGCTGCTGCTCCAGCGCTTGCAGCTAGACTGGGAGTGGTACTGGGCTGCAT
CTTGCTGAAGTATGTCACCGTCAGTATCCCAAGTCCCACGAGTCATCCTGTGGCTGATGGTGGAGTTG
GCTATCATCGGCTCAGACATGCAAGAAGTCATTGGCTCAGCCATTGCTATCAATCTTCTGTCTGTAGGA
AGAATTCCTCTGTGGGGTGGCGTTCTCATCACCATTGCAGATACTTTTGTATTTCTTCTTGGACAAA
TATGGCTTGGCGAAGCTAGAAGCATTTTTTGGCTTCTCATCACTATTATGGCCCTCACATTTGGATAT
GAGTATGTTACAGTGAACCCAGCCAGAGCCAGGTAAGGGCATGTTTCGTACCATCCTGTTCCAGGC
TGTGCACTCCACAGATTGAACAGGCTGTGGGCATCGTGGGAGCTGTCATCATGCCACACAACATGTAC
CTGCATTCTGCCTTAGTCAAGTCTAGACAGGTAACCCGGAACAATAAGCAGGAAGTTCGAGAAGCCAAT
AAGTACTTTTTTCATTGAATCCTGCATTGCACTCTTTGTTTCCTTCATCATCAATGTCTTTGTGTCTCA
GTCTTTGCTGAAGCATTTTTTGGGAAAACCAACGAGCAGGTGGTTGAAGTCTGTACAAATACCAGCAGT
CCTCATGCTGGCCTCTTTCCTAAAGATAACTCGACACTGGCTGTGGACATCTACAAGGGGGTGTGTG
CTGGGATGTTACTTTGGCCTGCTGCACCTACATTTGGCAGTGGGGATCCTGGCTGCAGGACAGAGC
TCCACCATGACAGGAACCTATTCTGGCCAGTTTGTGTCATGGAGGGATTCTGAACCTAAAGTGGTCACGC
TTTGCCCGAGTGGTCTGACTCGCTCTATTGCCATCATCCCCACTCTGTTGTTGCTGTCTTCCAAGAT
GTAGAGCATCTAACAGGGATGAATGACTTTCTGAATGTTCTACAGAGCTTACAGCTTCCCTTTGCTCTC
ATACCCATCCTCACATTTACGAGCTTGGCGCCAGTAATGAGTACTTTGCCAATGGACTAGGCTGGCGG
ATTGCAGGAGGAATCTTGGTCTTATCATCTGTTCCATCAATATGACTTTGTAGTGGTTTATGTCCGG
GACCTAGGGCATGTGGCATTATATGTGGTGGCTGCTGTGGTCAGCGTGGCTTATCTGGGCTTTGTGTTT
TACTTGGGTTGGCAATGTTTGATTGCACTGGGCATGTCCTTCTGGACTGTGGGCATACGGTAAGCATC
TCTAAAGGCCTGCTGACAGAAGAAGCCACCCGTGGCTACGTTAAATAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001174125

Insert Size: 1773 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001174125.1](#)

RefSeq Size: 4156 bp

RefSeq ORF: 1773 bp

Locus ID: 4891

UniProt ID: [P49281](#)

Cytogenetics: 12q13.12

Protein Families: Transmembrane

Protein Pathways: Lysosome

MW: 64.7 kDa

Gene Summary: This gene encodes a member of the solute carrier family 11 protein family. The product of this gene transports divalent metals and is involved in iron absorption. Mutations in this gene are associated with hypochromic microcytic anemia with iron overload. A related solute carrier family 11 protein gene is located on chromosome 2. Multiple transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Apr 2010]
Transcript Variant: This variant (1) encodes the longest isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.