

Product datasheet for **SC328957**

ODF2L (NM_001184766) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ODF2L (NM_001184766) Human Untagged Clone
Tag:	Tag Free
Symbol:	ODF2L
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC328957 representing NM_001184766.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGC
ATGGAGAAGGCTGTAAATGATGGAAGTCATTGAGAAGAACTCTTTGCCATCTTAAACTATATCAGAG
AAAGAAGATTTACCACGGTGTACCAAGTAAAGTCATCTCAGCTGCCTGAAGCAGGACATTCTAAATGAA
AAGACTGAATTGGAAGCAACACTTAAGGAAGCGGAGTTGGTAACCTATTCTGTAGAATTGCTTTTGCCA
CTATTTAAGGACACCATTTGAAAAGATTAATTTTAAAAATGCGAATCTTTCTGCATTGAATTTGAAGATT
TCTGAACAGAAGGAAATATTAATAAAAGAATTAGACACCTTCAAAAGTGTAAGCTTTAGAACAT
CTTCTTAGAAAGAGAGACTACAAACAAACAGGAGACAATTTATCCAGCATGTTGTAGAAAATCTAACT
GATAATGAAAGTAAAACTAATCTTAAGAAGAAGGTATTTGAAAAGGAGGCCCATATCCAAGAACTT
TCTTGTGTTGTTTCAGAGTGAAAAGGCAATACTTTGAAAGCAAACCGTTTTTCCCAATCAGTAAAAGTA
GTACATGAACGACTACAGATCCAAATTCATAAAAGGGAAGCCGAGAACGATAAGTTGAAAGAATATGTA
AAGAGCTTAGAAACCAAGATAGCCAAGTGAACCTGCAATCAAGAATGAATAAGAATGAGGCTATAGTG
ATGAAAGAAGCAAGTAGGCAAAAACTGTAGCTTTAAAAAAGGCATCTAAAGTTTACAAACAAAGGCTT
GACCATTTTACAGGAGCTATTGAAAAGCTTACTTCCCAATTAGAGATCAGGAAGCCAAGTTGTCTGAA
ACAATTTTCAGCTTCCAATGCCTGGAAAAGTCATTATGAGAAAATTGTGATAGAAAAAACGAATTGGAA
GTACAGATTGAAACAATGAAAAAGCAATCATTAACTTTTGAAGATCTGAAGAAAATGGAAGACCAT
GGAAAAAATTCATGTGAAGAAATCTTAGAAAAGTTCACTCAATTGAATATGAAATGAACTCTGAAT
CTTGAGAATACAAAATTAAGACTACACTTGCTGCTTTGAAGGATGAAGTTGTATCTGTTGAAAATGAA
CTCTCAGAATTGCAAGAAGTAGAAAAAACAGAAAACCTTATTGAAATGTATAAACTCAGGTACAA
AAGTTGCAAGAAGCAGCTGAAATAGTAAAAAGCAGATGTGAAAATTTGCTACATAAAAAATAACCAATA
ACAAAAACCAAAAAATAAAAAATGTAGAGAAGGTTGATGGAATCACAACTTCTGACAAAGCTTTCTCTG
GAAGAGGAAAACTGTCTTATTCAGCTGAAGTGTAACCTTCAACAAAAATTAGAACAGATGGACGCA
GAAAATAAAGAGCTTGAGAAGAAGCTGGCAACCAAGAAGATGTCTTAAGCACAGCAATCTTAAGTTT
AAAGAGAAATCTGCAGAGTATACAGCATTGGCCAGACAAGTGAAGCTGCTTTAGAAGAAGGAAGACAA
AAGGTTGCTGAAGAAATAGAGAAATGTCATCTAGAGAGAGTGCTTTACAAATTAATATAGATCTG
GAACTGAGCTTAGAAAGAAAAATGAAGAACAAATCACTTGTTGCAAAATGAACAGTAAAGCACAA
CATCAGGAAGTGTTTGAAGAAGTACAAATAGTCTTGAGAAGTCGAAAAATCAAAATGAGAGTATT
AAGAATTATTTACAGTTTCTTAAACTTCTTATGTAACAATGTTGAATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001184766

Insert Size: 1776 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001184766.1</u>
RefSeq Size:	4081 bp
RefSeq ORF:	1776 bp
Locus ID:	57489
UniProt ID:	<u>Q9ULJ1</u>
Cytogenetics:	1p22.3
MW:	68.6 kDa
Gene Summary:	Acts as a suppressor of ciliogenesis, specifically, the initiation of ciliogenesis. [UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (4) has multiple differences in the coding region, compared to variant 1, one of which results in a translational frameshift. The resulting protein (isoform d) has a distinct C-terminus and is shorter than isoform a. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.