

Product datasheet for **SC328974**

SLC26A11 (NM_001166348) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLC26A11 (NM_001166348) Human Untagged Clone
Tag:	Tag Free
Symbol:	SLC26A11
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001166348, the custom clone sequence may differ by one or more nucleotides

```

ATGCCTTCTTCGGTGACGGCGCTGGGTCAGGCCAGGTCCTCTGGCCCCGGGATGGCCCCG
AGCGCCTGCTGCTGCCCTGCGGCCCTGCAGAGGAGGCTGCCCATCTGGCGTGGCTG
CCCAGCTACTCCCTGCAGTGGCTGAAGATGGATTCGTCGCCGGCCTCTCAGTTGGCCTC
ACTGCCATTCCCCAGGCGCTGGCCTATGCTGAAGTGGCTGGACTCCCGCCCCAGTATGGC
CTCTACTCTGCCCTTCATGGGCTGCTTCGTGATTTCTTCTCTGGGCACCTCCCGGGATGTG
ACTCTGGGCCCCACCGCATTATGTCCCTCCTGGTCTCCTTCTACACCTTCCATGAGCCC
GCCTACGCTGTGCTGCTGGCCTTCTGTCCGGCTGCATCCAGCTGGCCATGGGGGTCCTG
CGTTTGGGGTTCTGCTGGACTTCATTTCTACCCCGTCATTAAGGCTTACCTCTGCT
GCTGCCGTACCATCGGCTTTGGACAGATCAAGAACCTGCTGGGACTACAGAACATCCCC
AGGCCGTTCTTCTGCAGGTGTACCACACCTTCTCAGGATTGCAGAGACCAGGGTAGGT
GACGCCGCTCCTGGGGCTGGTCTGCATGCTGCTGCTGCTGGTGTGAAGCTGATGCGGGAC
CACGTGCCTCCCGTCCACCCGAGATGCCCCCTGGTGTGCGGCTCAGCCGTGGGCTGGTC
TGGGCTGCCACGACAGCTCGCAACGCCCTGGTGGTCTCCTTCGCAGCCCTGGTTGCGTAC
TCCTTCGAGGTGACTGGATACCAGCCTTTCATCCTAACAGGGGAGACAGCTGAGGGGCTC
CCTCCAGTCCGGATCCCGCCCTTCTCAGTGACCACAGCCAACGGGACGATCTCCTTACC
GAGATGGTGCAGGACATGGGAGCCGGGCTGGCCGTGGTGGCCCTGATGGGCTCCTGGAG
AGCATTGCGGTGGCCAAAGCCTTCGCATCTCAGAATAATTACCGCATCGATGCCAACAG
GAGCTGCTGGCCATCGGTCTCACCAACATGTTGGGCTCCCTCGTCTCCTCCTACCCGGTC
ACAGGCAGCTTTGGACGGACAGCCGTGAACGCTCAGTCGGGGGTGTGCACCCCGGGGGG
GGCCTGGTGACGGGAGTCTGGTGTGCTGTCTCTGGACTACCTGACCTCACTGTTCTAC
TACATCCCCAAGTCTGCCCTGGCTGCCGTCATCATATGGCCGTGGCCCCGCTGTTTCGAC
ACCAAGATCTTTCAGGACGCTCTGGCGTGTAAAGAGGCTGGACCTGCTGCCCTGTGCGTG
ACCTTCTGCTGTGCTTCTGGGAGGTGCAGTACGGCATCCTGGCCGGGGCCCTGGTGTCT
CTGCTCATGCTCTGCACTCTGCAGCCAGGCCTGAGACCAAGGTGTGAGAGGGGCCGGTT
CTGGTCTGCAGCCGGCCAGCGGCTGTCTTCCCTGCCATGGAGGCTCTGCGGGAGGAG
ATCCTAAGCCGGGCCCTGGAAGTGTCCCCGCCACGCTGCCTGGTCTGGAGTGCACCCAT
GTCTGCAGCATCGACTACACTGTGGTGTGGGACTCGGCGAGCTCCTCCAGGACTCCAG
AAGCAGGGCGTCGCCCTGGCCTTTGTGGGCTGCAGGTCCCGTCTCCTGTCTCTGCTG
TCCGCTGACCTGAAGGGTTCCAGTACTTCTTACCCTGGAAGAAGCAGAGAAGCACCTG
AGGCAGGAGCCAGGGACCCAGCCCTACAACATCAGAGAAGACTCCATTCTGGACAAAAG
GTTGCCCTGCTCAAGGCATAA
    
```

Restriction Sites: Please inquire

ACCN: NM_001166348

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001166348.1](#), [NP_001159820.1](#)

RefSeq Size: 2860 bp

RefSeq ORF: 1821 bp

Locus ID: 284129

UniProt ID: [Q86WA9](#)

Cytogenetics: 17q25.3

Protein Families: Transmembrane

Gene Summary: This gene encodes a member of the solute linked carrier 26 family of anion exchangers. Members of this family of proteins are essential for numerous cellular functions including homeostasis and intracellular electrolyte balance. The encoded protein is a sodium independent sulfate transporter that is sensitive to the anion exchanger inhibitor 4,4'-diisothiocyanostilbene-2,2'-disulfonic acid. Alternate splicing results in multiple transcript variants.[provided by RefSeq, Oct 2009]
Transcript Variant: This variant (3) differs in the 5' UTR compared to variant 1. Variants 1, 2, 3 and 4 encode the same protein.