

Product datasheet for **SC329008**

ZNF185 (NM_001178110) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF185 (NM_001178110) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF185
Synonyms:	SCCELL
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC329008 representing NM_001178110.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGTATCTCAGCTCTTGGAGGCCACCAAAGGGAAGCCTCTGCCACCAGCGGAGGAGGAGCGCAAT
AACGTTCTCAAGCAGATGAAAGTGCGAACCACGCTGAAGGGGGACAAGAGCTGGATTACCAAGCAGGAT
GAATCGGAGGGTCGCACCATAGAGCTGCCCTCAGGCCGGAGTCGCACCACATCCTTTTCATCAGCTGGG
GAGGTTCCGAAGCCTAGGCCTCCGAGCACAAGGGCTCCCACTGGCTACATCATCCGGGGAGTGTTACCC
AAGCCATAGACAGCTCTCCAGCCCCAGCAGAGTTCCCAAGGCCAACGGGACTCCAAAAAGTGCT
GCCAGTCTGGTGAGAACAGCTAACGCTGGTCTCCCGCCCCCTCCTCTGGCTACAAGATGACCACT
GAGGATTACAAGAAGCTGGCACCCCTACAATATCAGGCGCAGCTCTACATCAGGGGACACCGAGGAGGAG
GAGGAGGAGGAGGTGGTCCATTCTCTCAGATGAACAGAAACGGAGGTGAGAGGCTGCAAGCGGTGTT
CTGAGGAGGACAGCTCCCGGGAGCACTCTACGTCCTGTGAGCGCCAGAAGAGCACTGGTCTACC
CAGGAGACACAGGCACCGTTTATCGCGAAGAGGGTGGAGGTGGTGAAGAGGACGGGCCCTCTGAGAAG
AGCCAGGACCCACCTGCTCTGGCAAGATCCACTCTGGCTCAAACAGGTCTCCCGAGGCAACAAAGAC
AAGGAGGCCCCCTGCTCCAGAGAGCTCCAGAGGGACTTGGCTGGTGAGGAGGCTTTCAGGGCCCCAAC
ACAGATGCTGCAAGGTCAAGTGACAGTTGAGTGATGGCAATGTGGGATCCGGAGCCACGGGCTCCCGG
CCTGAAGGCTTGGCTGCAGTAGACATCGGCTCCGAGAGAGGAAGTCCAGTGCCACTTCACTCTGCT
GTCCCTGCTGATAGGAAGAGCAACAGCAGCAGCCAGGAGGATGCAAAGGCAGACCCAAAGGGGGCC
TTGGCTGATTATGAGGGGAAGGATGTGGCCACCAGGGTCGGAGAGGCCTGGCAGGAGAGGCTGGAGCT
CCAAGAGGTGGCCAAGGAGACCCAGCTGTACCCGCTCAGCAACCTGCAGATCCCAGCACCCAGAGCGG
CAGAGCAGCCCCAGCGGATCTGAGCAACTTGTGACAGCAGAGAGTTGTGGCAGCAGCGTGTGACTGAT
TTTGAGGGGAAGGATGTGGCCACCAAGGTTCGGAGAGGCCTGGCAGGACAGGCCTGGAGCCCCAAGAGGT
GGCCAAGGAGACCCAGCTGTACCACTCAGCAACCTGCAGATCCCAGTACCCAGAACAGCAGAACAGC
CCCAGCGGATCTGAGCAATTCTGTCAGACGAGAGCTGCACCAGCAGGGTGAGGAGCCCCCTGAGCTGC
ATGGTCACTGTTACTGTCACTGCCACATCTGAGCAGCCTCACATTTATATCCAGCCCCGCAAGTGAA
TTGGACTCCAGCTCTACCACCAAAGGGATTCTCTTCTGTAAGGAGTACGTGAATGCTAGTGAAGTGTCT
TCTGGGAAGCCAGTATCTGCACGCTATAGCAACGTCAGCAGCATTGAGGACTCATTGCCATGGAGAAG
AAGCCTCCATGTGGCAGCACTCCATACTCTGAGAGGACAACCTGGAGGGATCTGTACTTACTGCAACCGT
GAGATCCGAGACTGTCCAAAGATTACCCTAGAACATCTTGGTATCTGCTGCCATGAATATTGCTTTAAG
TGTGGGATTTGCAGTAAACCGATGGGCGATCTCCTGGATCAGATCTTCATTCACCGTGACACCATTAC
TGTGGGAATGCTATGAGAAGCTCTTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:



ACCN: NM_001178110

Insert Size: 1893 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001178110.1](#)

RefSeq Size: 4167 bp

RefSeq ORF: 1893 bp

Locus ID: 7739

UniProt ID: [O15231](#)

Cytogenetics: Xq28

MW: 67.3 kDa

Gene Summary: Zinc-finger proteins bind nucleic acids and play important roles in various cellular functions, including cell proliferation, differentiation, and apoptosis. This gene encodes a LIM-domain zinc finger protein. The LIM domain is composed of two contiguous zinc finger domains, separated by a two-amino acid residue hydrophobic linker. The LIM domain mediates protein:protein interactions. Multiple alternatively spliced transcript variants encoding different isoforms have been identified.[provided by RefSeq, May 2010]
Transcript Variant: This variant (6) lacks two consecutive in-frame exons in one region and an in-frame exon in another region in the CDS, as compared to variant 1. The resulting isoform (6) lacks two internal segments, as compared to isoform 1.