

## Product datasheet for **SC329016**

### PREPL (NM\_001171613) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	PREPL (NM_001171613) Human Untagged Clone
Tag:	Tag Free
Symbol:	PREPL
Synonyms:	CMS22
Mammalian Cell Selection:	None
Vector:	<u><a href="#">pCMV6-XL5</a></u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >NCBI ORF sequence for NM\_001171613, the custom clone sequence may differ by one or more nucleotides

```

ATGGATGCATTTGAAAAAGTGAGAACAAAATTAGAAACACAGCCACAAGAAGAATATGAA
ATCATCAATGTGGAAGTTAAACATGGTGGTTTTGTTTATTACCAAGAAGTTGTTGCTTG
GTTTCGTTCCAAAGATGAAGAAGCAGACAATGATAATTATGAAGTTTTATTCAATTTGGAG
GAACTTAAGTTAGACCAGCCCTTCATTGATTGTATCAGAGTTGCTCCAGATGAAAAATAT
GTGGCTGCCAAGATAAGAACTGAAGATTCTGAAGCATCTACCTGTGTAATTATAAAAGCTC
AGCGATCAGCCCGTAATGGAAGCTTCTTTCCCGAATGTGTCCAGTTTTGAATGGGTAAAG
GACGAGGAAGATGAAGATGTTTTATTCTACACCTTCCAGAGGAACCTTCGCTGTCATGAC
GTATATCGAGCCACTTTTGGTGATAACAAACGTAATGAACGCTTTTACACAGAAAAAGAC
CCAAGCTACTTTGTTTCCTTTATCTTACAAAAGACAGTCGTTTCCTCACCATAAATATT
ATGAACAAGACTACTTCTGAAGTGTGGTTGATAGATGGCCTGAGCCCTTGGGACCCACCA
GTACTTATCCAGAAGCGAATACATGGGGTCTTTACTATGTTGAACACAGAGATGATGAA
TTATACATTCTACTAATGTTGGAGAACCTACAGAATTTAAGCTAATGAGAACAGCGGCT
GATACCCCTGCAATTATGAATTGGGATTTATTTTTACAATGAAGAGAAATACAAAAGTG
ATAGACTTGGACATGTTTAAAGGATCACTGTGTTCTATTTCTGAAGCACAGCAATCTCCTT
TATGTTAATGTGATTGGTCTGGCTGATGATTCAGTTCGGTCTCTAAAGCTCCCTCCTTGG
GCCTGTGGATTATAATGGATACAAATTCTGACCCAAAGAACTGCCCTTTCAACTTTGC
TCTCCAATACGTCCCCAAAATATTACACATACAAGTTTGCAGAAGGCCAACTGTTTGAG
GAAACTGGGCATGAAGACCCAATCACAAAGACTAGTCGCGTTTTACGTCTAGAAGCCAAA
AGCAAGGATGGAAAATTAGTGCCAATGACTGTTTTCCACAAAAGTACTCTGAGGACTTG
CAGAAGAAACCTCTCTTGGTACATGTATATGGAGCTTATGGAATGGATTTGAAAAATGAAT
TTCAGGCTGAGAGCGGGTCTTGGTGGATGATGGATGATGATGATGATGATGATGATGATG
CGAGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGT
AATGGCCTTGCTGATTTAGAGGCTTGCATTAAGACGCTTCATGGCCAAGGCTTTTCTCAG
CCAAGTCTAACAAACCCTGACTGCTTTCAGTGTGGAGGGGTGCTTGCAGGAGCATTGTGT
AATTCTAATCCAGAGCTGGTGGAGCGGTGACTTTGGAGGCACCTTTCTTGGATGTTCTC
AACACCATGATGGACACTACACTTCTGACATTAGAAGAATTAGAAGAATGGGGGAAT
CCTTCATCTGATGAAAAACACAAGAACTACATAAAACGTTACTGTCCCTATCAAAATATT
AAACCTCAGCATTATCCTTCAATTCACATAACGGCATATGAAAACGATGAACGGTACCT
CTGAAAGGAATTGTAAGTTATACTGAGAACTCAAGGAAGCCATCGCGGAGCATGCTAAG
GACACAGGTGAAGGCTATCAGACCCCTAATATTATTCTAGATATTCAGCCTGGAGGCAAT
CATGTAATTGAGGATTCTCACAAAAGATTACAGCCCAAATTAATTCCTGTACGAGGAA
CTTGGACTTGACAGCACCAGTGTTCGAGGATCTTAAGAAATACCTGAAATTCTGA

```

- Restriction Sites:** Please inquire
- ACCN:** NM\_001171613
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
- OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001171613.1, NP_001165084.1</u>
<b>RefSeq Size:</b>	5907 bp
<b>RefSeq ORF:</b>	1917 bp
<b>Locus ID:</b>	9581
<b>UniProt ID:</b>	<u>Q4I6C6</u>
<b>Cytogenetics:</b>	2p21
<b>Protein Families:</b>	Druggable Genome, Protease
<b>Gene Summary:</b>	<p>The protein encoded by this gene belongs to the prolyl oligopeptidase subfamily of serine peptidases. Mutations in this gene have been associated with hypotonia-cystinuria syndrome, also known as the 2p21 deletion syndrome. Several alternatively spliced transcript variants encoding either the same or different isoforms have been described for this gene.[provided by RefSeq, Jan 2010]</p> <p>Transcript Variant: This variant (6) contains an alternate exon at the 5' end compared to variant 1, resulting in translation initiation from an in-frame downstream AUG and a shorter isoform (4) compared to isoform 1. Variants 6 and 7 encode the same isoform. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>