

Product datasheet for **SC329020**

SIGLEC10 (NM_001171156) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SIGLEC10 (NM_001171156) Human Untagged Clone
Tag:	Tag Free
Symbol:	SIGLEC10
Synonyms:	PRO940; SIGLEC-10; SLG2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC329020 representing NM_001171156.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGCTACTGCCACTGCTGCTGTCTCGCTGTGGCGGGTCCCAGGCTATGGATGGGAGATTCTGGATA
CGAGTGCAGGAGTCAGTGTGGTCCCGAGGGCCTGTGCATCTCTGTGCCCTGCTTTTCTCTACCCC
CGACAAGACTGGACAGGGTCTACCCAGCTTATGGCTACTGGTTCAAAGCAGTACTGAGACAACCAAG
GGTGTCTCTGTGGCCACAAACCACAGAGTCGAGAGGTGAAATGAGCACCCGGGCGGATTCCAGCTC
ACTGGGGATCCCAGCAAGGGAACTGCTCCTTGGTGATCAGAGACGCGCAGATGCAGGATGAGTCACAG
TACTTCTTTGGGTGGAGAGAGGAAGCTATGTGAGATATAATTTTCATGAACGATGGGTTCTTTCTAAAA
GTAACAGTGTCTCAGTTCACGCCAGACCCAGGACCACAACACCGACCTCACCTGCCATGTGGACTTC
TCCAGAAAGGGTGTGAGCGCACAGAGGACCGTCCGACTCCGTGTGGCCTATGCCCCAGAGACCTTGT
ATCAGCATTTTACGTGACAACACGCCAGCCCTGGAGCCCAGCCCAGGGAAATGTCCCATACCTGGAA
GCCAAAAGGCCAGTTCCTGCGGCTCCTCTGTGTGCTGACAGCCAGCCCCTGCCACACTGAGCTGG
GTCTGTGAGAACAGAGTCTCTCCTCGTCCCATCCCTGGGGCCCTAGACCCTGGGGCTGGAGCTGCC
GGGGTGAAGGCTGGGGATTGAGGGCGCTACACCTGCCAGCGGAGAACAGGCTTGGCTCCAGCAGCGA
GCCCTGGACCTCTCTGTGAGTATCCTCCAGAGAACCTGAGAGTGATGGTTTCCAAGCAACAGGACA
GTCCTGGAAAACCTTGGGAACGGCACGTCTCTCCAGTACTGGAGGGCCAAAGCCTGTGCCTGGTCTGT
GTCACACACAGCAGCCCCCAGCCAGGCTGAGCTGGACCCAGAGGGGACAGGTTCTGAGCCCCTCCAG
CCCTCAGACCCCGGGTCTGGAGCTGCCTCGGGTTCAAGTGGAGCACGAAGGAGAGTTCACCTGCCAC
GCTCGGCACCCACTGGGCTCCAGCACGTCTCTCAGCCTCTCCGTGCACTACTCCCCGAAGCTGCTG
GGCCCCCTCCTGCTCCTGGGAGGCTGAGGGTCTGCACTGCAGCTGCTCCTCCAGGCCAGCCCAGCCCC
TCTCTGCCCTGGTGGCTTGGGGAGGAGCTGCTGGAGGGGAACAGCAGCCAGGACTCCTTCGAGGTCACC
CCCAGCTCAGCCGGGCCCTGGGCCAACAGCTCCCTGAGCCTCCATGGAGGGCTCAGCTCCGGCCTCAGG
CTCCGCTGTGAGGCTGGAACGTCCATGGGGCCAGAGTGATCCATCCTGCAGTGCCAGATAAGAAG
GGACTCATCTCAACGGCATTCTCCAACGGAGCGTTTCTGGGAATCGGCATCACGGCTCTTCTTTCTC
TGCCTGGCCCTGATCATCATGAAGATTCTACCGAAGAGACGGACTCAGACAGAAACCCGAGGCCAGG
TTCTCCCGGCACAGCACGATCCTGGATTACATCAATGTGGTCCCGACGGCTGGCCCCCTGGCTCAGAAG
CGGAATCAGAAAGCCACACCAACAGTCTCGGACCCCTCTCCACCAGGTGCTCCCTCCCAGAATCA
AAGAAGAACCAGAAAAGCAGTATCAGTTGCCAGTTTCCAGAACCCAAATCATCCACTCAAGCCCCA
GAATCCCAGGAGAGCCAAGAGGAGCTCCATTATGCCACGCTCAACTTCCCAGGCGTCAGACCCAGGCT
GAGGCCGGATGCCAAAGGGCACCCAGGCGGATTATGCAGAAGTCAAGTTCCAATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001171156

Insert Size: 1920 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001171156.1](#)

RefSeq Size: 3226 bp

RefSeq ORF: 1920 bp

Locus ID: 89790

UniProt ID: [Q96LC7](#)

Cytogenetics:	19q13.41
Protein Families:	Druggable Genome, Transmembrane
MW:	70.3 kDa
Gene Summary:	<p>SIGLECs are members of the immunoglobulin superfamily that are expressed on the cell surface. Most SIGLECs have 1 or more cytoplasmic immune receptor tyrosine-based inhibitory motifs, or ITIMs. SIGLECs are typically expressed on cells of the innate immune system, with the exception of the B-cell expressed SIGLEC6 (MIM 604405).[supplied by OMIM, Jul 2002]</p> <p>Transcript Variant: This variant (2) lacks an in-frame segment in the CDS, as compared to variant 1. The resulting isoform (2) lacks an internal segment, as compared to isoform 1.</p>