

## Product datasheet for **SC329028**

### SCMH1 (NM\_001172219) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	SCMH1 (NM_001172219) Human Untagged Clone
Tag:	Tag Free
Symbol:	SCMH1
Synonyms:	Scml3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC329028 representing NM\_001172219.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCAGCCTAATGTGATAGACTGGAGTGATGTTAGAAAACACAAATATGGTCACCTATCAGAGTCTGCA
TCCCAATATCAAGAAGCTGCTGACATCCTGGATCTAGGTCAATTTACCTGGGACAAATACCTAAAAGAA
ACATGTTCAAGTCCAGCGCTGTCCATTGCTTCAAGCAGTCTACACACCTCCAAGCAACGAGTTCAAG
ATCAGTATGAAATTGGAAGCACAGGACCCAGGAACACCACATCCACCTGTATTGCCACAGTAGTTGGA
CTGACAGGTGCCCGCCTTCGCCTGCGCCTTGATGGGAGCGACAACAAAAATGACTTCTGGCGGCTGGTT
GACTCAGTGAAATCCAGCCTATTGGGAACTGTGAAAAGAATGGGGGTATGCTACAGCCACCTCTTGGA
TTTCGGCTGAATGCGTCTTCTTGCCCATGTTCCTTTTGAAGACGCTAAATGGAGCAGAGATGGCTCCC
ATCAGGATTTTCCACAAGGAGCCACCATCGCCTTCCACAACCTTCTTCAAATGGGAATGAAGCTAGAA
GCTGTGGACAGGAAGAACCCTCATTTCATTGCCAGCCACTATTGGGGAGGTTCCGGGCTCAGAGGTG
CTTGTCACCTTTGATGGGTGGCAGGGGCTTTGACTACTGGTGGCCTTCGACTCCCGAGACATCTTC
CCTGTGGGCTGGTGTTCCTTGACTGGAGACAACCTGCAGCTCCTGGCACCAAAGTTGTGATTCCAAAG
AATCCCTATCCTGCCTCCGATGTGAATACTGAGAAGCCAGCATCCACAGCAGCACCAAAACTGTCTTG
GAACATCAACCAGGGCAGAGGGGGCGTAAACCAGGAAAGAAGCGGGGCCGACCCAAAGACCCTAATT
TCCCATCCCATCTCTGCCCATCCAAGACAGCTGAACCTTTGAAATCCCAAAGAAGAGAGGTCCCAAA
CCTGGCAGCAAGAGGAAACCTCGGACTTTGCTGAACCCACCACCTGCCTCACCAACAACCAGCACTCT
GAACCGGATACCAGCACTGTACCCAGGATGCTGCCACCATCCCAGCTCAGCCATGCAGGCCCAACA
GTTTGTATCTACTTGAACAAGAATGGCAGCACAGGCCCCCACTTAGATAAGAAGAAGGTCCAGCAACTC
CCTGACCATTTTGGACCAGCCCGTGCCTCTGTGGTGTTCAGCAGGCTGTCCAGGCTGTATCGACTGT
GCTTATACCAGAAAACCGTCTTCAGCTTCTCAAGCAAGGCCATGGTGGTGAGGTTATCTCAGCGGTG
TTTGACCGGGAACAGCATACCCTCAACCTCCCAGCAGTCAACAGCATCACCTACGTCCTCCGCTTCTG
GAGAAACTCTGCCACAACCTTCGTAGTGACAATCTGTTTGGCAACCAGCCCTTTACACAGACTCACTTG
TCACTCACTGCCATAGAGTACAGCCACAGCCACGACAGGTACCTACCAGGTGAAACCTTTGTCTGGGG
AATAGTCTGGCCCGCTCCTTGAACCACTCAGACTCAATGGACTCTGCCTCAAATCCCACCAACCTT
GTCAGCACCTCCCAAAGGCACCGGCCCTTGCTTTCATCCTGTGGCCTCCCACCAAGCACTGCCTCAGCT
GTGCGCAGGCTATGCTCCAGGGGTCCGACCGATACCTGGAGAGCCGCGATGCCTCTCGACTGAGTGGC
CGGGACCCCTCCTCATGGACAGTCGAGGATGTGATGCAGTTTGTCCGGGAAGCTGATCCTCAGCTTGA
CCCCAGCTGACCTGTTTCGAAACACGAGATCGATGGCAAGGCCCTGCTGCTGCTGCGCAGTGACATG
ATGATGAAGTACATGGGCTGAAGCTGGGGCTGCACTCAAGCTCTCTACCACATTGACCGGCTGAAG
CAGGGCAAGTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_001172219

**Insert Size:** 1947 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001172219.1](#)

**RefSeq Size:** 3155 bp

**RefSeq ORF:** 1947 bp

**Locus ID:** 22955

**UniProt ID:** [Q96GD3](#)

**Cytogenetics:** 1p34.2

**Protein Families:** Transcription Factors

**MW:** 72.1 kDa

**Gene Summary:** Associates with Polycomb group (PcG) multiprotein complexes; the complex class is required to maintain the transcriptionally repressive state of some genes.[UniProtKB/Swiss-Prot Function]  
Transcript Variant: This variant (4) lacks an alternate exon and uses an upstream start codon, compared to variant 3. The resulting isoform (d) has a longer N-terminus, compared to isoform c.