

Product datasheet for **SC329031**

IGSF10 (NM_001178145) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	IGSF10 (NM_001178145) Human Untagged Clone
Tag:	Tag Free
Symbol:	IGSF10
Synonyms:	CMF608
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF:	>NCBI ORF sequence for NM_001178145, the custom clone sequence may differ by one or more nucleotides <pre> ATGGTGGAAAGCGTAAGGCTTGAGAACCAGCCCTGTGACCTAAGAGTGGGCAGCTGGATC CACGTCTACCCTAATGGATCCCTGTTTATTGGATCAGTAACAGAAAAAGACAGTGGTGTC TACTTGTGTGTGGCAAGAAACAAAATGGGGGATGATCTGATACTGATGCATGTTAGCCTA AGACTGAAACCTGCCAAAATTGACCACAAGCAGTATTTTAGAAAGCAAGTGCTCCATGGG AAAGATTTCCAAGTAGATTGCAAAGCTTCCGGCTCCCCAGTGCCAGAGATATCTTGGAGT TTGCCTGATGGAACCATGATCAACAATGCAATGCAAGCCGATGACAGTGGCCACAGGACT AGGAGATATACCTTTTTCAACAATGGAACTTTATACTTCAACAAAGTTGGGGTAGCGGAG GAAGGAGATTATACTTGCTATGCCCAGAACCCCTAGGGAAAGATGAAATGAAGGTCCAC TTAACAGTTATAACAGCTGCTCCCGGATAAGGCAGAGTAACAAAACCAACAAGAGAATC AAAGCTGGAGACACAGCTGTCTTGACTGTGAGGTCACTGGGGATCCCAAACCAAAAATA TTTTGGTTGCTGCCTTCCAATGACATGATTTCTTCTCCATTGATAGGTACACATTTTCT GCCAATGGGTCTTTGACCATCAACAAAGTGAACTGCTCGATTCTGGAGAGTACGTATGT GTAGCCCGAAATCCCAGTGGGGATGACACCAAAATGTACAACTGGATGTGGTCTCTAAA CCTCCATTAATCAATGGTCTGTATACAAACAGAACTGTTATTAAGCCACAGCTGTGAGA CATTCCAAAAACACTTTGACTGCAGAGCTGAAGGGACACCATCTCTGAAGTCATGTGG ATCATGCCAGACAATATTTTCTCACAGCCCCATACTATGGAAGCAGAATCACAGTCCAT AAAAATGGAACCTTGGAATTAGGAATGTGAGGCTTTCAGATTAGCCGACTTTTATCTGT GTGGCCCGAAATGAAGGTGGAGAGAGCGTGTGGTAGTACAGTTAGAAGTACTGGAATG CTGAGAAGACCGACATTTAGAAATCCATTTAATGAAAAATAGTTGCCAGCTGGGAAAG TCCACAGCATTGAATTGCTCTGTTGATGGTAACCCACCACCTGAAATAATCTGGATTTTA CCAAATGGCACACGATTTTCCAATGGACCACAAAGTTATCAGTATCTGATAGCAAGCAAT GGTCTTTTATCATTTTCAAAACAACTCGGGAGGATGCAGGAAAAATATCGCTGTGCACT AGGAATAAAGTTGGCTATATTGAGAAATTAGTCATATTAGAAATTGGCCAGAAGCCAGTT ATTCTTACCTATGCACCAGGGACAGTAAAAGGCATCAGTGGAGAATCTCTATCACTGCAT TGTGTGTCTGATGGAATCCCTAAGCCAAATATCAAATGGACTATGCCAAGTGGTTATGTA GTAGACAGGCCTCAAATTAATGGGAAATACATATTGCATGACAATGGCACCTTAGTCATT AAAGAAGCAACAGCTTATGACAGAGGAACTATATCTGTAAGGCTCAAAATAGTGTGGT CATACACTGATTACTGTTCCAGTAATGATTGTAGCCTACCCTCCCGAATTACAAATCGT CCACCCAGGAGTATTGTACCAGGACAGGGGCAGCCTTTCAGCTCCACTGTGTGGCCTTG GGAGTTCCCAAGCCAGAAATCACATGGGAGATGCCTGACCACTCCCTTCTCTCAACGGCA AGTAAAGAGAGGACACATGGAAGTGAGCAGCTTCACTTACAAGGTACCCTAGTCATTGAG AATCCCCAAACCTCCGATTCTGGGATATACAAATGCACAGCAAAGAACCCACTTGGTAGT GATTATGCAGCAACGTATATTCAAGTAATCTGA </pre>
Restriction Sites:	Please inquire
ACCN:	NM_001178145
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001178145.1, NP_001171616.1</u>
RefSeq Size:	2824 bp
RefSeq ORF:	1953 bp
Locus ID:	285313
UniProt ID:	<u>Q6WRI0</u>
Cytogenetics:	3q25.1
Protein Families:	Druggable Genome, Transmembrane
Gene Summary:	Involved in the control of early migration of neurons expressing gonadotropin-releasing hormone (GNRH neurons) (By similarity). May be involved in the maintenance of osteochondroprogenitor cells pool (By similarity).[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) differs in the 5' UTR and coding sequence compared to variant 1. The resulting isoform (2) has a shorter and distinct N-terminus compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.