

Product datasheet for **SC329034**

LARP4 (NM_001170808) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LARP4 (NM_001170808) Human Untagged Clone
Tag:	Tag Free
Symbol:	LARP4
Synonyms:	PP13296
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC329034 representing NM_001170808.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGTTGCTTTTCGTGGAGCAGGTAGCATCTAAAGGAAGTGGTTTAAATCCTAATGCCAAAGTATGGCAA
GAAATGCTCCTGGAAATCTGATGCCACCCAGTAAGTCACTGGAAGTCTTGGCATGAAATA
GCAGCTACATCAGGTGCTCATCCTGAGGGTAATGCAGAGCTCTCAGAAGATATATGTAAAGAATATGAA
GTAATGTATTCTTCATCTTGTGAAACCACAAGAAATACTACAGGCATTGAAGAATCAACTGATGGGATG
ATTTTAGGACCAGAAGATCTGAGTTACCAATATATGATGTTTCCGGAGAAAGCAATTCAGCAGTTTCT
ACAGAAGACCTAAAAGAATGTCTGAAGAAACAATTAGAATTCTGTTTTCCAGAGAAAATTTGTCAAAG
GATCTTTACTTGATATCTCAAATGGATAGTGATCAGTTTCCCAATTTGGACAGTTGCCAACATGGAA
GAAATAAAAAAGTTGACTACAGACCCTGATCTAATTCTTGAAGTGTAAAGATCTTCTCCCATGGTACAA
GTTGATGAGAAGGTGAGAAAGTGAAGCAAGTCAATAAGCGTTGATTGTAATCTTAGAGAGATTCCT
GAAACAACACCAATAGAGGAAGTGAAGGTTTGTCAAAGTGAAGTCAAGTCAAGTCAAGTCAAGTCAAGT
GAGTTTGCACACAATAGCAACTGGTATATCACTTTCCAGTCAGACACAGATGCACAACAGGCTTTTAAA
TACTTAAGAGAAGAAGTTAAAACATTTACAGGGCAAGCCAATTATGGCAAGGATAAAAAGCCATCAATACA
TTTTTGTAGAATGGTTATCGATTAATGGATTCTAGTATCTATAGTCACCCCAATTCAAACTCAAGCA
CAGTATGCCTCCCGAGTCTTATGCAGCCTGTATATAATCCTCACCAACAGTACTCGGTCTATAGTATT
GTGCCTCAGTCTGGTCTCCAAATCCTACACCTTACTTTGAAACACCACTGGCTCCCTTTCCCAATGGT
AGTTTTGTGAATGGCTTAAATTCGCCAGGATCTTATAAAACAAATGCTGCTGCTATGAATATGGGTCGA
CCATTCAAAAAATCGGAGAACTCTCTCAGAGGTCGAAGACGACGAGAAGATGACAGGATCTCAAGA
CCTCATCTTCAACAGCTGAATCAAAGGCTCCAACACCAAGTTGACTTATTAGCCTCAATTTTCCA
CCTTTACCTGGAAGTTCATCAAGAATGCCAGGTGAAGTCAAGTCAAGTCAAGTCAAGTCAAGTCAAGT
AAAGGTGTCTACAAAGAAAAGGATAATGAAGAGTTGACAATTAGTTGCCAGTGCCTGCAGATGAGCAG
ACAGAATGCATTCTGCCAGCAACTCAATATGAGTACCAGTTCTCCATGTGCTGCTGAGCTTACTGCA
TTAAGCACAACCTCAGCAAGAAAAGGATCTAATAGAAGATTCCTCTGTTTCAAGGATGGTCTCAATCAG
ACAATAATACCAGTTTCTCCTCAAGTACTACAAAGCCATCGAGGGCAAGTACTGCTTACCATGTAAT
AATAACATAAATGCAGCTACAGCTGTGGTCTACAGGAACCCGAAAGTTAAGTTATGCTGAAGTGTGC
CAGAAGCCCCCTAAAGAGCCATCTCAGTCTTGTGACGCACTACGGGAACTTCGCTCCAATGTGGTG
TCTCCACCAAAAATGAAGACAATGGAGCTCTGAGAAGTCCGTTGAGAAACCACATGAGAAGCCAGAA
GCAAGGGCTAGTAAGGATTATTCTGGCTCCGAGGCAATAAATCCCGAGGGAGCAGCAGGAAAAATC
AGGGAACAGAGACGCCAGTTTAGCCATAGGGCTATACCTCAGGGAGTGACTCGACGTAATGGCAAAGAG
CAATATGTGCCACCCAGATCACCAAAGTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001170808

Insert Size: 1962 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001170808.1](#)

RefSeq Size: 6324 bp

RefSeq ORF: 1962 bp

Locus ID: 113251

UniProt ID: [Q71RC2](#)

Cytogenetics: 12q13.12

Protein Families: Transcription Factors

MW: 72.7 kDa

Gene Summary: RNA binding protein that binds to the poly-A tract of mRNA molecules (PubMed:21098120). Associates with the 40S ribosomal subunit and with polysomes (PubMed:21098120). Plays a role in the regulation of mRNA translation (PubMed:21098120). Plays a role in the regulation of cell morphology and cytoskeletal organization (PubMed:21834987, PubMed:27615744). [UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (6) lacks an alternate in-frame exon, compared to variant 1, resulting in a shorter protein (isoform f), compared to isoform a.