

## Product datasheet for **SC329051**

### **TTC39B (NM\_001168340) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	TTC39B (NM_001168340) Human Untagged Clone
Tag:	Tag Free
Symbol:	TTC39B
Synonyms:	C9orf52
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC329051 representing NM\_001168340.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGATGCTGTGCTGGCCTGCCGGCTCCGGGGCGGGGAAATCGAGTAGCTGCACTCAGGCCCGCCCGG
AGACCCGGAGGGAGCGCCGGCCGCTCCCTTTTGGCCTCCTTTGGCTGGGCTGAGCCAGAGCCGAGA
GCAGGAGTCGGCTCTGAGTTCCTGCTTGGTTTTTGGGTGGCAGCAGCCAGAGGAGGAATATGGCGCTC
TTAGGCAGCCGAGCGGAGCTGGAAGCGGACGAGGACGTTTTTCAAGATGCCTTGAAACCATCTCAATA
TCTTCTCACTCAGATATGGCAACAAGCAGCCTTCACTTTGCATCATGTGATACACAACAGGCGCCGAGA
CAGCGCGGAGCGTCTACTGTCAGCAGTTCATCATCAACCAAGGTGGATCTCAAGAGTGGCCTCGAAGAA
TGTGCTGTGGCATTGAATTAATTTCTAAGCAACAAATTTACAGACGCCTTAGAATTGCTTCGCCCTGG
GCTAAGGAGAGTATGTACCATGCCTGGGCTACAGTACCATTGTGGTGTGCAGGCTGCCTGACCTTC
GAGCAACAGGACATCCAGAACGGCATTCTGCCATGAAGGACGCTTTACAAACCTGCCAAAAATACAGG
AAAAAATACACAGTTGTAGAATCTTTCTCAAGTCTTCTTCTAGAGGATCCCTGGAGCAACTGAGTGAA
GAGGAAATGCATGCTGAAATCTGTTATGCCGAGTGTCTCTACAGAAAGCTGCACTCACGTTTGTGCAG
GATGAAAATATGATCAACTTCATCAAAGGTGGACTCAAAATTAGAACAAGTTACCAAATATATAAGAA
TGTCTTTCTATCCTGCATGAAATTCAGAAGAACAACCTCAGCAGGAATCTTCTATGAGTTGAAGGG
GGTGTCAAGCTTGGCAGTGGGCCCTTAATTTGATGCTGTCCCTTTTACCAGCAAGGATTATCAGACTA
CTAGAATTTATTGGATTTTCTGGGAATAGGGAGTTGGGCTCCTGCAGTTACGTGAAGGTGCTTCAGGA
AGGAGCATGAGGTCTGCACTGTGCTGCTTAACCATCTTAGCTTTTCACTTACATCTCCCTAATACTT
GGTACAGGAGAAGTGAATGTTGCGGAAGCAGAGCGTCTCCTGGCACCTTCTCCAGCAGTTTCCAAT
GTACGCTCAGTATTTCAAAGTGCACAAGAAGATTTCAAAAATGCATTTCAAGTTCAAGAAGAATGGAAA
CAGTTTACCATCTCTGCTACTGGGAGCTAATGTGGATTAATGTTTTCCAACAAAACCTGGATGCAGGCA
TATTACTATTCAGATCTGCTTTGCAAAGAGAGTAAATGGTCCAAGGCAACTTATGTGTTCTTGAAAGCA
GCAATTTTGAATGCTTCCAGAGGAGGATGTAGTAGCAACGAATGAGAATGTGGTAACTTTATTCAGA
CAGGTGGACAGCTTGAAGCAGAGAATTGCCGGGAAATCTATCCCTACTGAGAAGTTGCTGTGAGGAAG
GCTCGGCGTTACTCCGCCTCCTTACCTGCACCTGTGAAGCTCATCTTACCTGCCCTGGAATGATGTAT
GTCTGGAATGGTTTTTCAATAGTGAGCAAAAGAAAAGACCTTTCTGAAAATCTGTTAGTAACTGTTGAA
AAAGCTGAGGCAGCTTTACAAAGTCAAAATTTTAAACAGCTTCTCTGTGGATGATGAGTGCCTTAGTGAAG
TACTTAAAGGATGTTGCCTCAAGAACTTACAGCGGCCCTTGCAAGCTGAACTATGTTACAATCATGTT
GTGGAAAGTGAAAAGCTACTGAAGTATGACCACTACCTAGTGCCGTTCACTCTATTTGAATTGGCATCT
TTGTATAAAAGCCAGGGGAAATTTGACAAGGCCATAAAGTTCTAGAAAAGTCAAGGAACTACAAA
GATTACTCCCTGGAGTCCAGACTACACTTCAGAATTCAGGCAGCTCTTACCTCTGGAGAAAACCTTCC
TCAGATTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul

**Plasmid Map:**


**ACCN:** NM\_001168340

**Insert Size:** 2010 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001168340.1](#)

**RefSeq Size:** 3237 bp

**RefSeq ORF:** 2010 bp

**Locus ID:** 158219

**UniProt ID:** [Q5VTQ0](#)

**Cytogenetics:** 9p22.3

**MW:** 75.5 kDa

**Gene Summary:** Regulates high density lipoprotein (HDL) cholesterol metabolism by promoting the ubiquitination and degradation of the oxysterols receptors LXR (NR1H2 and NR1H3).  
[UniProtKB/Swiss-Prot Function]  
Transcript Variant: This variant (3) has multiple differences in the central coding region but maintains the reading frame, compared to variant 1. This results in a shorter protein (isoform 3), compared to isoform 1.