

Product datasheet for **SC329053**

Granuphilin (SYTL4) (NM_001174068) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Granuphilin (SYTL4) (NM_001174068) Human Untagged Clone
Tag:	Tag Free
Symbol:	Granuphilin
Synonyms:	SLP4
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF:	>NCBI ORF sequence for NM_001174068, the custom clone sequence may differ by one or more nucleotides <pre> ATGTCCGAGTTACTGGACCTTTCTTTCTGTCTGAGGAGGAAAAGGATTTGATTCTCAGT GTTCTACAGCGAGATGAAGAGGTCCGAAAGCAGATGAGAAAAGGATTAGGCGACTAAAG AATGAGTTACTGGAGATAAAAAGGAAAGGGCCAAAGAGGGGCAGCCAACACTACAGTGAT CGGACCTGTGCCCGGTGCCAGGAGAGCCTGGGCCGTTTGTAGTCCCAAAACCAATACTTGT CGGGGTGTGAATCACCTGGTGTGTGCGGACTGCCGCATACAGGAAAGCAATGGTACCTGG AGGTGCAAGGTGTGCGCCAAGGAAATAGAGTTGAAGAAAGCAACTGGGGACTGGTTTTAT GACCAGAAAGTGAATCGCTTTGCTTACCGCACAGGTAGTGAGATAATCAGGATGTCCTG CGCCACAAACCTGCAGTGAGTAAAAGAGAGACAGTGGGACAGTCCCTCCTTCATCAGACA CAGATGGGTGACATCTGGCCAGGAAGAAAGATCATTGAGAGCGGCAGAAAGGAGCCAGT GTGCTATTTGAAGTGCCAAAGCTGAAAAGTGAAAGAGTGCATTGGAAGCTGAGAGTGAG AGTCTGGATAGCTTCACAGCTGACTCGGATAGCACCTCCAGGAGAGACTCTCTGGATAAA TCTGGCCTCTTCCAGAATGGAAGAAGATGTCTGCTCCCAATCTCAAGTAGAAAAGGAA ACTCAGCCTGGAGGTCAAATGTGGTATTTGTGGATGAGGGTGAGATGATATTTAAGAAG AACACCAGAAAAATCCTCAGGCCTTCAGAGTAACTAAATCTGTGATAGATCTTCGCCCA GAAGATGTGGTACATGAAAGTGGCTCCTTGGGAGACAGAAGCAAAATCCGTCAGGCCTC AATGTGGATATGGAAGAGGAAGAAGAAGAAGACATTGACCACCTAGTGAAGTTACAT CGCCAGAAGCTAGCCAGAAGCAGCATGCAAAGTGGCTCCTCCATGAGTACGATCGGCAGC ATGATGAGCATCTACAGTGAAGCTGGTGATTTCCGGGAACATCTTTGTGACTGGCAGGATT GCCTTTTCCCTGAAGTATGAGCAGCAAAACCCAGAGTCTGGTTGTCCATGTGAAGGAGTGC CATCAGCTGGCCTATGCTGATGAAGCCAAGAAGCGCTCTAACCCTATGTGAAGACTTAC CTTCTGCCTGACAAGTCCCGCCAAGGAAAAAGAAAAACCAGCATCAAGCGGGACACTATT AATCCACTATATGATGAGACGCTGAGGTATGAGATCCCAAGTCTCTCCTGGCCAGAGG ACCCTGCAGTTCTCAGTTTGGCATCATGGTCGTTTTGGCAGAAACACTTTCTTGGAGAG GCAGAGATCCAGATGGATTCTGGAAGCTTGATAAGAACTGGATCATTGCCTCCCTTTA CATGGAAGATCAGTGCTGAGTCCCGACTGGCTTGCCATCACACAAAGGCGAGTTGGTG GTTTCATTGAAATACATCCAGCCTCCAAACCCCTGTTGGAGGTGACCGGAAAAAGAGT AAAGGTGGGGAAGGGGGAGAGCTCCAGGTGTGGATCAAAGAAGCCAAGAAGTACGGCT GCCAAAGCAGGAGGGACTTCAGACAGCTTTGTCAAGGGATACCTCCTCCATGAGGAAC AAGGCCAGTAAACGTAAACTCCTGTGATGAAGAAGACCCTGAATCCTCACTACAACCAT ACATTTGTCTACAATGGTGTGAGGCTGGAAGATCTACAGCATATGTGCCTGGAACAGT GTGTGGGACCGGGAGCCCTGGCCAGCAATGACTTCTGGGAGGGGTGAGGCTGGGTGTT GGCACTGGGATCAGTAATGGGGAAGTGGTGGACTGGATGGACTCGACTGGGGAAGAAGTG AGCCTGTGGCAGAAGATGCGACAGTACCCAGGGTCTTGGGCAGAAGGGACTCTGCAGCTC CGTTCCTCAATGGCCAAGCAGAAGCTGGGTTTATGA </pre>
Restriction Sites:	Please inquire
ACCN:	NM_001174068
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001174068.1](#), [NP_001167539.1](#)

RefSeq Size: 3767 bp

RefSeq ORF: 2016 bp

Locus ID: 94121

UniProt ID: [Q96C24](#)

Cytogenetics: Xq22.1

Gene Summary: This gene encodes a member of the synaptotagmin like protein family. Members of this family are characterized by an N-terminal Rab27 binding domain and C-terminal tandem C2 domains. The encoded protein binds specific small Rab GTPases and is involved in intracellular membrane trafficking. This protein binds Rab27 and may be involved in inhibiting dense core vesicle exocytosis. Alternate splicing results in multiple transcript variants that encode the same protein. [provided by RefSeq, Mar 2010]
 Transcript Variant: This variant (3) differs in the 5' UTR, compared to variant 1. Variants 1, 2 and 3 encode the same protein.