

Product datasheet for **SC329060**

PCCA (NM_001178004) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PCCA (NM_001178004) Human Untagged Clone
Tag:	Tag Free
Symbol:	PCCA
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC329060 representing NM_001178004.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGGGGTTCTGGGTCGGGACAGCACCGCTGGTCGCTGCCGGACGGCGTGGGCGGTGGCCGCCGAC
CAGCTGATGCTGAGCGCGCGCTGCCGACCCCTGAAGCATGTTCTGTACTATTCAAGACAGTGCTTAATG
GTGTCCCGTAATCTTGGTTCAGTGGGATATGATCCTAATGAAAAAATTTTGATAAAATCTTGTGCT
AATAGAGGAGAAAATTGCATGTCGGGTTATTAGAACTTGCAAGAAGATGGGCATTAAGACAGTTGCCATC
CACAGTGATGTTGATGCTAGTTCGTTTCATGTGAAAATGGCGGATGAGGCTGTCTGTGTTGGCCAGCT
CCCACCAGTAAAAGCTACCTCAACATGGATGCCATCATGGAAGCCATTAAGAAAACCAGGCCCAAGCT
GTACATCCAGGTTATGGATTCCTTTTCAGAAAACAAGAATTTGCCAGATGTTTGGCAGCAGAAGATGTC
GTTTTTCATTGGACCTGACACACATGCTATTCAAGCCATGGCGACAAGATTGAAAGCAAATATTAGCT
AAGAAAGCAGAGGTTAATAACAATCCCTGGCTTTGATGGAGTAGTCAAGGATGCAGAAGAAGCTGTCAGA
ATTGCAAGGGAAAATTGGCTACCTGTGATCAAGGCCTCAGCAGGTGGTGGTGGGAAAGGCATGCCG
ATTGCTTGGGATGATGAAGAGACCAGGGATGGTTTTAGATTGTCATCTCAAGAAGCTGCTCTAGTTTT
GGCGATGATAGACTACTAATAGAAAAATTTATTGATAATCCTCGTCATATAGAAATCCAGGTTCTAGGT
GATAAACATGGGAATGCTTTATGGCTTAATGAAAGAGAGTGCTCAATTCAGAGAAGAAATCAGAAGGTG
GTGGAGGAAGCACCAAGCATTTTTTTGGATGCCGAGACTCGAAGAGCGATGGGAGAAACAGCTGTAGCT
CTTGCCAGAGCAGTAAAAATTCCTCTGCTGGGACCGTGGAGTTCCTTGTGGACTCTAAGAAGAATTTT
TATTTCTTGAAATGAATACAAGACTCCAGTTGAGCATCCTGTACAGAATGCATTACTGGCCTGGAC
CTAGTCCAGGAAATGATCCGTGTTGCTAAGGGCTACCTCTCAGGCACAAAACAAGCTGATATTCGCATC
AACGGCTGGGCACTTGAATGTCGGGTTTATGCTGAGGACCCCTACAAGTCTTTTGGTTTACCATCTATT
GGGAGATTGTCTCAGTACCAAGAACCCTTACATCTACCTGGTGTCCGAGTGGACAGTGGCATCCAACCA
GGAAGTGATATTAGCATTATTATGATCCTATGATTTCAAAAATAATCACATATGGCTCTGATAGAAT
GAGGCACTGAAGAGAATGGCAGATGCACTGGATAACTATGTTATTCGAGGTGTTACACATAATATTGCA
TTACTTCGAGAGGTGATAATCAACTCACGCTTTGTAAAAGGAGACATCAGCACTAAATTTCTCTCCGAT
GTGTATCCTGATGGCTTCAAAGGACACATGCTAACCAAGAGTGAGAAGAACCAGTTATTGGCAATAGCA
TCATCATTGTTTGTGGCATTCCAGTTAAGAGCACAAACATTTTCAAGAAAATTCAGAATGCCTGTTATT
AAACCAGACATAGCCAACTGGGAGCTCTCAGTAAAATTGCATGATAAAGTTCATACCGTAGTAGCATCA
AACAAATGGGTCAAGTCTCGGTGGAAGTTGATGGGTGAAACTAAATGTGACCAGCACGTGGAACCTG
GCTTCGCCCTTATTGTCTGTGAGCGTTGATGGCACTCAGAGGACTGTCCAGTGTCTTTCTCGAGAAGCA
GGTGGAAACATGAGCATTGATTTCTTGGTACAGTGGTAGCAGAAGGTCAAGAAATTTGTGTGATTGAA
GCCATGAAAATGCAGAATAGTATGACAGCTGGGAAAACGGCACGGTGAATCTGTGCACTGTCAAGCT
GGAGACACAGTTGGAGAAGGGGATCTGCTCGTGGAGCTGGAAATGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:



ACCN: NM_001178004

Insert Size: 2046 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001178004.1](#)

RefSeq Size: 2436 bp

RefSeq ORF: 2046 bp

Locus ID: 5095

UniProt ID: [P05165](#)

Cytogenetics: 13q32.3

Protein Families:	Druggable Genome
Protein Pathways:	Metabolic pathways, Propanoate metabolism, Valine, leucine and isoleucine degradation
MW:	75 kDa
Gene Summary:	<p>The protein encoded by this gene is the alpha subunit of the heterodimeric mitochondrial enzyme Propionyl-CoA carboxylase. PCCA encodes the biotin-binding region of this enzyme. Mutations in either PCCA or PCCB (encoding the beta subunit) lead to an enzyme deficiency resulting in propionic acidemia. Multiple transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, May 2010]</p> <p>Transcript Variant: This variant (3) lacks an alternate in-frame exon in the 3' coding region compared to variant 1. The resulting protein (isoform c) is shorter compared to isoform a.</p> <p>Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>