

Product datasheet for **SC329063**

C18orf63 (NM_001174123) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C18orf63 (NM_001174123) Human Untagged Clone
Tag:	Tag Free
Symbol:	C18orf63
Synonyms:	DKFZP781G0119
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC329063 representing NM_001174123.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGAATGATTCTAGGCAGCAGTCTCTGTTCTTCATCACACTTCCAGATTTAAACAACTCTGTGCTGTC
AGAATAAATACTGAGTAATAAGGTGGCAGATACTGAGATTAGGACTATACAAATGAAGATGTGCAGGCAA
TTGCTGTTTTTGCATCAAGATATTCTTACATCACCTGTGTCTGGAATATTAACCAAATCTGGGTGGTG
ATGGCGATACCATTTTATAAAGCAAGAAAACTTAATGCCTATGTTGAAAAATATGGAGCTAAGATGGAG
GCTCCACAAAGAGTAATTCCTGTAATTCTTCAGAACTGCCTGTCATATTCATTTCATGGCTAGACTTGCT
CCAGCCTGGAATAGAAGTGGTCATCTCTTGATACAAGGAAGAGATTTCTTTCTCAGATGGGAAAAACAA
AGTGCTGTTGTATTAAACATCAATGTAACAGAACTCAAGTTTGTCTAAGTATAGAAGCTGCACAATC
AGACTGCCAGCACCTGAGCTAAAAGAGTTTGAGATTTCCAGAGTATTATAAAGGATTTTCATGCTAAC
AAGCATGCTGTCATTGAGAGACATTCCATTTAAGCAACTGGTGTACGTTTTGCCAAGTATGAAATG
GGCAAAATATAAATATTTTTCATGCCATCCCCGCTGCCTGTCTTTTCACTCATATGGAGATTTCCAG
AGACACTGGGATGCCCTGTATGGCTATAAACTTCAGGTGATTGTGAAAAAATAAAAATATACTGCAAC
ATCTATTTCAAAATGCTCGGAGAAAGGACCTTCACATATCCTCTCAGTTGCATCAGAAAGTCAGCCCATG
CAGTTCTTTCCAAGGTAGATTCCGAAGTTGTGTTGAAAAGTTTTCTTTCAGATTTAAAACTAAACTG
CCCCATATATGTGGATTTCCCATAAAGATGACAAGTAAACCATGCTACTACACACAAGAACTAACTAAA
CCTAATATACAGGAACACAAAGTCAAGCCACCCCAATTTGACCACTAAAAAATGCTTAGGGCATCTCTG
ACTCAAGCCACTTCCAGAAAGCCTGCCTGTGCTCAAGTCTTCTACCATGTTTCAGTAGCAGTGGACCAC
AAGGTGGAGCTTTCAGTCAGCCAGCCAACATCAGGCATTTTCTCAGCTTTCATCTCCAGCCAGAGTCT
GTTTCAGGTAGAAAGAAATCCCTGTCTATCAGGGCTCCACAAGTACATTGAGAAGTATTAATGCCAAC
AGAGGAAATACTCAAGTTCAGCACACAAATCTTAGCTCCCAAAGCAACATCACCCCTAAGTTTGTACCA
GTTTTCAAAATAGATTGTTACAAATGAACAAAAATACCTCAGTACTTGGCAGCCAAAAAGAAAAACAG
CATGATGTGACACAATCTAAATTGTTTTCACTCAAACTAGTATGATCCAGCATGACAACTGAATTTA
GGTCCAGCTATAAAAAACCGTTATAGTAGTAACATTGAGATGCAGGCTGTAACAATTTAAATCAGGAG
AATTCAGACCTCTGCAAGAAAAAATACAGAGTCTTCTGAAAAATGACAAAAATTTCCCTCTTCTCGT
GGAAAAATCGACTGTGAGTTTAAACAAAAATAAGCAACTATCTAATAGTGCAGTATTTGTGGTGTCAAT
AACAAATTTAGGGTGGTAAAAAGTGTGTTGACTTCCAAATGAAAGGAAAAAGAAATTTAACAGGCAAA
GGTATAACACAAATTTTAGGGAAAGCCATGGGTCACTAAAAGTAAAAGACAGCCACACATTTTGAA
TCAGATGGAGAAACCGAAGATCCAGCACTGCTACAGCAGCAATCAGAAAAATCAAGCTAAAGAAAGTTGGT
ACAAGTGACCACAGGTTGATAGTAAGCAAAATAGCCACAGGTCTAAAAGAAAAATATGTCCAGAGTCT
TCCAAAATTCAAAGAAGCATATTCCGATACTGTGCACTATGGCCAATCCAGTTCTTCTAAGAAGCAG
ATACTTGACTCGGATAAATCAAACTTAAGAAATCTCTCATCATTGCTTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001174123

Insert Size: 2058 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001174123.1</u>
RefSeq Size:	5065 bp
RefSeq ORF:	2058 bp
Locus ID:	644041
UniProt ID:	<u>Q68DL7</u>
Cytogenetics:	18q22.3
MW:	77.2 kDa