

Product datasheet for **SC329066**

EPB41L5 (NM_001184937) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	EPB41L5 (NM_001184937) Human Untagged Clone
Tag:	Tag Free
Symbol:	EPB41L5
Synonyms:	BE37; LULU; LULU1; YMO1; YRT
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC329066 representing NM_001184937.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCTGAGTTTCTCCGTAGAACAACACTAGGGCGTCGGTCTATGCGTAAACATGCAGAGAAGGAACGACTC
CGAGAAGCACAACGCGCCGCCACACATATTCTGCAGCTGGAGATTCTAAGTCCATCATCACGTGTCGG
GTGTCCCTTCTGGATGGTACTGATGTTAGTGTGGACTTGCCAAAAAAGCCAAAGGACAAGAGTTGTTT
GATCAGATTATGTACCACCTGGACCTGATTGAAAGCGACTATTTTGGTCTGAGATTTATGGATTCAGCA
CAAGTAGCACATTGGTTGGATGGTACAAAAAGCATCAAAAAGCAAGTAAAAATTGGTTCCACCCTATTGT
CTGCATCTTCGAGTTAAGTTTTATTCTCAGAACCAAATAACCTTCGTGAGGAGCTAACCCGGTATTTA
TTTGTCTTCAGTTAAAACAAGATATTCTCAGTGGAAAATTAGACTGTCCCTTTGATACAGCAGTGCAA
TTGGCAGCTTATAATCTGCAAGCTGAACCTGGTGACTATGATCTTGCTGAGCATAGTCTGAACTTGTC
TCAGAGTTCAGATTTCGTGCCTATTCAGACTGAAGAGATGGAAGTGGCTATTTTTGAGAAATGGAAGGAA
TACAGAGGTCAAACACCAGCACAGGCTGAAACCAATTATCTGAATAAAGCCAAATGGCTAGAAATGTAT
GGGTTGATATGCATGTGGTCAAGGCTAGAGATGGGAATGACTATAGTTGGGACTAACCCAACAGGA
GTCCTTGTTTTTGAAGGAGATACAAAATTGGCTTATTTTTTGGCCGAAGATAACCAGATTGGATTTT
AAGAAGATAAATTAACCTTGGTGGTGTAGAAGATGATGATCAGGGCAAAGAACAGGAACATACATTT
GTCTTTAGACTGGATCATCCAAAAGCATGCAAACATTTATGGAATGTGCTGTGGAGCATCATGCTTTC
TTCCGCTTCGAGGCCCGTCCAAAAGATTCTCATCGATCAGGATTTATTCGACTAGGATCACGATTT
AGATATAGTGGGAAAACAGAGTATCAGACCACAAAACCAATAAAGCAAGAAGATCAACATCCTTTGAA
AGAAGGCCAGCAAACGATATTCTAGACGAACTCTACAAATGAAAGCATGTGCTACAAAACCTGAAGAA
CTTAGTGTTCACAATAATGTTTCGACCCAAAGTAATGGCTCCCAACAGGCTTGGGGGATGAGATCTGCT
CTGCCTGTGAGTCTTCCATTTCTCTGCTCCTGTGCCAGTGGAGATAGAGAATCTTCCACAGAGTCTT
GGAACAGACCAGCATGACAGGAAATGCATTCTCTGAATATTGATTTGCTGAATAGCCCAGACTTATTG
GAAGCAACGATTGGTGTGTAATTGGGGCATCTGACACTATGGAACATCCCAAGCACTGAATGACGTT
AATGTAGCCACCAGGCTTCCGGGATTAGGGGAACCTGAAGTTGAATATGAGACATTAAGACACCTCA
GAGAAGCTCAAACAGCTTGAGATGGAGAACAGTCTTTGCTGTCCCCTCGATCCAACATCGATGTTAAC
ATAAACAGCCAGGAGGAAGTGGTGAAGTTGACTGAGAAATGCCTTAATAATGTCATTGAGAGCCAGGA
TTGAATGTCATGAGAGTCTCCTGACTTCAAGAGTAACATTTTGAAGGCTCAAGTAGAAGCAGTGCAT
AAGGTTACAAAAGAAGATAGCTTATTAAGTCATAAAAATGCCAATGTTTCAGGATGCTGCCACAAACAGT
GCTGTGTTAAATGAGAATAATGTGCCCTCCCAAGAGTCTCTTGAAGTCTGATGCTTATCACACCT
GCCGACAGTGGTCTGTTCTAAAGGAAGCTACAGATGAATTGGATGCCTTGTGCTGATCTCTAACTGAG
AATCTAATTGATCACACAGTTGCACCTCAGGTGCTTCCACATCCATGATCACACCCCGGTGGATTGTT
CCGCTCTGGTCCCATTTTGGCAGAAGAAGCTGTCTGAAGCAGAAGTGTCTACTGACCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:


ACCN: NM_001184937

Insert Size: 2064 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001184937.1](#)

RefSeq Size: 6490 bp

RefSeq ORF: 2064 bp

Locus ID: 57669

UniProt ID: [Q9HCM4](#)

Cytogenetics: 2q14.2

MW: 77.5 kDa

Gene Summary: May contribute to the correct positioning of tight junctions during the establishment of polarity in epithelial cells.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (2) lacks an exon in the 3' coding region compared to variant 1. The resulting protein (isoform 2) has a shorter, distinct C-terminus compared to isoform 1.
Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.