

Product datasheet for **SC329071**

HDX (NM_001177479) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	HDX
Synonyms:	CXorf43; D030011N01Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



Fully Sequenced ORF: >SC329071 representing NM_001177479.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAAATTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAATCTACGTTCTGTATTTACTGTAGAACAACAAAGGATTTACAGCGTTATTATGAAAATGGAATG
ACAAATCAAAGTAAAAATGCTTTAGCTCATATTACAGTGTGCACAGGAGACTAAGCTGGACTTCAGT
GTAGTCAGGACGTGGGTTGGCAATAAGAGAAGAAAGATGAGTAGTAAGAACTCTGAATCTGGAACAGCA
ACAACAGGAACCTCTTTGTAGCTCCAGACATCACAGTCAGAAATGTGGTTAATATTGCTCGACCTCA
AGCCAGCAGTCTTCTTGGACATCTGCCAATAATGATGTCTTGAAGTATATACAGTCCAGCCAGT
TCATCAAGTAGGCAAGGAACAACAACATACAGACACACAATTACAGAAGCACATAAAATCCCTATT
CAGAAAACAGCCACTAAAAATGATACTGAGTTTCAGTTACACATTCTGTCCAAAGACAAGTAGCACAC
TGTAATAATGCTTCCCTACTCTAGGTGAAAAACAATTATTTTGTCAAGACAGACAAGTGTCTAAAT
GCTGGAACTCAGTATTCAATCACGCAAGAAAACTATGGAACTCTTCACTACAAGCTTCTGAAATG
ACAGTACCTCAAAAGCCTTCTGTGTGCCACCGACCTTGTAATTTGAACAGTTGGGATTCAAAGGTCA
TATAAGCCTGAACACACAGGCCAGCATTACATAACTTATGTGGGCAAAAGCCAATATTAGAGACCT
TACTGTAGAACAACAACTTGGAAATCCGTGAAGTGTTTCATTGGCAGTTAGCGATTACCCCCAGAGA
ATTCTGGGAGGAAATGCCCCACAGAAGCCTAGCTCAGCAGAAGGAAATGTTTGTCCATTGCAATGGAG
ACTGGAGATGCTGAGGATGAATATGCCAGAGAGAGAGCTGGCATCGATGAGAGCAGACATACCAAGC
TATTCGAGATTTTATGAAAGTGGCAGTTCCCTTCGAGCTGAGAACCAGTACAACCTTGCCCGGACCA
GGAAAGAAATATGCCAAATTCACAAATGGTGAATATTAGAGATATGTCAGACAATGTACTGTATCAAAAC
AGAACTACCATTTGACACCACGGACCTCATTACATACAGCATCTAGTACAATGTACAGTAATACCAAT
CCATTACGGAGTAATTTTCTCCTCATTTGCATCATCAACCAATTGAGATTATCAAAAACCAAAAC
AATTACAGATTTTCAAGAAACCTTACTGTGCCTTGGATTACAGGGTGTTCTAGAAAAAGAGCACTACAG
GACCGCACTCAGTTCAGTGACCGAGACTTAGCCACCCTTAAGAAGTATTGGGACAATGGCATGACCAAGC
CTGGGCTCTGTTTGTAGAGAGAAAATGAAGCTGTGGCACTGAATTAATGTTGACTGTGAAATAGTT
CGGACTTGGATTGGGAATCGAAGAAGGAATATCGTTAATGGGGATTGAAGTTCCACCTCCAAGAGGA
GGCCCTGCTGATTCTCTGAGCAGCCTGAGTCTGGTTCTTTATCTGCACTCACACCAGGAGAGGAAGCT
GGGCTGAAGTAGGAGAGGATAATGACAGAAATGATGAAGTATCCATCTGTTTGTCTGAAGGAAGCTCT
CAAGAAGAGCCCAATGAAGTTGTTCCGAATGATGCAAGGGCTCATAAGGAAGAGGACCACCATGCAGTA
ACCACAGATAATGTGAAAATAGAAATATTGATGATGAAGAAAGTACATGATAAGTAATTCTGAAGTA
GAACAAGTAACTCTTTCTGGATTATAAGAATGAAGAAGTCAAATTCATTGAAAATGAGCTCGAGATT
CAAAAGCAAAATACTTTAACTTCAGACTTTTGTAGAAGCTTGATATTAGCAATGAAAGCTGATGAT
AAGGAACAACAGCAGGCACTGCTGTGAGATTTACCTCTGAATTAGAGGAAATGGATTTCAATCATGCC
TCACTGGAGCCTGATGATACCTCATTCACTGTATCTTCTTGTGAGAGAAAAATGTCTCAGAAAGTTG
TGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGCGCCGGC
```

Restriction Sites: SgfI-MluI

ACCN: NM_001177479

Insert Size: 2073 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_001177479.1
RefSeq Size:	6335 bp
RefSeq ORF:	2073 bp
Locus ID:	139324
UniProt ID:	Q7Z353
Cytogenetics:	Xq21.1
Protein Families:	Transcription Factors
MW:	77.2 kDa