

Product datasheet for **SC329103**

LDB3 (NM_001171610) Human Untagged Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | LDB3 (NM_001171610) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | LDB3 |
| Synonyms: | CMD1C; CMH24; CMPD3; CYPHER; LDB3Z1; LDB3Z4; LVNC3; MFM4; ORACLE; PDLIM6; ZASP |
| Mammalian Cell Selection: | None |
| Vector: | <u>pCMV6-XL5</u> |
| E. coli Selection: | Ampicillin (100 ug/mL) |



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001171610, the custom clone sequence may differ by one or more nucleotides

```

ATGTCTTACAGTGTGACCCTGACTGGGCCCGGGCCCTGGGGCTTCCGTCTGCAGGGGGC
AAGGACTTCAACATGCCCTCACTATCTCCCGGATCACACCAGGCAGCAAGGCAGCCCAG
TCCCAGCTCAGCCAGGGTGACCTCGTGGTGGCCATTGACGGCGTCAACACAGACCCATG
ACCCACCTGGAAGCCCAGAACAAGATCAAGTCTGCCAGCTACAATTGAGCCTCACCTG
CAGAAATCAAAGGTCCCATTCCATCTCCACGACAGCACCTCCAGTCCAGACCCCTCTG
CCGGTGATCCCTACCAGAAGGACCCCGCTCTGGACACGAACGGCAGCCTGGTGGCACCC
AGCCCCAGCCCTGAGGCGAGGGCCAGCCAGGCACCCAGGCACCCGGAGCTCAGGCC
ACCTTTAGCCCTGCCTTCTCCCGCCCTCCGCCTTCTCCTCACTCGCCGAGGCCTCTGAC
CCTGGCCCTCCGCGGGCCAGCCTGAGGGCCAAGACCAGCCAGAGGGGGCCCGGGACCTA
CTCGGCCAAAAGCCCTGCCGGCTCGAGCCAGCCGAGGCAATATAACAACCCCATTTGGC
CTGTACTCGGCAGAGACCCTGAGGGAGATGGCTCAGATGTACCAGATGAGCCTCCGAGGG
AAGGCCTCGGGTGTCCGACTCCCAGGAGCGCCGACTACCAGGAACGCTTCAACCCAGT
GCCCTGAAGGACTCGGCCCTGTCCACCCACAAGCCATCGAGGTGAAGGGGCTGGCGGC
AAGGCCACCATCATCCATGCGCAGTACAACACGCCCATCAGCATGTATCCAGGATGCC
ATCATGGATGCCATCGCTGGGCAGGCCAAGCCAAAGGCAGTGACTTCAGTGGGAGCCTC
CCTATTAAGGACCTTGCCGTAGACAGCGCCTCTCCCGTCTACCAGGCTGTGATTAAGAGC
CAGAACAAGCCAGAAGATGAGGCTGACGAGTGGGCACGCCGTTCTCCAACCTGCAGTCT
CGCTCCTTCCGCATCCTGGCCCAGATGACGGGGACAGAATTATGCAAGACCCTGATGAA
GAAGCTCTGCGAAGGTCAAGGCCCCAGGCCTTCTCCTACAGCCCCGAGTGGCCGCTCT
TCAGCACCTGCCACCCACACCAGCTACAGTGAAGGGCCCCGCCGCCCTGCACCCAAGCC
CGGGTTGTCAACACTGCCAGCATCCGGCCTTCTGTCTACCAGCCAGTGCCTGCATCTACC
TACAGCCCGTCCCAGGGGCAATTACAGTCCCCTCCCTACACCCCTCCCTGCCCT
GCCTACACCCCTCCCTGCCCTGCCTACACCCCTCACCTGTCCCCACCTACACTCCA
TCCCCAGCACCAGCCTATACCCCTCACCTGCCCCAACTATAACCCTGCACCCTCGGTG
GCCTACAGCGGGGCCCTGCGGAGCCTGCCAGCCGTCCACCCTGGGTGACAGATGATAGC
TTCTCCAGAAAGTTTCCCCGGGCAAGAGCACCACTCCATCAGCAAGCAGACCTGCC
CGGGGAGGCCAGCCTACACCCAGCGGTCTCAGGTGCCACCCTTCCAGGGGGACC
GTCCAGAGGGCTGAGCGATTCCAGCCAGCAGCCGACTCCACTCTCGGTCACTGCAAC
AATGTCATCCGGGCCCATTTCTGGTAGCCATGGGCCGTTCTTGGCACCTGAAGAGTTC
ACCTGTGCTACTGCAAGACTTCCCTGGCAGATGTGTGCTTTGTGAAGAGCAGAACAAC
GTTTACTGTGAGCGATGTTATGAGCAATCTTTGCCCCGCTGTGTGCCAAGTGCAACACC
AAAATTATGGGGGAAGTAATGCATGCCTTGAGACAGACATGGCACACCACCTGCTTCGTC
TGTGCGGCCTGCAAGAAGCCTTTTGGGAACAGCCTCTTCCACATGGAAGACGGGGAGCCC
TACTGCGAAGAAAGACTACATCAATCTGTTTCAGCACCAAGTGCCATGGCTGCGATTTCCC
GTGGAGGCTGGCACAAGTTTATCGAAGCCCTGGGCCACACTTGGCACGACACCTGCTTC
ATTTGCGCAGTCTGCCATGTGAATCTGGAGGGGACCCGTTCTACTCCAAGAAGGACAGA
CCCTGTGCAAGAAGCACGCACACCATCAACTTGATG

```

Restriction Sites: Please inquire

ACCN: NM_001171610

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

| | |
|-------------------------------|--|
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>NM_001171610.1, NP_001165081.1</u> |
| RefSeq Size: | 5436 bp |
| RefSeq ORF: | 2199 bp |
| Locus ID: | 11155 |
| UniProt ID: | <u>O75112</u> |
| Cytogenetics: | 10q23.2 |
| Gene Summary: | <p>This gene encodes a PDZ domain-containing protein. PDZ motifs are modular protein-protein interaction domains consisting of 80-120 amino acid residues. PDZ domain-containing proteins interact with each other in cytoskeletal assembly or with other proteins involved in targeting and clustering of membrane proteins. The protein encoded by this gene interacts with alpha-actinin-2 through its N-terminal PDZ domain and with protein kinase C via its C-terminal LIM domains. The LIM domain is a cysteine-rich motif defined by 50-60 amino acids containing two zinc-binding modules. This protein also interacts with all three members of the myozenin family. Mutations in this gene have been associated with myofibrillar myopathy and dilated cardiomyopathy. Alternatively spliced transcript variants encoding different isoforms have been identified; all isoforms have N-terminal PDZ domains while only longer isoforms (1, 2 and 5) have C-terminal LIM domains. [provided by RefSeq, Jan 2010]</p> <p>Transcript Variant: This variant (5) has multiple differences in the coding region, compared to variant 1. The encoded isoform (5) is longer than isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p> |