

Product datasheet for **SC329110**

NUTM2G (NM_001170741) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NUTM2G (NM_001170741) Human Untagged Clone
Tag:	Tag Free
Symbol:	NUTM2G
Synonyms:	FAM22G; NUTMG
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC329110 representing NM_001170741.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCTTCAAATGGAGCATACCCAGTGCTGGGACCCGGCGTGACCGTGAACCTGGCACCTCCCTGTCT
GTGTTACAGGCTCTGCCCTTTGCCACACCTCTCCCGGCCAACACACAGGCCGCCCTCGTGACTGCA
GTGGTTCTCCAGCCGGCCCTCTGGTGCTCTCTGCCTTCCCAGCACCCCTCTAGTGGCAGGACAGGAT
GGCCGCGGCCCAAGTGGGCTGGGCTTCCAACGCTTTGTCCAGATGAGGACAGAAGTGGGCTGTG
AAGCCCCCTCAGGCACAGACCTTGATCCTAACTCAGGCCCCCTCGTCTGGCAGGCTCCAGGCACCTC
TGTGGAGGTGTATGTGTCCACCTCCCCTACTCCTGGCAGCTGCTCTGGGGTGGCGTTACCTCTGCC
CAGGTGGTTGGGGCACCCAGGCCTGTGAGGGAGGCTGGTCCCATGGCCTTCTCTTCCACCACCACCA
CCGGCTGCCAGGTGGCCCCATCGTGTCCCAGGGAACGCTGGGCCATGGCCACAAGGGGCTCATGGA
GAGGGCAGCCTGGCTCCCTCCCAGGCCAAGGCCCGGCCGACGACTCCTGTAACCCAAGAGTGTCTAT
GAGAACTCCGACTCTGGCAGCACTACAAGCCCCTGGCCCGGAGGCACCTTCCCAGAGTCTGACACT
GAAGCGCTTTCCTGCTTCTCATCCAGTTCTCCGATCCCTGGCCCGCGGAAGCCACCATGACGCTG
GAGGAGGACTGTGGCGGCCATGCGGGAATGGCAGCACACGAGCAACTTGACCGGATGATTTTCTAC
GAGATGGCGGCAAGTTCTTGAGTTTGTAGGCTGAGGAGGAGATGCAGATTGAGAAATCGCAGTGGATG
AAGGGGCCCCAGAGCCTGCCTCTCCAGCCCCCGGAGGCTTGAACCTCGAGGACCCCTGCCCTGAG
GTGGTCAAGCAGCCAGTGATCTTCCAGCAAGGATGGCCCCAAGGGCCCGACTGCCTGCCACCA
CCCAGGCCCCAGAGGCCAGCGGAGACCAAGGCCACCTGCCACCACCCAGGCCCCGAGGCCAGCAGAG
ACCAAGGTCCCTGAGGAGATCCCCCTGAAGTGGTGCAGGAGTATGTGGACATCATGGAGGAGCTGCTG
GGGTCTCACCCTGGGGACACAGGGGAGCCTGAGGGACAACGGGAAAAGGGCAAAGTGGAGCAGCCGAC
GAAGAGGACGGGATGACCTCAGACCCGGGCTCCTGAGCTACATTGACAAGCTGTGTTCCAGGAAGAC
TTTGTCACCAAGGTGGAGGCCGTATTACCCCGATTCTGGAAGAATTGCTTTCCCAGATCCACAG
ATGGATTTCTTGCCCTAAGCCAGGAGCTGGAGCAGGAGGAAGGACTCACCTTGCCCAGCTAGTGGAG
AAGCGCTCCTGTCTTGAAGGAGAAAGGGTGGGGAGGGCAGCCCTCGACATGGCAGGCCAGGTTG
GACTCAAGTCTTCTGAGTTTGCAGCTGGCCAAGAAGCAGCGAGAGAGTCCCTGACCCCAACAAAGG
GTCAGCGTGGAACCTCCCCACCCAGACGGCTGCCAGGACCTCAGGGACAGGCAGAGTGCGCACT
GGCATGGCCAGGTCCGAAGACCTGCTGTGCTTTTGGGATGTGAGGACTCCCCAGGCTGAAGGCTGTC
CGGCCAACCTCTCTCCAGGACCACAGACCCACCTGCCCGGCTGGGGACCAAGGACGCTTGGGT
CTCCCTGGAGAGTCTCTGTCAAGGAGTACATGGGCTGGCTAAGGGGTCAAGTGAAGGAGACGGAATC
CCTGGCATGGTCTATGTCTGGGTTCCACACAGGCTGAGGCCCTGGAGGCTGTCCAGAGCCCTGTC
CCTTCTCGGGCCTTCTCAGCCAGGAGGGAGAGGACCCAGGGAGCTTCAATCTCCATCTGCTCAG
AAAAGAGGCTCAGCCATCACCTTCTCTGCCAGCAAGTCCAAGAAGCGACCTCTCTTTGGAAGCCCA
TCCCCTGCTGAAAGACACCGCACCCAGGGCTGGGCTCAGGGTCTCTGGGGAGCAATCCCTGGCTTG
GGGCTGGGTGGCCCTCACAGTCTAAAAGAGAAAGGGTGACCCCTTGGCTCCAGGAGGAAGAAGAAG
CGGCATTGTAGCCAGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001170741

Insert Size: 2226 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001170741.1</u>
RefSeq Size:	2516 bp
RefSeq ORF:	2226 bp
Locus ID:	441457
UniProt ID:	<u>Q5VZR2</u>
Cytogenetics:	9q22.33
Protein Families:	Secreted Protein
MW:	79 kDa