

Product datasheet for **SC329121**

TRPV4 (NM_001177433) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TRPV4 (NM_001177433) Human Untagged Clone
Tag:	Tag Free
Symbol:	TRPV4
Synonyms:	BCYM3; CMT2C; HMSN2C; OTRPC4; SMAL; SPSMA; SSQTL1; TRP12; VRL2; VROAC
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001177433, the custom clone sequence may differ by one or more nucleotides

```

ATGGCGGATTCCAGCGAAGGCCCGCGCGGGCCCGGGGAGGTGGCTGAGCTCCCCGGG
GATGAGAGTGGCACCCCAGGTGGGGAGGCTTTTCTCTCTCCTCCCTGGCCAATCTGTTT
GAGGGGGAGGATGGCTCCCTTTTCGCCCTCACCGGCTGATGCCAGTCGCCCTGCTGGCCCA
GGCGATGGGCGACCAAATCTGCGCATGAAGTTCAGGGCGCCTCCGCAAGGGGGTGCC
AACCCCATCGATCTGCTGGAGTCCACCCATATAGAGTCCCTCGGTGGTGCCTGGCCCAAG
AAAGCACCCATGGACTACTGTTTGACTACGGCACCTATCGTCACCACTCCAGTGACAAC
AAGAGGTGGAGGAAGAAGATCATAGAGAAGCAGCCGAGAGCCCCAAAGCCCTGCCCT
CAGCCGCCCCCATCCTCAAAGTCTTCAACCGGCTATCCTCTTTGACATCGTGTCCCGG
GGCTCCACTGCTGACCTGGACGGGCTGCTCCATTCTTGCTGACCCACAAGAAACGCCTA
ACTGATGAGGAGTTTCGAGAGCCATCTACGGGGAAGACCTGCCTGCCAAGGCCTTGCTG
AACCTGAGCAATGGCCGCAACGACACCATCCCTGTGCTGCTGGACATCGCGGAGCGCACC
GGCAACATGAGGGAGTTCATTAACCTGCCCTTCCGTGACATCTACTATCGAGGGGAGCTG
CCCCTGTGCTGGCTGCCTGCACCAACCAGCCCCACATTGTCAACTACCTGACGGAGAAC
CCCCACAAGAAGGCGGACATGCGGCGCCAGGACTCGCGAGGCAACACAGTCTGCATGCG
CTGGTGGCCATTGCTGACAACACCCGTGAGAACACCAAGTTTGTACCAAGATGTACGAC
CTGCTGCTGCTCAAGTGTGCCCGCTCTTCCCCGACAGCAACCTGGAGGCCGTGCTCAAC
AACGACGGCCTCTCGCCCTCATGATGGCTGCCAAGACGGGCAAGATTGGGAACGCCAC
GAGATGCTGGCTGAGAGCCATCAATGAACTGCTGCGGGAAGTGGCGCAAGTTCCGGG
GCCGTCTCCTTACATCAACGTGGTCTCTACCTGTGTGCCATGGTCTCTTACTCTC
ACCGCCTACTACCAGCCGCTGGAGGGCACACCGCCGTACCCTTACCGCACACGGTGGAC
TACCTCGGGTGGCTGGCGAGGTCAATTACGCTTCACTGGGGTCTGTTCTTCTTCCACC
AACATCAAAGACTTGTTCATGAAGAAATGCCCTGGAGTGAATTCTCTTCTTATTGATGGC
TCTTCCAGTGTCTACTTCTACTCTGCTGCTGGTATCGTCTCAGCAGCCCTTAC
CTGGCAGGGATCGAGGCCTACCTGGCCGTGATGGTCTTTGCCCTGGTCTGGGCTGGATG
AATGCCCTTACTTACCCCGTGGGCTGAAGCTGACGGGGACCTATAGCATCATGATCCAG
AAGATTCTTCAAGGACCTTTCCGATTCTGCTGCTACTTGTCTTCTCATGATCGGC
TACGCTTACGCCCTGGTCTCCCTCCTGAACCCGTGTGCCAACATGAAGGTGTGAATGAG
GACCAGACCAACTGCACAGTGCCCACTTACCCCTCGTGCCGTGACAGCGAGACCTTACG
ACCTTCTCCTGGACCTGTTAAGCTGACCATCGGCATGGGCGACCTGGAGATGCTGAGC
AGCACCAGTACCCCGTGGTCTTATCATCCTGCTGGTACCTACATCATCCTCACCTTT
GTGCTGCTCCTCAACATGCTCATTGCCCTCATGGGCGAGACAGTGGGCCAGGTCTCCAAG
GAGAGCAAGCACATCTGGAAGCTGACAGTGGGCCACCACCATCCTGGACATTGAGCGCTCC
TTCCCCGATTCTGAGGAAGGCCTTCCGCTCTGGGGAGATGGTACCCTGGGCAAGAGC
TCGGACGGCACTCCTGACCGCAGGTGGTCTCAGGGTGGATGAGGTGAACTGGTCTCAC
TGGAACCAGAACTTGGGCATCATCAACGAGGACCCGGGCAAGAATGAGACCTACCAAGTAT
TATGGCTTCTCGCATACCGTGGGCCCTCCGCAGGGATCGTGGTCTCCTCGTGGTACCC
CGCGTGGTGGAACTGAACAAGAACTCGAACCCGGACGAGGTGGTGGTGCCTCTGGACAGC
ATGGGGAACCCCGCTGCGATGGCCACCAGCAGGTTACCCCGCAAGTGGAGGACTGAT
GACGCCCGCTCTAG

```

Restriction Sites: Please inquire

ACCN: NM_001177433

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001177433.1</u> , <u>NP_001170904.1</u>
RefSeq Size:	2871 bp
RefSeq ORF:	2295 bp
Locus ID:	59341
UniProt ID:	<u>Q9HBA0</u>
Cytogenetics:	12q24.11
Protein Families:	Druggable Genome, Ion Channels: Transient receptor potential, Transmembrane
Gene Summary:	<p>This gene encodes a member of the OSM9-like transient receptor potential channel (OTRPC) subfamily in the transient receptor potential (TRP) superfamily of ion channels. The encoded protein is a Ca²⁺-permeable, nonselective cation channel that is thought to be involved in the regulation of systemic osmotic pressure. Mutations in this gene are the cause of spondylometaphyseal and metatropic dysplasia and hereditary motor and sensory neuropathy type IIC. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Apr 2010]</p> <p>Transcript Variant: This variant (3) lacks two in-frame exons in the coding region, compared to variant 1. The encoded isoform (c) is shorter, compared to isoform a.</p>