

## Product datasheet for **SC329129**

### ZDHC8 (NM\_001185024) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ZDHC8 (NM_001185024) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZDHC8
Synonyms:	DHC8; ZDHC1; ZNF378
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

**Fully Sequenced ORF:** >SC329129 representing NM\_001185024.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGCCCCGACGCCCGGGACGCGCCTCAAACCCGCCAAGTACATCCCGGTGGCCACGCCCCCGCGCTG
TTCCCGCTCACAATGGCATCATCTTCTCTTTGCTGGCCAACTTCAGCATGGCCACTTTCATGGAC
CCTGGTGTTCCTCCCGAGCGGATGAGGATGAGGACAAGGAGGACGACTTCCGGCTCCGCTGTACAAG
AACGTGGATGTGCGAGGTATCCAGGTCCGCATGAAGTGGTGTGCCACGTGCCACTTCTACCGCCCGCG
CGCTGCTCCCACTGCAGCGTCTGTGACAACTGTGTAGAGGACTTTGACCACCACTGCCCTGGGTCAAC
AACTGCATCGGGCGTCAAACATCGCTACTTCTCTGTTCTGCTGTCACTCAGTGCACACATGGTG
GGCGTCGTGGCCTTCGGCCTGGTCTACGTGCTGAACCACGCTGAGGGGCTGGGAGCCGCGCACACCACC
ATCACCATGGCTGTCATGTGTGGCCGGCCTCTTCTTCCATCCCTGTATTGGCCTCACTGGCTTCCAT
GTGGTGTGGTCACTCGGGGGCGCACCAACGAGCAGGTGACTGGGAAGTTCCCGGGGGTGTGAAC
CCTTTCACCCGAGGCTGCTGTGGGAATGTGGAGCACGTGCTGTGTAGCCCTGGCGCCCGGTAGGTG
GTGGAGCCACCCCGCTGCCGCTCGCGGTGAGTTTGAAGCCGCTTTCCTTAGGCCCTGAACTCCTGGAC
CGAGTGCACCGCTCAAGGTCAAGCTTAGTGACAACGGGCTGAAGGCTGGCCTGGGCCGTAGCAAGTCC
AAGGGCAGCCTGGACCGCTGGATGAGAAGCCACTGGACTTGGGGCCACCACTGCCCCCAAGATAGAG
GCTGGCACGTTACAGCAGTGACCTGCAGACCCCGCGCCAGGCAGTGTGAGAGTGCCCTGTCCGTGCAG
AGGACCAGCCCCGACACCTGCCATGTACAAGTTTAGCCGGCTTTCACCGGTCCCAAGGTGCC
TTCTGTGACAGCGCAGGTTCCAGGCCCTGATTCCCTGACCCTGGGGGACGACAGCATCCGTAGC
CTGGACTTGTGTCCGAGCCGAGCCTGGACCTCCCTGACTATGGGCCAGGGGGCTGCATGCAGCATTAC
CCGCCATCCCCACCGCTCAGCCCTCTGATGCCTTCTCGGGCGCTTTCGCTCCCTGAGCCTCAAGGCC
TCGAGCCGGCGGGGGGATCATGTGGCCTGCAGCCCTGCGCTCTGAGGGGGGGCCCCACGCC
CACCGTAGCATTTTTCGCCCATGCACTGCCAACCGCAACGGCAGCCTGTCTATGACAGCCTGCTC
AATCCTGGCTCGCCTGGTGGCCACGCTGCCCTGCCACCCAGCAGTTGGCGTGGCCGATACCACTCA
CCCTACCTGCATCCTGGGGCAACGGGGACCCGCCACGGCCCTACCCCGCAGCTTCAGCCCGTGTG
GGCCCCGCCCCGGGAGCCCTCGCTGTGCGCTACGACAACCTGTCCAGGACCATATGGCATCCATC
CAGGAGCGCAAGGACAGGAGGAGCGTGAGCGCCTGCTGCGCTCCAGGCCGACTCACTTTCGGCGAC
TCAGGCGTCTATGACGCTCCAGCTCCTACAGCCTGCAGCAGGCCAGTGTGCTGTCCGAGGGCCCCGA
GGTCCCGCTGCGCTATGGCTCCAGAGACGACCTTGTGGCTGGGCCCGCTTCGGTGGCCCGCAAC
CCTGCCCTGCAGACGTCCTGCTCGCTGTCCAGCTCCGTGAGCCGTGACCGGGACGTGCTCCTCC
TCCCTGCAGGCTGATCAGGCCAGCAGCAACGCCCGGGGGCCCCGGCCAGCAGTGGCTCACACAGGTCA
CCTGCACGCCAGGGCCTGCCCTCCCCGCCGCACTCCCACTACCATCCTACGCGGGCCCCAAAGCT
GTGCTTTCATCCACACGGACCTCCAGAGCCACCGCCTCGCTGACCGTGCAGAGGGGGCGGATTGGC
ACCTGCACCCGTGGATGGGGGCGCGTGGCCAGCCTTGGGTGCCTCCTGGGCTGCACCTGTGCCACCTT
GGCCGCCGAGGAGCCGCCACCACTGCGGGCCCCCTGGAGCCAGGCCCGGGGACCCCCACGCCGGG
GCCATGTGCCGCTGCATTGGCTCCAGTCTTTTCCCGAGCCTCTCGGGGCCTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_001185024

**Insert Size:** 2337 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001185024.1](#)

**RefSeq Size:** 3237 bp

**RefSeq ORF:** 2337 bp

**Locus ID:** 29801

**UniProt ID:** [Q9ULC8](#)

**Cytogenetics:** 22q11.21

<b>Protein Families:</b>	Transmembrane
<b>MW:</b>	83.2 kDa
<b>Gene Summary:</b>	<p>This gene encodes a four transmembrane protein that is a member of the zinc finger DHHC domain-containing protein family. The encoded protein may function as a palmitoyltransferase. Defects in this gene may be associated with a susceptibility to schizophrenia. Alternate splicing of this gene results in multiple transcript variants. A pseudogene of this gene is found on chromosome 22.[provided by RefSeq, May 2010]</p> <p>Transcript Variant: This variant (1) encodes the longer isoform (1).</p>