

Product datasheet for **SC329140**

TRP 7 (TRPC7) (NM_001167577) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TRP 7 (TRPC7) (NM_001167577) Human Untagged Clone
Tag:	Tag Free
Symbol:	TRPC7
Synonyms:	TRP7
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC329140 representing NM_001167577.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTTGAGGAACAGCACCTTCAAAAACATGCAGCGCCGGCACACAACGCTGAGGGAGAAGGGCCGTCGC
CAGGCCATCCGGGTCCCGCCTACATGTTCAACGAGAAGGGCACCAGTCTGACGCCCGAGGAGGAGCGC
TTCTGGACTCGGCTGAGTATGGCAACATCCCGTGGTCCGGAATAATGCTGGAGGAGTCCAAGACCTT
AACTTCAACTGTGTGGACTACATGGGGCAGAACGCTCTGCAGCTGGCCGTGGGCAACGAGCACCTAGAG
GTCACGGAGCTGCTGCTGAAGAAGGAGAACCTGGCACGGGTGGGGGACGCGCTGCTGCTGGCCATCAGC
AAGGGCTATGTGCGCATCGTGGAGGCCATCCTCAACCACCGGCCTTCGCGCAGGGCCAGCGCTGACG
CTCAGCCCCTGGAACAGGAGCTGCGGACGACGACTTCTATGCCTACGACGAGGACGGCAGCGCTTC
TCCCACGACATCAGCCCATCATCTGGCGGCGCACTGCCAGGAGTATGAGATCGTGCACATCTGCTG
CTCAAGGGCGCCCGCATCGAGCGGCCACGACTACTTCTGCAAGTGAATGAGTGCACCGAAGAACAG
CGGAAAGACTCCTTCAGCCACTCGGCTCGCGCATGAACGCCTACAAAGGACTGGCGAGTCTGCCTAC
TTGTCCCTGTCCAGCAAGACCTGTCCCTACCGCCCTGGAGCTCAGCAACGAGTTAGCCAGACTAGCC
AACATTGAGACTGAATTAAGTTCGTTGCTCATCTAACTGTGACGAGCAATTGCTTACCATGTGGTAT
GAAAATCTCTCAGGCTTACGTCAACAGTCTATCGCTGTGAAATTCCTGGCTGCTTTGGAGTCTCCATA
GGCCTCCCTTTCTCGCCATAGCCTATTGGATTGCTCCGTGCAGCAAGCTAGGACGAACCTGAGGAGC
CCTTTCATGAAGTTTGTAGCTCATGCAGTTTCTTTTACAATCTTCTTGGGATTATTAGTTGTGAATGCA
TCTGACCGATTTGAAGGTGTTAAAACCTGCCAAACGAAACCTTACAGACTACCCAAAACAAATCTTC
AGAGTAAAACACACAGTCTCCTGGACAGAAATGCTCATTATGAAGTGGTCTTAGGAATGATTTGG
TCCGAATGCAAGGAAATCTGGGAGGAGGGCCACGGGAGTACGTGCTGCACCTTGTGAACCTGCTGAT
TTCGGGATGCTGCCATCTTCGTGGCCTCTTCCAGCAGCGCTTATGGCCTTCTGAAGGCCACGGAG
GCACAGCTGTACGTGGACCAGCAGTGCAGGACGACACGCTGCACAATGTCTCGTTCCGCCGGAAGTG
GCATACTCACCTACGCCAGGACAAGTGGTGGCCTTCAGACCCTCAGATCATATCGGAAGGGCTCTAC
GCGATAGCCGTCGTGCTGAGCTTCTCTCGCATTGCATACATTCTGCCAGCCAACGAGAGTTTTGGGCC
CTGCAGATCTCGCTAGGGAGAACTGTGAAAGATATCTTCAAGTTCATGGTCATTTTCATCATGGTATTT
GTGGCCTTCATGATTGGGATGTTCAACCTGTACTCTTACTACCGAGGTGCCAAATACAACCCAGCGTTT
ACAACGGTTGAAGAAAGTTTTAAACTTTGTTTTGGTCCATATTCGGCTTATCTGAAGTAATCTCAGTG
GTGCTGAAATACGACCACAAATTCATCGAGAACATTGGCTACGTTCTCTACGGCCTTATAACGTCACC
ATGGTGGTAGTGTGCTCAACATGCTAATAGCCATGATAAACAACCTCTATCAGGAAATTGAGGAGGAT
GCAGATGTGGAATGGAAGTTCGCCCGAGCAAACTCTGGCTGTCTTACTTTGATGAAGGAAGAACTCTA
CCTGCTCCTTTAATCTAGTGCCAAGTCTAAATCATTTTATTATCTCATAATGAGAATCAAGATGTGC
CTCATAAACTCTGCAAACTAAGGCCAAAAGCTGTGAAAATGACCTTGAATGGGCATGCTGAATTCC
AAATTCAGAAGACTCGTACCAGGCTGGCATGAGGAATTCTGAAAATCTGACAGCAAATAACACTTTG
AGCAAGCCACCAGATACCAGAAAATCATGAAACGGCTCATAAAAAGATACGTCCTGAAAGCCAGGTG
GACAGAGAAAATGACGAAGTCAATGAAGGCGAGCTGAAGGAAATCAAGCAAGATATCTCCAGCCTGCGC
TATGAGCTTCTTGAGGAAAAATCTCAAGCTACTGGTGGAGCTGGCAGACCTGATTCAACAACCTCAGCGAG
AAGTTTGAAGAACTTAAACAAAGACCCTGAGGGTGAACAAGGGCAAAGACATTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001167577

Insert Size: 2406 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001167577.1](#)

RefSeq Size: 2752 bp

RefSeq ORF: 2406 bp

Locus ID: 57113

UniProt ID: [Q9HCX4](#)

Cytogenetics: 5q31.1

Protein Families:	Druggable Genome, Ion Channels: Transient receptor potential, Transmembrane
MW:	92.4 kDa
Gene Summary:	<p>Thought to form a receptor-activated non-selective calcium permeant cation channel. Probably is operated by a phosphatidylinositol second messenger system activated by receptor tyrosine kinases or G-protein coupled receptors. Activated by diacylglycerol (DAG) (By similarity). May also be activated by intracellular calcium store depletion.</p> <p>[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (2) lacks an alternate in-frame exon in the central coding region, compared to variant 1, resulting in an isoform (2) that is shorter than isoform 1.</p>