

Product datasheet for **SC329143**

SCUBE2 (NM_001170690) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SCUBE2 (NM_001170690) Human Untagged Clone
Tag:	Tag Free
Symbol:	SCUBE2
Synonyms:	CEGB1; CEGF1; CEGP1; scube/You
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC329143 representing NM_001170690.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGGGGTGCGGGGCCGAACCGTCCCGGGCGGCCTGGGCGGTGCTGCTGCTGCTGCTGCTGCCG
CCACTGCTGCTGCTGCGGGGGCCGTCCCGCGGGTGGGGCCGTGCCCGGGGCCGAGGAGATGTA
GATGAGTGTGCCAAGGGCTAGATGACTGCCATGCCGACGCCCTGTGTCAGAACACACCCACCTCTAC
AAGTGCTCCTGCAAGCCTGGCTACCAAGGGGAAGGCAGGCAGTGTGAGGACATCGATGAATGTGAAAT
GAGCTCAATGGAGGCTGTGTCCATGACTGTTTGAATATTCAGGCAATTATCGTTGCACTTGTTTTGAT
GGCTTCATGTTGGCTCATGACGGTCATAATTGTCTTGATGTGGACGAGTGCCTGGAGAACAATGGCGGC
TGCCAGCATACCTGTGTCAACGTCATGGGAGCTATGAGTGTGCTGCAAGGAGGGGTTTTCTGAGT
GACAATCAGCACACCTGCATTACCCTCGGAAGAGGGCCTGAGCTGCATGAATAAGGATCACGGCTGT
AGTCACATCTGCAAGGAGGCCCAAGGGGCAGCGTCGCCTGTGAGTGCAGGCCTGGTTTTGAGCTGGCC
AAGAACCAGAGAGACTGCATCTTGACCTGTAACCATGGGAACGGTGGGTGCCAGCACTCCTGTGACGAT
ACAGCCGATGGCCAGAGTGCAGCTGCCATCCACAGTACAAGATGCACACAGATGGGAGGAGCTGCCTT
GAGCGAGAGGACTGTCTGGAGGTGACAGAGAGCAACACCACATCAGTGGTGGATGGGGATAAACGG
GTGAAACGGCGGCTGCTCATGGAACGTGTGCTGTCAACAATGGAGGCTGTGACCGCACCTGTAAGGAT
ACTTCGACAGGTGCCACTGCAGTTGCTCTGTTGGATTCACTCTCCAGTTGGATGGGAAGACATGTAAG
GATATTGATGAGTGCCAGACCCGCAATGGAGGTTGTGATCATTCTGCAAAAACATCGTGGGCAGTTTT
GACTGCGGCTGCAAGAAAGGATTTAAATTAACAGATGAGAAGTCTTGCCAAGATGTGGATGAGTGC
TCTTTGGATAGGACCTGTGACCACAGCTGCATCAACCACCTGGCACATTTGCTGTGCTGCAACCGA
GGGTACACCCTGTATGGCTTACCACCTGTGGAGACCAATGAGTGCAGCATCAACAACGGAGGCTGT
CAGCAGGTCTGTGTGAACACAGTGGGCAGCTATGAATGCCAGTGCCACCCTGGGTACAAGCTCCACTGG
AATAAAAAAGACTGTGTGGCTTCTGTGACCTGAGCTGCATCGTAAAGCGAACCGAGAAGCGGCTCCGT
AAAGCCATCCGCACGCTCAGAAAGGCCCTCCACAGGGAGCAGTTTCACTCCAGCTCTCAGGCATGAAC
CTCGACGTGGCTAAAAAGCCTCCAGAACATCTGAACGCCAGGCAGAGTCTGTGGAGTGGGCCAGGGT
CATGCAGAAAACCAATGTGTGAGTTGCAGGGCTGGGACCTATTATGATGGAGCACGAGAACGCTGCATT
TTATGTCCAAATGGAACCTTCCAAAATGAGGAAGGACAAATGACTTGTGAACCATGCCCAAGACCAGGA
AATTCTGGGGCCCTGAAGACCCAGAAGCTTGAATATGTCTGAATGTGGAGTCTGTGTCAACCTGGT
GAATATTCTGAGATGGCTTTGCACCTTCCAGCTCTGTGCCCTGGGCAGTTCCAGCCTGAAGCTGGT
CGAACTTCTGCTTCCCCTGTGGAGGAGCCTTGCCACCAACATCAGGGAGCTACTTCTTTCCAGGAC
TGTGAAACCAGAGTTCAATGTTACCTGGACATTTCTACAACACCACCACTCACCGATGTATTCTGTGC
CCAGTGGGAACATACCAGCCTGAATTTGGAAAAATAATTGTGTTTCTTGCCAGGAAATACTACGACT
GACTTTGATGGCTCCACAACATAAACCAGTGTAAAACCTCATCCAATTCTGTGACAACATATGAAACC
TGCCAGACCTACGAACGCCCATCGCCTTACCTCCAGGTCAAAGAAGCTGTGGATTCAAGTCC
AATGAAGGAACAGCGCTAGAGGGTCCAGGTCCCATACGTGACATATGATGAGGACTACCAGGAACCTC
ATTGAAGCATAGTTGAGATGGCAGGCTCTATGCATCTGAGAACCATCAGGAAATACTTAAGGATAAG
AACTTATCAAGGCTCTGTTTGTGCTTGGCCATCCCAGAACTATTTCAAGTACACAGCCAGGAG
TCCCAGAGATGTTTCCAAGATCGTTCATCCGATTGCTACGTTCCAAGTGTCCAGGTTTTGAGACCT
TACAAATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001170690

Insert Size: 2424 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001170690.1](#)

RefSeq Size: 4058 bp

RefSeq ORF: 2424 bp

Locus ID: 57758

UniProt ID: [Q9NQ36](#)

Cytogenetics: 11p15.4

Protein Families: Transmembrane

MW: 88.9 kDa

Gene Summary: Lipid-binding protein required for SHH long-range signaling by binding to the dually lipid-modified SHH (ShhNp) and by promoting ShhNp mobilization, solubilization and release from the cell membrane (PubMed:22902404, PubMed:22677548). Acts by enhancing the proteolytic processing (shedding) of the lipid-modified N- and C- terminal of ShhNp at the cell surface (PubMed:24522195). Synergizes with DISP1 to increase SHH secretion (PubMed:22902404). Probable cell surface coreceptor for VEGFR2 involved in VEGFR2-mediated angiogenesis (PubMed:27834687).[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (2) has multiple differences in the coding region but maintains the reading frame, compared to variant 1. This variant encodes isoform 2, which is shorter than isoform 1.