

Product datasheet for **SC329156**

CAGE1 (NM_001170693) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CAGE1 (NM_001170693) Human Untagged Clone
Tag:	Tag Free
Symbol:	CAGE1
Synonyms:	bA69L16.7; CT3; CT95; CTAG3
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001170693, the custom clone sequence may differ by one or more nucleotides

```

ATGAACAAGGACTATCAAAAATTTTGGTCATCACCTTCAGATCCTGTACATTTTGAAGTG
GATACTTCTCATGAAAAAGTAGAAAGCATGTCAGAATCGGATACCATGAATGTCAGCAAT
CTTTCTCAAGGTGAATGCTTTTCACATTCTCCAATCTGTATGGAAACACCGGCACCACT
TGTGACTTGCCTCAGAACGAAATAAAGAATTTTGAAGGGAAAAATGAGTATGAATCCACA
CTTTGTGAAGATGCTTATGGCACACTAGACAATTTGTTAAATGATAACAACATTGAAAAAT
TACTCAACGAATGCACTAATTCAGCCAGTTGACACCATCAGCATATCTTCCTTGAGACAA
TTTGAAACCGTTTGCAAATTTCACTGGGTAGAAGCATTGATGATGAAATGACAGAGAAG
CCAGAATTTCAAAGTCAAGTGTATAATTATGCAAAAGACAACAATATAAAGCAAGACTCA
TTTAAGGAAGAAAAATCCAATGGAACTAGTGTTTCTGCAATACAGACCAACTTGGTAAC
GAGTATTTTAGACAGCCTCCTCTAGAAGCCCGCCTCTAATCCACTGCAGTGAGAGATG
CTGAAATTTACAGAAAAGTCACTGGCTAAAAGTATTGCCAAGGAATCTGCCCTCAACCCT
AGCCAACCTCCAAGCTTCTTATGTAAGACTGCAGTTCTTCAAAGAAAATACAGAATTAT
GGGGAGATTCTGAGATGTCAGTCAGTTATGAAAAGGAAGTCACAGCAGAGGGGTGTGGAG
AGGCCAGAAAATTGCTCAACTTGGTCTTCGGCAGGCATTTCTTGAGGAGTGAAGCATGT
CGGGAGAACTGTGAGATGCCTGACTGGGAGCAAAGTGCTGAAAGCTTACAACCTGTTCAA
GAGGACATGGCTTTAAATGAAGTCTTGCAAGAAATTAACATACTAACAGAAAGCAGGAA
GTGCGAATCCAAGAACTCCAGTGTAGTAACCTGTATTTAGAGAAGAGGGTTAAAGAACTA
CAGATGAAGATTACCAACAGCAAGTGTTTCATTGATGTCATCAATAAGCTAAAGGAGAAT
GTTGAAGAATTAATTGAAGACAAATACAAAATAATCCTAGAGAAGAATGATACTAAAAAG
ACATTGCAGAATTTGGAAGAGGTTTGTAGCTAACACGCAAAAACATCTTCAGGAATCCAGG
AATGACAAGGAAATGTTACAGCTTCAATTTAAGAAGATCAAGGCTAATTATGTGTGTTTA
CAGGAAAAGGTACATGACTGAAATGCAACAAAAAATAAATCTGTAAAGTCAGTATTAGAG
ATGGACAAAACCTTAAGCAAGAAAGAAGAGGTAGAGAGACTACAACAACCAAAAAA
GAACTGGAAGGACACAGCTTCTGCTTTGGACTTGTTGAAACGGGAAAAAGAGGCCAA
GAACAAGAGTTCTGTCTTTACAGGAGGAATCCAGAACTTGAAAGGAAAAACCTGGAA
GAAAGACAGAACTGAAATCTAGACTTGAGAAATTGCTCACTCAAGTTAGAAATTTGCAA
TTTATGTCTGAAATGAAAGAACGAAGAATATAAACTTCAGCAGCAAAATCAACGAAGTA
AAAAATGAGAATGCAAACTTAAACAGCAGGTTGCAAGGAGTGAAGAGCAAAATATGTC
CCTAAATTTGAGACAGCTCAGTTAAAGGATCAATTAGAGGAAGTCTTGAAGTCAGATATT
ACCAAGGATACAAAAACGACACATTCTAATCTGCTCCCGGATTGCTCACCTTGTGAAGAG
AGGCTGAATCCTGCAGATATAAAAAGAGCTTCTCAGCTGGCCTCCAAAATGCACAGTCTT
CTGGCTCTGATGGTGGGACTTCTCACATGCCAGGACATCATCAATTCTGATGCTGAACAT
TTCAAAGAGAGTGAGAAGGTTAGTGATATAATGCTGCAAAAACGAAGAGCCTCCATCTT
AAAAAGAAAACCTTTAGATAAAGAGCTACTGAAACATAAAGATAGAATCACAACCTTTAGA
GAGTTAATTGCTAAGGAAAAAGCATTTCAAGATCATGCTATTAAGGTCATAGACTGTGAT
TCAGATGAAGCCAAGAGTATCAGAGATGTACCTACCTTCTGGGAGCCAACTGGATAAG
TACCACAGTCTAAATGAGGAGCTTGATTTCTTGGTTACATCATATGAAGAAATCATTGAA
TGTGCTGACCAAAAGGCTTGCAATATCCCACTCCAGATTGCACATTTGGAAGAGAGAAAT
AAACATTTAGAGGATTTAATTAGAAAGCCAGAGAAAAAGCCAGAAAACCAAGATCAAAA
AGCTTAGAAAAATCATCCGAAGTCCATGACCATGTTGGCATATACTTGTACATATTTATA
GGATACAGGCTGGGCGCGGTGGCTGACGCCTGTAATCCAGCACTTTGGGAGGCCGAGGC
GGGTGGATCACTTGA

```

Restriction Sites: Please inquire
ACCN: NM_001170693

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001170693.1, NP_001164164.1</u>
RefSeq Size:	3074 bp
RefSeq ORF:	2475 bp
Locus ID:	285782
UniProt ID:	<u>Q8TC20</u>
Cytogenetics:	6p24.3