

Product datasheet for **SC329165**

GPRASP2 (NM_001184876) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GPRASP2 (NM_001184876) Human Untagged Clone
Tag:	Tag Free
Symbol:	GPRASP2
Synonyms:	DFNX7; GASP2
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001184876, the custom clone sequence may differ by one or more nucleotides

```

ATGACTGGGGCAGAGATTGAGCCTAGTGCCAGGCCAAGCCTGAAAAGAAGGCTGGGGAA
GAGGTTATCGCTGGGCTGAGAGAGAGAATGATGTCCCTCTGGTGGTCAGACCCAAGTT
AGGACCCAGGCAACTACTGGGGCAAGGCCAAAAGTGAAGCAAGTCTGTGCCTGCGGCA
AGGCCAAAAGTGAAGCCAAAGCAATGTCTGGGGCAAGGCCAAAAGTGAAGTCCAAGTA
ATGGGTGGTGAAGCCAAAAGCGAGGCTCAAGGAATCACAGGGGCCAGGCCAAAAGC
GATGCCAGGGCAGTAGGTGGCGCTCGTTCTAAAAGTGAAGCAAGCAATCCCTGGAGCA
AGGCCAAGGATGAGGCCAGGCATGGGCCAGAGTGAATTTGGGACTGAAGCAGTGTCA
CAGGCAGAAGGAGTGTCCAGACTAATGCCGTTGCTTGGCCACTGGCCACTGCTGAGTCT
GGATCAGTTACTAACTAAGGGCTGTCTATGGATAGAGAACTAGTCAATGTGGATGCT
GAAACCTTCTGGCACCCAGGGTCAGAAAGGAATCCAGCCCTGGTTGGACCAGGGGAG
GAGACTAATATGGGGTCTTGGTGTATTCCAGGCCAGGGCCAGAGAGGAGGCCTCTAAT
GAGTCTGGGTTCTGGTCAGCAGATGAGACCTCTACAGCGTCTTCTTTCTGGACTGGAGAA
GAGACAAGTGTGAGTATGAGCCAGGGAAGAGTCCAATACCAGGTCCAGGCACAGGGCT
AAACATCAGACTAATCCAGGTCCAGGCCAGATCCAAGCAAGAAGCCTATGTTGATTCC
TGGTCTGGATCTGAGGATGAGGCCAGCAACCCATTCTCCTTCTGGGTTGGAGAAAATACC
AATAACTTGTTCAGGCCAGAGTCAGGGAGGAGGCAAATATCAGGTCCAAGCTCAGGACA
AATAGAGAAGATTGTTTTGAATCTGAGTCTGAAGATGAGTTCTATAAGCAGTCTGGGTT
TTGCTGGAGAAGAGGCCAATAGTAGATTCAGGCACAGAGACAAAGAAGATCCTAATACT
GCCTTGAAACTCAGGGCCAGAAAGATGTTGACAGTGTAGGGTCAAACAAGAACCAGG
TTTGAGGAGGAAGTCATTATTGGGTCTGGTCTGGGCAGAAAAGAGGCCAGTTGGAG
GGTGGAGCTTCAGCAATCTGGAATCTGAGCCAGGAACTGAGGAGGGGGCCATTGGCGGA
TCCGCGTACTGGGCTGAGGAAAAGTCCAGTTTGGGGGCTGTGGCCAGAGAAGAGGCCAAG
CCGGAGTCTGAAGAAGAGGCCATATTTGGTCTGGTCTGGGACAGAGATGAGGCCTGC
TTTGACCTAAATCCCTGTCTGTGTACAAGTCAAGTGTAGGTTTCAAGATGCAGCTGAG
GAGCTTAATGCATCCTCCAGGCCCAAACCTGGGACGAGGTCACTGTTGAATTCAAACCT
GGTCTTTTTATGGGGTTGGCTCCGATCCACAAGCCCTTTGGAATCCCGAAGAGGCT
TCTGAAATGCTTGAAGCAAAGCCAAAGAACCTGGAACCTAGCCAGAGGAGAAGAGCAG
GAATCTTTGCTTCAAGCTGATCAGCCTAGTCTGAGTTTCAATTTCAAGTATGATCCTTCC
TACCGGTCAAGTCCGGGAAATTCGAGAGCATCTTAGGGCCAGGGAGAGTGCAGAGTCTGAG
AGTTGGTCTGCAGCTGCATACAATGTGAGCTGAAAATGGTTCTGAAGAGTTTGAAGAA
TTCTTTTATTAATGACAAAATTCGGGATCCTTTTATTCATGAAATATCTAAAATTGCA
ATGGGTATGAGAAGTGTCTCAATTTACCCGAGATTTTCAATTCGAGATTCAGGTGTTGTC
TCATTTATGAAACCTTGCTTAATATCCATCCTCTAGAGTTAGGACAAGTTTTTTGGAA
AATATGATTCACATGGCTCCACCTTATCCAAATCTAAACATGATTGAGACATTCATATGT
CAAGTGTGTGAGGAAACCTTGCACATAGTGTGGATTCCCTTGGAGCAGTCACTGGAATA
AGGATGCTTAGACACCTCACTATGACTATTGACTATCACACACTGATTGCCAACTATATG
TCCGGGTTTCTCTCTTATTAACCACAGCCAATGCGGAGAACGAAGTTTCAAGTTCTGAAA
ATGCTATTGAATTTGTCTGAAAATCCTGTGTGGCAAAAAAACTATTCAGTGCCAAAGCT
CTTTCAATATTTGTGGTCTCTTTAACATAGAAGAGACAAATGATAATATTCAAATTTGTT
ATTAATAATGTTTCAGAATATCAGTAACATTATAAAAAGTGGAAAGATGCTCCTAATTGAT
GATGATTTCAAGTCTTGGCCGCTTATTTCTGCATTTCTGTAATTTGAGGAGTTAGCTAAG
CAACTACAAGCCCAAATAGACAACCAAATGATCCTGAGGTGGGACAACAAAGTTAA

```

Restriction Sites: Please inquire
ACCN: NM_001184876
Insert Size: 3857 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001184876.1, NP_001171805.1</u>
RefSeq Size:	3857 bp
RefSeq ORF:	2517 bp
Locus ID:	114928
UniProt ID:	<u>Q96D09</u>
Cytogenetics:	Xq22.1
Protein Families:	Druggable Genome
Gene Summary:	<p>The protein encoded by this gene is a member of a family that regulates the activity of G protein-coupled receptors (GPCRs). The encoded protein has been shown to be capable of interacting with several GPCRs, including the M1 muscarinic acetylcholine receptor and the calcitonin receptor. Several transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, May 2010]</p> <p>Transcript Variant: This variant (5) differs in the 5' UTR compared to variant 3. All five variants encode the same protein.</p>