

Product datasheet for **SC329166**

GPRASP2 (NM_001184874) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GPRASP2 (NM_001184874) Human Untagged Clone
Tag:	Tag Free
Symbol:	GPRASP2
Synonyms:	DFNX7; GASP2
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001184874, the custom clone sequence may differ by one or more nucleotides

```

ATGACTGGGGCAGAGATTGAGCCTAGTGCCAGGCCAAGCCTGAAAAGAAGGCTGGGGAA
GAGGTTATCGCTGGGCTGAGAGAGAGAATGATGTCCCTCTGGTGGTCAGACCCAAGGTT
AGGACCCAGGCCAACTACTGGGGCAAGGCCAAAAGTCTGAGACCAAGTCTGTGCCTGCGGCA
AGGCCAAAAGTCTGAGGCCAAGCAATGTCTGGGGCAAGGCCAAAAGTCTGAGGTCCAAGTA
ATGGGTGGTGAAGCAAGCAAGCAAGCAAGCAAGCAAGCAAGCAAGCAAGCAAGCAAGCA
GATGCCAGGGCAGTAGGTGGCGCTCGTTCTAAAAGTCTGATGCCAAGCAATCCCTGGAGCA
AGGCCAAGGATGAGGCCAGGCATGGGCCAGAGTGAATTTGGGACTGAAGCAGTGTCA
CAGGCAGAAGGAGTGTCCAGACTAATGCCGTTGCTTGGCCACTGGCCACTGCTGAGTCT
GGATCAGTTACTAACTAAGGGCCTGTCTATGGATAGAGAACTAGTCAATGTGGATGCT
GAAACCTTCTGGCACCCAGGGTCAGAAAGGAATCCAGCCCTGGTTTGGACCAGGGGAG
GAGACTAATATGGGGTCTTGGTGTCTATCCAGGCCAGGGCCAGAGAGGAGGCCTCTAAT
GAGTCTGGGTTCTGGTCAGCAGATGAGACCTCTACAGCGTCTTCTTTCTGGACTGGAGAA
GAGACAAGTGTGAGTATGAGCCAGGGAAAGAGTCCAATACCAGGTCCAGGCACAGGGCT
AAACATCAGACTAATCCAGGTCCAGGCCAGATCCAAGCAAGAAAGCCTATGTTGATTCC
TGGTCTGGATCTGAGGATGAGGCCAGCAACCCATTCTCCTTCTGGGTTGGAGAAAATACC
AATAACTTGTTCAGGCCAGAGTCAGGGAGGAGGCAAATATCAGGTCCAAGCTCAGGACA
AATAGAGAAGATTGTTTTGAATCTGAGTCTGAAGATGAGTTCTATAAGCAGTCTGGGTT
TTGCTGGAGAAGAGGCCAATAGTAGATTCAGGCACAGAGACAAAGAAGATCCTAATACT
GCCTTGAAACTCAGGGCCAGAAAGATGTTGACAGTGTAGGGTCAAACAAGAACCAGG
TTTGAGGAGGAAGTCATTATTGGGTCTGGTCTGGGCAGAAAAGAGGCCAGTTGGAG
GGTGGAGCTTCAGCAATCTGGAATCTGAGCCAGGAACTGAGGAGGGGGCCATTGGCGGA
TCCGCGTACTGGGCTGAGGAAAAGTCCAGTTTGGGGGCTGTGGCCAGAGAAGAGGCCAAG
CCGGAGTCTGAAGAAGAGGCCATATTTGGTCTGGTCTGGGACAGAGATGAGGCCTGC
TTTGACCTAAATCCCTGTCTGTGTACAAGGTCAGTGTAGGTTTCCAGAGATGCAGCTGAG
GAGCTTAATGCATCCTCCAGGCCCAAACCTGGGACGAGGTCAGTGTGAATTCAAACCT
GGTCTTTTTATGGGGTTGGCTCCGATCCACAAGCCCTTTGGAATCCCGAAGAGGCT
TCTGAAATGCTTGGGCAAGGCCAAGAACCTGGAACCTAGCCAGAGGAGAAGAGCAG
GAATCTTTGCTTCCAGCTGATCAGCCTAGTCTGAGTTTACATTTTCAAGTATGATCCTTCC
TACCGGTGAGTCCGGGAAATTCGAGAGCATCTTAGGGCCAGGGAGAGTGCAGAGTCTGAG
AGTTGGTCTGCAGCTGCATACAATGTGAGCTGAAAATGGTTCTGAAGAGTTTGAAGAA
TTCTTTTATTAATGGACAAAATTCGGGATCCTTTTATTCATGAAATATCTAAAATTGCA
ATGGGTATGAGAAGTGTCTCAATTTACCCGAGATTTTATTTCGAGATTCAGGTGTTGTC
TCATTTATGAAACCTTGCTTAATATCCATCCTCTAGAGTTAGGACAAGTTTTTTGGAA
AATATGATTCACATGGCTCCACCTTATCCAAATCTAAACATGATTGAGACATTCATATGT
CAAGTGTGTGAGGAAACCTTGCACATAGTGTGGATTCCCTTGGAGCAGTACTGGAATA
AGGATGCTTAGACACCTACTATGACTATTGACTATCACACACTGATTGCCAACTATATG
TCCGGGTTTCTCTCTTATTAACCACAGCCAATGCGGAGAACGAAGTTTTCAGTTCTGAAA
ATGCTATTGAATTTGTCTGAAAATCCTGTGTGGCAAAAAAACTATTCAGTGCCAAAGCT
CTTTCAATATTTGTGGTCTCTTTAACATAGAAGAGACAAATGATAATATTCAAATTTGTT
ATTAATAATGTTTCAGAATATCAGTAACATTATAAAAAGTGGAAAGATGTCCTTAATTGAT
GATGATTTCAAGTCTTGGCCGCTTATTTCTGCATTTCTGTAATTTGAGGAGTTAGCTAAG
CAACTACAAGCCCAAATAGACAACCAAATGATCCTGAGGTGGGACAACAAAGTTAA

```

Restriction Sites: Please inquire
ACCN: NM_001184874
Insert Size: 3946 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001184874.1</u> , <u>NP_001171803.1</u>
RefSeq Size:	3946 bp
RefSeq ORF:	2517 bp
Locus ID:	114928
UniProt ID:	<u>Q96D09</u>
Cytogenetics:	Xq22.1
Protein Families:	Druggable Genome
Gene Summary:	<p>The protein encoded by this gene is a member of a family that regulates the activity of G protein-coupled receptors (GPCRs). The encoded protein has been shown to be capable of interacting with several GPCRs, including the M1 muscarinic acetylcholine receptor and the calcitonin receptor. Several transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, May 2010]</p> <p>Transcript Variant: This variant (3) represents the longest transcript. All five variants encode the same protein.</p>