

Product datasheet for **SC329188**

BCAR1 (NM_001170719) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	BCAR1 (NM_001170719) Human Untagged Clone
Tag:	Tag Free
Symbol:	BCAR1
Synonyms:	CAS; CAS1; CASS1; CRKAS; P130Cas
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001170719, the custom clone sequence may differ by one or more nucleotides

```

ATGGAGAACGTGCTGGCCAAAGCGCTCTATGACAATGTGGCCGAGTCCCCGGATGAGCTC
TCCTTCCGCAAGGGTGACATCATGACGGTGTGGAGCAGGACACGCAGGGCCTGGACGGC
TGGTGGCTCTGCTCGCTGCATGGGCGCCAGGGCATCGTGCCTGGGAACCGCCTCAAGATC
TTGGTGGGCATGTATGATAAGAAGCCAGCAGGGCCTGGCCCCGGCCCTCCCGCCACCCCG
GCCCCAGCTCAGCTGGCTCCATGCCCCAGCGCTCCGGCTCCCACTACACGCCCATG
CTCCCCAACACCTACCAGCCCCAGCCAGACAGCGTCTACCTGGTGCCCACTCCCAGCAAG
GCTCAGCAAGGCCTCTACCAAGTCCCGGGTCCCAGCCCTCAGTTCAGTCTCCCCAGCC
AAGCAGACATCCACCTTCTCGAAGCAGACACCCCATCACCGTTTCCCAGCCCGGCCACA
GACCTGTACCAGGTGCCCCAGGGCCTGGAGGCCCTGCCAGGATATTTACCAGGTGCCA
CCTTCTGCCGGGATGGGGCATGACATCTACCAGTCCCCCGTCCATGGACACACGCAGC
TGGGAGGGCACGAAGCCCCGGCAAAGGTGGTGGTGCCACCCCGTGGGCAGGGCTAT
GTATACGAGGCCCGCCAGCCGGAGCAGGACGAGTACGACATCCCAGCACCTGTGGCC
CCGGGGCCACAGGACATCTATGATGTGCCCGGTTCCGGGGCTGCTTCCAGCCAGTAT
GGCCAGGAGGTGTATGACACACCCCTTGGCTGTCAAGGGTCCCAATGGCCGAGACCCG
TTGCTGGAGGTGTATGAGTGGCCCCAGTGTGGAGAAGGGCTGCCACCGTCCAACCAC
CACGCAGTCTACGACGTTCTCCATCGGTGAGCAAGGATGTGCCCGATGGCCCACTGCTG
CGTGAGGAGACCTACGATGTGCCCCCGCCTTCGCCAAGGCCAAGCCCTTGACCCGGCC
CGCACCCCACTGGTACTGGTGCGCCCCCTCCAGACTCCCCGCCGGCCGAGGACGTGAT
GACGTGCCGCCCCCGCTCCTGACCTCTACGACGTGCCCCCTGGCTTGCGGCGCCTGGC
CCGGGCACCTGTACGATGTGCCCGTGAACGGGTGCTTCTCCTGAGGTGGCTGATGGT
GGCGTGGTGCACAGTGGTGTATGCGGTGCCTCCCCAGCTGAACGTGAAGCCCCGCA
GAGGGCAAAGCGCCTGTCCGCCCTCAGCACCGGCAGCACAGCAGCAGCCAGTCTGCGTCC
TCCTTGGAGGTGGCAGGGCCGGCCGGGAACCCCTGGAGCTGGAAGTTGCTGTGGAGGCC
CTGGCACGGTGCAGCAGGGTGTGAGCGCCACCGTTGCCACCTTCTGGACCTGGCAGGC
AGCGCCGGTGCAGTGGGAGTGGCGTAGCCCTCTGAGCCACAGGAGCCGCTGGTGCAG
GACCTGCAGGCTGCTGTGGCCGCTGTCCAGAGTGGCGTCCAGAGCTGTTGGAGTTTGGC
CGCAGCGCGTGGCAATGCTGCCACACATCTGACCGTGCCTGCATGCCAAGCTTAGC
CGGCAGCTGCAGAAGATGGAGGACGTGCACCAGACGCTGGTGGCACATGGTCAGGCCCTC
GACGCTGGCCGGGAGGCTCTGGAGCCACCTTGGAGACCTGGACCGCTGGTGGCTGC
TCGCGGGCTGTGCCGAGGACGCCAAGCAGCTGGCCTCCTTCTGCACGGCAATGCCTCA
CTGCTCTTCAGACGGACCAAGGCCACTGCCCGGGGCTGAGGGGGTGGCACCCCTGCAC
CCCAACCCCACTGACAAGACCAGCAGCATCCAGTACGACCCCTGCCCTACCCCTAAG
TTCACCTCCCAGGACTCGCCAGATGGGCAGTACGAGAACAGCGAGGGGGGCTGGATGGAG
GACTATGACTACGTCCACCTACAGGGGAAGGAGGAGTTTGAAGAAGCCAGAAGGAGCTG
CTGGAAAAGGGCAGCATCACGGCAGGGCAAGAGCCAGCTGGAGTTGCAGCAGCTGAAG
CAGTTTGAACGACTGGAACAGGAGGTGTACGGCCATAGACCAGACCTGGCCAAGTGG
ACGCCAGCCCAACCCCTGGCCCCGGGGCAACAGGGCCCTGGGGCCCTCGGACCGGCAG
CTGCTGCTCTTCTACCTGGAGCAGTGTGAGGCCAACCTGACCACACTGACCAACGCCGTG
GACGCCTTCTTTACCGCGTGGCCACCAACCAGCCGCCAAGATCTTTGTGGCGCACAGC
AAGTTCGTATCCTCAGCGCCACAAGCTGGTGTTCATCGGGGACACACTGTACGGCAG
GCCAAGGTGCTGACGTGCGCAGCCAGGTGACCCACTACAGCAACCTGCTGTGCGACCTC
CTGCGCGGCATCGTGGCCACCACCAAGGCCGCTGCCTTGCAGTACCCATCGCCTTCCGG
GCCAGGACATGGTGGAGAGGGTCAAGGAGCTGGGCCACAGCACCCAGCAGTTCGCCCG
GTCCTAGGCCAGCTGGCAGCCGCTGA

```

Restriction Sites: Please inquire
ACCN: NM_001170719

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001170719.1</u> , <u>NP_001164190.1</u>
RefSeq Size:	3372 bp
RefSeq ORF:	2607 bp
Locus ID:	9564
UniProt ID:	<u>P56945</u>
Cytogenetics:	16q23.1
Protein Families:	Druggable Genome
Protein Pathways:	Chemokine signaling pathway, Focal adhesion, Leukocyte transendothelial migration, Regulation of actin cytoskeleton

Gene Summary:

The protein encoded by this gene is a member of the Crk-associated substrate (CAS) family of scaffold proteins, characterized by the presence of multiple protein-protein interaction domains and many serine and tyrosine phosphorylation sites. The encoded protein contains a Src-homology 3 (SH3) domain, a proline-rich domain, a substrate domain which contains 15 repeat of the YxxP consensus phosphorylation motif for Src family kinases, a serine-rich domain, and a bipartite Src-binding domain, which can bind both SH2 and SH3 domains. This adaptor protein functions in multiple cellular pathways, including in cell motility, apoptosis and cell cycle control. Dysregulation of this gene can have a wide range of effects, affecting different pathways, including cardiac development, vascular smooth muscle cells, liver and kidney function, endothelial migration, and cancer. [provided by RefSeq, Sep 2017]

Transcript Variant: This variant (7) has an alternate 5' sequence, as compared to variant 1. The resulting isoform (7) is shorter and has a different N-terminus, as compared to isoform 1.

Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.