

## Product datasheet for **SC329208**

### Alpha Dystroglycan (DAG1) (NM\_001177642) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Alpha Dystroglycan (DAG1) (NM_001177642) Human Untagged Clone
Tag:	Tag Free
Symbol:	DAG1
Synonyms:	156DAG; A3a; AGRNR; DAG; LGMDR16; MDDGA9; MDDGC7; MDDGC9
Mammalian Cell Selection:	None
Vector:	<u><a href="#">pCMV6-XL5</a></u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >NCBI ORF sequence for NM\_001177642, the custom clone sequence may differ by one or more nucleotides

```

ATGAGGATGTCTGTGGGCCTCTCGCTGCTGCTGCCCTCTGGGGAGGACCTTTCTCCTC
CTGCTCTCTGTGGTTATGGCTCAGTCCCCTGGCCAGTGAACCCCTCAGAGGCTGTCAGG
GACTGGGAAAACAGCTTGGGCATCCATGCACTCAGTGTCTCAGACCTCCACGAGGCT
GTTCCACAGTGGTTGGCATTCTGATGGCACGGCTGTCGTGGGCGCTCATTTCGAGTG
ACCATTCCAACAGATTTGATTGCCTCCAGTGGAGATATCATCAAGGTATCAGCGGCAGGG
AAGGAGGCTTTGCCATCTTGGCTGCACTGGGACTCACAGAGCCACACCCTGGAGGGCCTC
CCCCTTGACACTGATAAGGGTGTGCATTACATTTAGTGAGCGCTACACGGCTGGGGGCC
AACGGGAGCCACATCCCCAGACCTCCAGTGTGTTCTCCATCGAGGTCTACCCTGAAGAC
CACAGTGAGTGCAGTGGTGGAGACAGCCTCCCAGACCCTGGTGGAGTGGTATCATCT
GCCTGTGCTGCGGATGAACCTGTGACTGTTTTGACGGTATTTGGATGCCGACCTACC
AAGATGACCCCAAAGCAAAGGATTGACCTCCTGCACAGGATGCGGAGCTTCTCAGAAGTA
GAGCTTACAACATGAAATAGTGCCGGTGGTGAATAACAGACTATTTGACATGTCGGCC
TTCATGGCTGGCCCGGAAATGCAAAAAGGTGGTGGAGAAATGGGGCCCTTCTCTCTGG
AAGCTGGGCTGCTCCCTGAACCAGAACAGTGTGCCTGACATTCATGGTGTAGAGGGCCCT
GCCAGGGAGGGCGCAATGTCTGCTCAGCTTGGCTACCCTGTGGTGGGTGGCACATCGCC
AATAAGAAGCCCCCTCTCCCAAACGCGTCCGGAGGCAGATCCATGCTACACCCACACCT
GTCAGTGCATTGGGCCCCAACACGGCTATCCAGGAGCCCCATCCAGGATCGTGCCA
ACCCCCACATCTCAGCCATTGCTCCTCCAACAGAGACCATGGCTCCTCCAGTACAGGAT
CCTGTTCTGGGAAACCCAGGTACCATCCGGACTCGAGGCGCCATTATTCAAACCCCA
ACCTAGGCCCATCCAGCCTACTCGGGTGTGAGAAGTGGCACACAGTTCTGGCCAG
ATTCGCCAACGATGACCATTCCTGGCTATGTGGAGCCTACTGCAGTTGCTACCCCTCCC
ACAACCACCAACAAGAAGCCACGAGTATCCACACCAAAACAGCAACGCCTTCAACTGAC
TCCACCACCAACGACTCGCAGGCCAACCAAGAAACCGGACACCCCGGCCAGTGCC
CGGGTACCACCAAAAGTTTCCATCACCAGATTGAAAAGTGCCTCACCGCCTACTCGTATT
CGCACCACCAAGTGGAGTGCCCGTGGCGGAGAACCAACAGCGCCAGAGCTCAAG
AACCATATTGACAGGGTAGATGCCTGGGTGGCACCTACTTTGAGGTGAAGATCCCGTCA
GACACTTTCTATGACCATGAGGACACCACCTGACAAGCTGAAGCTGACCTGAAACTG
CGGAGCAGCAGCTGGTGGGCGAGAAGTCTGGGTACAGTTCAACAGCAACAGCCAGCTC
ATGTATGGCTTCCCACAGCAGCCACGTGGGCAACACAGATATTTGCATGCCACA
GACAAGGGGGGCTGTGCGCTGTGGATGCCTTCGAGATCCACGTCCACAGGCGCCCCAA
GGGATAGGGCTCCTGCAAGGTTCAAGGCCAAGTTTGTGGGTGACCCGGCACTGGTGTG
AATGACATCCACAAGAAGATTGCCTTGGTAAAGAAACTGGCCTTCGCCTTTGGAGACCGA
AACTGTAGCACCATCACCTGCAGAATATCACCCGGGGTCCATCGTGGTGGAAATGGACC
ACAACACACTGCCCTTGGAGCCCTGCCCAAGGAGCAGATCGTGGGCTGAGCCGCGG
ATCGCTGAGGATGATGAAAACTCGGCCTGCCTTCTCCAACGCCCTAGAGCCTGACTTT
AAGGCCACAAGCATCACTGTGACGGGCTCTGGCAGTTGTGGCACCTACAGTTTATCCCT
GTGGTACCACCCAGGAGAGTGCCCTCAGAGGCGCCGCCACAGAAGTGCCTGACAGGGAC
CCTGAGAAGAGCAGTGAGGATGATGTCTACCTGCACACAGTATTCCGGCCGTGGTGGTC
GCAGCCATCCTGCTCATTGCTGGCATCATTGCCATGATCTGCTACCGCAAGAAGCGGAAG
GGCAAGCTTACCCTTGGAGACCAGGCCACCTTCATCAAGAAGGGGGTGCCTATCATCTTT
GCAGACGAACTGGACGACTCCAAGCCCCACCCTCCTCCAGCATGCCACTATTCTGCAG
GAGGAGAAGGCTCCCCTACCCCTCCTGAGTACCCCAACCAGAGTGTGCCGAGACCACT
CCTCTGAACCAGGACCCATGGGAGAGTACACGCCCTGCGGGATGAGGATCCCAATGCG
CCTCCCTACCAGCCCCACCGCCTTACAGCACCCATGGAGGGCAAGGGCTCCCGTCCC
AAGAACATGACCCCATACCGGTACCTCCTCCCTATGTCCCACCTTAA

```

**Restriction Sites:** Please inquire  
**ACCN:** NM\_001177642  
**Insert Size:** 5595 bp

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001177642.1</a></u> , <u><a href="#">NP_001171113.1</a></u>
<b>RefSeq Size:</b>	5595 bp
<b>RefSeq ORF:</b>	5595 bp
<b>Locus ID:</b>	1605
<b>UniProt ID:</b>	<u><a href="#">Q14118</a></u>
<b>Cytogenetics:</b>	3p21.31
<b>Protein Families:</b>	Druggable Genome, Secreted Protein, Transmembrane
<b>Protein Pathways:</b>	Arrhythmogenic right ventricular cardiomyopathy (ARVC), Dilated cardiomyopathy, ECM-receptor interaction, Hypertrophic cardiomyopathy (HCM), Viral myocarditis
<b>Gene Summary:</b>	<p>This gene encodes dystroglycan, a central component of dystrophin-glycoprotein complex that links the extracellular matrix and the cytoskeleton in the skeletal muscle. The encoded preproprotein undergoes O- and N-glycosylation, and proteolytic processing to generate alpha and beta subunits. Certain mutations in this gene are known to cause distinct forms of muscular dystrophy. Alternative splicing results in multiple transcript variants, all encoding the same protein. [provided by RefSeq, Nov 2015]</p> <p>Transcript Variant: This variant (11) differs in the 5' UTR compared to variant 1. All variants encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>