

Product datasheet for **SC329217**

DFNB31 (WHRN) (NM_001173425) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DFNB31 (WHRN) (NM_001173425) Human Untagged Clone
Tag:	Tag Free
Symbol:	DFNB31
Synonyms:	CIP98; DFNB31; PDZD7B; USH2D; WI
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001173425, the custom clone sequence may differ by one or more nucleotides

```

ATGAACGCGCCGCTGGACGGCCTGTCGGTGAGCTCGTCTCCACCGGCTCGCTGGGCTCG
GCGGCCGGGGCGGGCGGGCGGGGGCGGGGGCTGCGGTTACTGTCTGCCAACGTGCGC
CAGCTGCACCAAGCGCTGACCGCGCTGTGAGCGAGGCGGAGCGGGAGCAGTTACCCAC
TGCTTGAACGCTTACCACGCGCGCCGCAACGTCTTCGACCTGGTGCACCCCTGCGCGT
CTGCTGGACAGTCCGGTCAAGCGCGCCTGTGCCCATGCTTCGTCTGGTCATCCCGCG
TCCGACCAGTGTCTTTCGACCAATACACGGCCGAGGGCCTCTACCTGCCCGCCACCACC
CCCTACAGGCAGCCCGCTGGGGCGGCCCGACAGCGGGGCCAGGGGAGGTGCGCCTG
GTGAGTTTTCGGCGTGCCAAGGCCACGAGGGCTTGGGCTTACGATCCGTGGGGCTCG
GAGCACGGCGTGGGCATCTACGTGTCTCTGGTGAACAGGCTCTCTAGCTGAGAAGGAA
GGACTGCGGGTGGGGACAGATTCTGCGCGTCAACGACAAATCCCTGGCCCGGGTGACC
CACGCGGAGGCCGTCAAGGCTCTGAAGGGCTCCAAGAAGCTGGTGTGTGTACTCA
GCAGGGCGCATCCCTGGGGGCTACGTACCAACCACATCTACACCTGGGTGGACCCGAG
GGCCGACGATCTCCCCACCTCGGGCCTGCCACGCCCACGGTGGTGCCTGAGGCGAG
CAGGAGGGTGACCGGAGGAGCACCTGCACCTCCTGCAAGGAGGGGATGAGAAAAAGGTG
AACCTGGTGTGGGGGACGGCCGGTCCCTGGGCCTCACGATCCGTGGGGGAGCTGAGTAC
GGCCTTGGCATTACATCACTGGCGTGGACCCAGGCTCTGAAGCAGAAGGCAGCGGGCTC
AAGTTGGGGACAGATTCTAGAAGTGAATGGCGGAGCTTTCTCAACATCCTACACGAC
GAGGCTGTGAGGCTGCTTAAGTCATCTCGGCACCTCATCTGACAGTGAAGGACGTGGG
AGGCTGCCCATGCCCGCACCCTGTGGACGAGACCAAGTGGATCGCCAGTTCGGGATC
AGGGAGACCATGGCGAACTCGGCAGGGTTTCTTGGCGATCTCACAAAGAAAGGAATAAC
AAGCCAGGATTTACAAGGGCCAGCCGGCTCCAGGTGACCCTGAGCAGCTGGGGAAAC
CAGACACGAGTGTCTGGAGGAGCAGGCTCGGCACCTGCTGAACGAGCAGGAACACGCC
ACCATGGCCTACTACCTGGATGAGTACCGTGGTGGCAGCGTCTCTGTGGAGGCCCTCGT
ATGGCCCTGTTCAAGCTGCTCAACACCCACGCCAAGTTCTCACTCCTCTGAGGTGAGA
GGCACCATTTCCCGCAAGACCTAGAACGCTTCGACCACCTGGTGTGAGGCGTGAGATT
GAGTCCATGAAGGCGGGCAGCCCCAGGCCCGGGGCTGGGGACACCTACTCCATGGTC
TCCTACAGTGACACGGTTTATCCACAGGCAGCCACGGCACCTCCACCACCGTCAGTCTG
GCCAGGAACACTCTGGACCTGGAGGAACTGGCGAGGCTGTCCAGGGCAATATCAACGCC
CTCCCAGATGTGTCCGTGGATGATGTGAGTCCACCTCCAGGGGCTGTCAAGCTTCAAG
CCTACTGCCTCGCCACCACCTCTGGCCCAAGGCAACGACCTCCCACTAGGCCAGCCAAGG
AAGCTGGGGAGAGAGGACCTCCAGCCACCTTCTCCATGCCTTCTGCTCGGGCACTGTC
TTCTCGGCTCCACAGAACCAGCCAGCCCGCAGCGGGCACCCACCCAGGGACCTCC
TCTGCACAGGACTTGCCCTTTCCTCCCATCTATGCCTCCGTCTCCCTGCCAACCCACG
TCCAAGAGGCGCTGGACGCCATCTGGCCCTGGTCAACCAACACCCATCGGCCCTTC
CCACGGGTCCAGTCAACCCCGCACCTGAAAAGCCCTCTGCAGAGGCCACAGTGGCTGG
GGCTGCCCTTCTGCCCCATCACCTCTGGCCACCAGACCAGACAGGCACAAACCAGCAC
TTTGTGATGGTGGAGGTCCACCGCCCGACAGCGAGCCAGACGTCAATGAAGTGAAGGCG
CTGCCCCAGACGCGCACCTCTACGCTCTCCAGCTCTCGGACAGCGGGCAGACTCTAAGC
GAGGACAGTGGTGTGGATGCTGGCGAGGAGAGGCCAGCGCCCAAGCCGAGGAAGGCGAG
TCGGTGTCCACCAAGAGCAGGAGTAGCAAGGAGTGCCTCGGAACGAGAGGCCACAGAT
GGGGCCAAACAACCGCTGGACTTCTGGAGCCACGTCCACTCTGGTCCGTGTGAAGAAA
AGTGCGGCCACCCTGGGCATCGCCATCGAGGGTGGCGCCAACACCCGCCAGCCCCGCT
AGGATTGTCACTATTAGAGAGGGCGCTCAGCTCACAACGTGGGCAGCTCAAGTGGGC
CACGTGATTCTGGAAGTGAATGGGCTGACGCTTCGGGGCAAGGAGCACCGGGAGGCCGCC
CGCATTATCGCCGAGGCCTTCAAGACTAAGGACCGTGACTACATTGACTTCTGGTCACT
GAGTTCAATGTGATGCTCTAG

```

Restriction Sites: Please inquire
ACCN: NM_001173425

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001173425.1, NP_001166896.1</u>
RefSeq Size:	4076 bp
RefSeq ORF:	2721 bp
Locus ID:	25861
UniProt ID:	<u>Q9P202</u>
Cytogenetics:	9q32
Gene Summary:	<p>This gene is thought to function in the organization and stabilization of stereocilia elongation and actin cytoskeletal assembly, based on studies of the related mouse gene. Mutations in this gene have been associated with autosomal recessive non-syndromic deafness and Usher Syndrome. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2016]</p> <p>Transcript Variant: This variant (3) contains an alternate, in-frame splice site in the 3' coding region, compared to variant 1. The encoded protein (isoform 3) is shorter than isoform 1.</p>