

Product datasheet for **SC329266**

NALP2 (NLRP2) (NM_001174082) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NALP2 (NLRP2) (NM_001174082) Human Untagged Clone
Tag:	Tag Free
Symbol:	NALP2
Synonyms:	CLR19.9; NALP2; NBS1; PAN1; PYPAF2
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001174082, the custom clone sequence may differ by one or more nucleotides

```

ATGGTGTCTTCGCGCAGATGGGCTTCAACCTGCAGGCTCTCCTGGAGCAGCTCAGCCAG
GATGAGTTGAGCAAGTTCAAGTATCTGATCACGACCTTCTCCCTGGCACACGAGCTCCAG
AAGATCCCCCACAAAGGAGGTAGACAAGGCTGATGGGAAGCAACTGGTAGAAATCCTCACC
ACCCATTGTGACAGCTACTGGGTGGAGATGGCGAGCCTCCAGGCTTTGAAAAGATGCAC
CGAATGGATCTGTCTGAGAGAGCAAAGGATGAAGTCAGAGAAGCAGCTTTGAAATCCTTT
AATAAAGGAAGCCTCTATCATTAGGGATAACACGAAAGAACGACCACCTTAGACGTG
GACGAAATGCTGGAGCGCTTCAAACAGAAAGCACAAGACAATAGGTGCAGGTAT
ATATTGAAGACGAAGTTCGGGAGATGTGAAGAGCTGGCCTGGAGATAGCAAAGAGGTC
CAGGTTATGGCTGAGAGATACAAGATGCTGATCCATTACGCAACCCAGGGTGCTTCCC
GGGCCCTTCTCATACACGGTGGTGTGTATGGTCTGCAGGCCTTGGGAAAACACGCTG
GCCAGAAACTAATGCTAGACTGGGCAGAGGACAACCTCATCCAAAATCAAATATGCG
TTCTACCTCAGCTGCAGGGAGCTCAGCCGCTGGGCCCGTGCAAGTTTGCAGAGCTGGTC
TTCAGGGACTGGCCTGAATTGCAGGATGACATTCCACACATCCTAGCCCAAGCACGGAAA
ATCTTGTTCGTGATTGACGGCTTTGATGAGCTGGGAGCCGACCTGGGGCGCTGATCGAG
GACATCTCGGGGACTGGGAGAAGAAGAAGCCGGTGCCCGTCTCCTGGGGAGTTTGCTG
AACAGGGTGATGTTACCCAAGGCCGCCCTGCTGGTCACCACGCGGCCAGGGCCCTGAGG
GACCTCCGGATCCTGGCGGAGGACCGATCTACATAAGGGTGGAGGGCTTCTGGAGGAG
GACAGGAGGGCCTATTTCTGAGACACTTTGGAGACGAGGACCAAGCCATGCGTGCCTTT
GAGCTAATGAGGAGCAACGCGGCCCTGTTCCAGCTGGGCTCGGCCCCCGGGTGTGCTGG
ATCGTGTGCACGACTCTGAAGCTGCAGATGGAGAAGGGGGAGGCCCGTCCCCACCTGC
CTACCCGCACGGGGCTGTTCTGCGTTTCTCTGCAGCCGGTCCCGCAGGGCGCACAG
CTGCGGGGCGGCTGCGGACGCTGAGCCTCCTGGCCGCGCAGGGCCTGTGGGCGCAGACG
TCCGTGCTTACCAGAGGATCTGGAAAGGCTCGGGGTGCAGGAGTCCGACCTCCGTCTG
TTCCTGGACGGAGACATCCTCCGCCAGGACAGAGTCTCAAAGGCTGCTACTCCTTCATC
CACCTCAGCTTCCAGCAGTTTCTACTGCCCTGTTCTACACCCTGGAGAAGGAGGAGGAA

```



[View online »](#)

GAGGATAGGGACGGCCACACCTGGGACATTGGGGACGTACAGAAGCTGCTTTCCGGAGTA
 GAAAGACTCAGGAACCCCGACCTGATCCAAGCAGGCTACTACTCCTTTGGCCTCGCTAAC
 GAGAAGAGAGCCAAGGAGTTGGAGGCCACTTTTGGCTGCCGGATGTCACCCGGACATCAA
 CAGGAATTGCTGCGATGCGACATAAGTTGTAAGGGTGGACATTCAACGGTGACAGACCTG
 CAGGAGCTCCTCGGCTGTCTGTACGAGTCTCAGGAGGAGGAGCTGGTGAAGGAGGTGATG
 GCTCAGTTCAAAGAAATATCCCTGCACTTAAATGCAGTAGACGTTGTGCCATCTTCATTC
 TGGCTCAAGCACTGTCGAAACCTGCAGAAAATGCACTGCAGGTAATAAAGGAGAATCTC
 CCGGAGAATGTCAGTGCCTGAATCAGACGCCGAGGTTGAGAGATCCCAGGATGATCAG
 CACATGCTTCCTTTCTGGACGGACCTTTGTTCCATATTTGGATCAAATAAGGATCTGATG
 GGTCTAGCAATCAATGATAGCTTCTCAGTGCCTCCCTAGTAAGGATCCTGTGTGAACAA
 ATAGCCTCTGACACCTGTCATCTCCAGAGAGTGGTGTCAAACATTTCCCGAGCTGAT
 GCTCATCGGAACCTCTGCCTAGCTCTTCGAGGTCACAAGACTGTAACTGTGACCCCTT
 CAAGGCAATGACCAGGATGATATGTTTCCCGCATTGTGTGAGGTCTTGAGACATCCAGAA
 TGTAACCTGCGATATCTCGGTTGGTGTCTTGTCCGCTACCACTCAGCAGTGGGCTGAT
 CTCTCCTTGGCCCTGAAGTCAACCAGTCCCTGACGTGCGTAAACCTCTCCGACAATGAG
 CTCTGGATGAGGGTGCTAAGTTGCTGTACACAACCTTTGAGACACCCCAAGTGCTTTCTG
 CAGAGGTTGTCGTTGGAAAACCTGCACCTTACAGAAGCCAATTGCAAGGACCTTGCTGCT
 GTGTTGGTTGTCAGCCGGGAGCTGACACACCTGTGCTTGGCCAAGAACCCATTGGGAAT
 ACAGGGGTGAAGTTTCTGTGTGAGGGCTTGAGGTACCCCGAGTGTAACTGCAGACCTTG
 GTGCTTTGGAACCTGCGACATAACTAGCGATGGCTGCTGCGATCTCACAAGCTTCTCCAA
 GAAAAATCAAGCCTGTTGTGTTGGATCTGGGGCTGAATCACATAGGAGTTAAGGGAATG
 AAGTTCCTGTGTGAGGCTTTGAGGAAACCACTGTGCAACTTGAGATGTCTGTGGTTGTGG
 GGATGTTCCATCCCTCCGTTTCAGTTGTGAAGACCTCTGCTTGCCTCAGCTGCAACCCAG
 AGCCTCGTCACTCTGGACCTGGGTGAGAAATCCCTTGGGGTCTAGTGGAGTGAAGATGCTG
 TTTGAAACCTTGACATGTTCCAGTGGCACCCCTCCGGACACTCAGGTTGAAAATCGATGAC
 TTTAATGATGAACTCAATAAGCTGCTGGAAGAAAATAGAAGAAAAAACCACAACCTGATT
 ATTGATACTGAGAAACATCATCCCTGGGCAGAAAGGCCTTCTTCTCATGACTTCATGATC
 TGA

- Restriction Sites:** Please inquire
- ACCN:** NM_001174082
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
- OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001174082.1](#), [NP_001167553.1](#)

RefSeq Size: 3524 bp

RefSeq ORF: 3123 bp

Locus ID: 55655

UniProt ID: [Q9NX02](#)

Cytogenetics: 19q13.42

Gene Summary: This gene is a member of the nucleotide-binding and leucine-rich repeat receptor (NLR) family, and is predicted to contain an N-terminal pyrin effector domain (PYD), a centrally-located nucleotide-binding and oligomerization domain (NACHT) and C-terminal leucine-rich repeats (LRR). Members of this gene family are thought to be important regulators of immune responses. This gene product interacts with components of the I κ B kinase (IKK) complex, and can regulate both caspase-1 and NF- κ B (nuclear factor kappa-light-chain-enhancer of activated B cells) activity. The pyrin domain is necessary and sufficient for suppression of NF- κ B activity. An allelic variant (rs147585490) has been found that is incapable of blocking the transcriptional activity of NF- κ B. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Dec 2016]
Transcript Variant: This variant (3) lacks an alternate in-frame exon in the 5' coding region compared to variant 1. The encoded isoform (2) has the same N- and C-termini, but is shorter than isoform 1.