

Product datasheet for **SC329359**

TUT7 (NM_001185059) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TUT7 (NM_001185059) Human Untagged Clone
Tag:	Tag Free
Symbol:	TUT7
Synonyms:	PAPD6; TENT3B; ZCCHC6
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001185059, the custom clone sequence may differ by one or more nucleotides

```

ATGGGAGATACAGCAAAACCTTATTTTCGTGAAGCGCACTAAAGACCGGGGGACTATGGAT
GATGATGACTTCAGAAGGGGTCACCCCCAACAAAGATTATTTAATAATAGATGACCATGCT
AAAGGCCATGGCAGTAAATGGAAAAGGGCCTTCAAAAAAAGAAGATAACACCAGGGAAC
TATGGGAATACCCCCAGAAAAGGACCATGTGCTGTTTCCAGCAATCCATATGCATTTAA
AACCCAATCTACAGTCAACCCGCTTGGATGAATGACAGCCACAAAGATCAGAGTAAGAGA
TGGCTGTCTGATGAACATACTGGTAATTCAGACAACCTGGAGAGAATTCAAACCTGGACCT
AGAATTCCTGTTATAAACCGACAAAGAAAAGACTCCTTTCAAGAAAATGAAGATGGTTAT
AGGTGGCAAGACACAAGAGGCTGCAGAACTGTAAGACGACTGTTTCATAAAGACCTAACA
AGCCTAGAAACACGTCAGAAATGGAAGCAGGAAGTCCTGAAAACAAGAAGCAGAGGTCC
AGACCTAGGAAGCCACGGAAGACTAGAAATGAGGAAAATGAGCAGGATGGAGACTTGGAA
GGCCCTGTGATCGATGAGTCTGTACTTTCAACGAAGGAGCTGCTAGGCTTACAGCAGGCT
GAGGAGAGACTGAAGAGAGACTGCATTGACAGGCTAAAAAGGCGACACGAAACTACCT
ACAGCAAAGTACACCTGCAGACTCTGTGATGTTTTAATTGAATCCATTGCATTTGCCCAT
AAGCATATCAAGGAAAAGAGGCACAAGAAAAACATTAAGGAGAAGCAAGAGGAAGAGTTG
CTCACTACGTTACCCCCACCAACACCTCCCAGATAAATGCAGTTGGCATTGCCATTGAC
AAAGTGGTACAGGAATTTGGCTTACACAATGAGAACTTGAACAGAGGCTGGAAATTA
CGTATCATGGAAAATGTGTTCCAACACAAGTTACCAGATTGTTCCCTAAGATTATATGGG
TCATCCTGTAGCAGATTGGGTTTCAAAAAATTCGGATGTAACATTGACATCCAGTTTCCA
GCCATTATGTCTCAGCCAGATGTCCTCTTACTTGTTCAGAAATGTTTAAAGAACAGTGAC
TCCTTTATTGATGTTGATGCAGACTTCCATGCTAGGGTGCCAGTGGTGGTGTGCAGAGAA
AAGCAAAGTGGTCTTCTGTGTAAGTGAGCGCAGGAAATGAAAATGCTTGTCTGACAACA
AAGCATTTAACTGCCCTTGAAAACTAGAACCAAAGCTGGTTCCTTTGGTGATTGCATTT
AGGTACTGGGCAAAGCTTTGCAGTATAGATCGCCCTGAAGAAGGAGGTCTGCCACCTTAT
GTGTTTGCCTGATGGCCATTTTCTTTCTCAGCAGAGGAAAGAACCCCTTTTGCCTGTA
TATCTAGGATCATGGATTGAAGGATTCTCATTAAAGCAAACCTAGGGAATTTCAACCTTCAA
  
```


[View online »](#)

GACATTGAAAAAGATGTTGTGATCTGGGAACATACTGACAGTGCTGCAGGGGACACAGGC
 ATAAACAAAAGAGAGGCCACCAAGAGAAACGCCGATTAAAAGGGGACAGGTGTCTTAATA
 TTGGATGTGAAACACCAGCCTTCAGTACCAGTTGGGCAGCTCTGGGTGGAATTGCTGCGG
 TTCTATGCTTTAGAAATTAATTTGGCTGATTTAGTGATAAGTATTCGTGTCAAAGAATTG
 GTATCTCGGGAATTGAAGGATTGGCCAAAAAGCGCATTGCCATTGAAGATCCCTACTCT
 GTTAAAAAGAAATGTGGCAAGAACCCTAAATAGTCAACCTGTGTTGAATATATACTTCAT
 TGTTTAAGGACAACATACAAGTATTTTGCTCTTCCACACAAAATTACAAAATCCAGCCTT
 CTTAAAGCCTCTGAATGCAATTACATGTATTTTCAGAACATTCTAAAGAAGTAATAAATCAT
 CATCCAGATGTACAAACAAAAGATGATAAGCTCAAAAACCTCAGTTTTGGCCCAAGGTCTT
 GGTGCTACCAGTTCAGCTGCAAATACCTGTAAGGTACAGCCACTTACTCTTAAAGAGACT
 GCTGAAAGTTTTGGAAGCCACCAAAAGAAGAAATGGGAAATGAACACATCAGTGTCCAC
 CCTGAAAACCTCAGACTGTATCCAAGCAGATGTTAACTCTGATGATTACAAGGGTGATAAA
 GTATACCATCCAGAAACAGGAAGGAAAAACGAGAAAGAGAAAGTTGGAAGGAAGGCAAG
 CATCTGTTGACTGTTGATCAGAAACGTGGAGAGCATGTTGTCTGTGGCAGCACACGTAAT
 AATGAGTCAGAGAGCACTTTGGATTTAGAAGGCTTCCAAAATCCCACAGCTAAAGAGTGT
 GAGGGACTTGCCACTTTAGATAACCAAGGCTGATCTTGATGGAGAAAGTACAGAAGGTACT
 GAGGAACTAGAAGACTCTCTAAACCACTTTACCCACTCAGTACAGGGCCAGACATCAGAA
 ATGATTCCCTCTGATGAAGAGGAGGAGGACGACGAAGAAGAGGAGGAGGAAGAAGAACCT
 AGGCTCACCATTAAACCAAGGGAAGATGAAGATGGCATGGCTAATGAAGATGAGTTAGAC
 AACACCTACACTGGGTACAGGGGATGAGGACGCCCTATCTGAAGAGGATGATGAGTTAGGC
 GAAGCTGCTAAGTATGAAGACGTGAAAGAATGTGGAAAACATGTAGAAAGAGCTCTCCTA
 GTGGAACCTTAATAAAATAAGTCTCAAGGAAGAAAATGTATGTGAAGAAAAAATTCACCT
 GTGGATCAGTCTGATTTTTTTTATGAATTCAGTAACTTATCTTCACCAAGGCAAGTCT
 CCTACGCTAGTGTGCAGCTTATGCAAAACGAGAGGGTCATCTAAAGAAGGACTGTCCTGAA
 GACTTCAAAAGAATCCAGCTAGAACCTCTGCCACCATTACACCCAAGTTTTTAAATATC
 TTAGATCAAGTCTGTATCCAGTGTATAAGGATTTTTCTCCAACAATTATAGAAGATCAG
 GCTCGTGAACATATTGGCAAAACCTAGAAAGTTTCATAAGACAGGACTTTCAGGAACT
 AAATTGAGCCTGTTTGGCTCCTCCAAAATGGATTTGGGTTCAAACAGAGTGACCTTGAC
 GTCTGTATGACAATTAATGGACTTGAACTGCTGAGGGATTGGACTGTGTGAGAATATT
 GAAGAATTAGCAAGAGTCTCAGAAAACATTGAGTCTGAGAAACATCTTACCTATTACA
 ACAGCAAAGGTGCAATTTGTGAAGTTCTTCCATTTGAGAAGTGGTCTGGAAGTAGATATC
 AGTTTGTATAACACATTGGCCCTTCATAACACAAGGCTTTTATCTGCTTATTCGCCATT
 GATCCCAGAGTGAAGTATTTGTGCTATACCATGAAAGTATTTACAAAGATGTGTGATATT
 GGTGATGCATCTAGAGGCAGCTTATCATCGTATGCATATACTCTTATGGTGCTATATTTT
 CTCCAGCAGAGGAATCCACCAGTCATTCTGTCTTCAAGAGATATACAAAGGTGAAAAG
 AAACCTGAAATATTTGTTGATGGCTGGAATATTTATTTTTTTGATCAAATAGATGAAGTG
 CCTACCTATTGGTCAGAAATGTGGAAAAATACAGAATCTGTTGGGCAGTTATGGTTGGGC
 CTTCTTCGTTTCTACACAGAGGAATTTGATTTTAAAGAACATGTTATTAGCATCAGGAGA
 AAAAGTCTGCTTACAACTTTTAAGAAACAGTGGACCTCAAAATACATTGTTATTGAAGAT
 CCTTTGATTTGAATCATAATCTTGGAGCTGGATTATCAAGGAAAATGACAAATTTTATA
 ATGAAGGCTTTTATCAATGGTAGAAGATTTTGGTATTCTGTCAAGGGATTTCCAAAG
 GACTACCCCTCAAAAATGGAATACTTTTTTGTATCCAGATGTGTTAACTGAAGGAGAGCTG
 GCCCAAATGATAGATGTTGTGCAATTTGTGGAAAAATCGGACACTTCATGAAGGACTGT
 CCTATGAGGAGAAAAGTAAGACGGCGGCGAGATCAGGAAGATGCCCTGAACCAAAGATAC
 CCTGAGAACAAAGGAAAAAGAAGCAAGAGGACAAAGAAATTCACAACAAGTACACAGAA
 AGGGAGGTGTCAACAAAAGAAGATAAGCCCATACAGTGCACACCTCAGAAAGCCAAGCCA
 ATGCGGGCAGCTGCTGACCTGGGGAGGGAGAAGATCCTCAGGCCACCAGTAGAAAAATGG
 AAGAGACAGGATGACAAAGACTTAAGAGAAAAACGTTGTTTTATTGTGGAAGAGAAGGG
 CACATTAAGAAAGGAATGCCACAGTTTAAAGGCTCTTCAGGTAGCCTTTCCAGTAAATAT
 ATGACTCAGGGAAAAGCCTCAGCGAAGAGGACCCAGCAGGAATCATGA

Restriction Sites:

Please inquire

ACCN:	NM_001185059
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001185059.1, NP_001171988.1</u>
RefSeq Size:	5586 bp
RefSeq ORF:	4488 bp
Locus ID:	79670
UniProt ID:	<u>Q5VYS8</u>
Cytogenetics:	9q21.33

Gene Summary:

Uridyltransferase that mediates the terminal uridylation of mRNAs with short (less than 25 nucleotides) poly(A) tails, hence facilitating global mRNA decay (PubMed:19703396, PubMed:25480299). Essential for both oocyte maturation and fertility. Through 3' terminal uridylation of mRNA, sculpts, with TUT7, the maternal transcriptome by eliminating transcripts during oocyte growth (By similarity). Involved in microRNA (miRNA)-induced gene silencing through uridylation of deadenylated miRNA targets (PubMed:25480299). Also functions as an integral regulator of microRNA biogenesis using 3 different uridylation mechanisms (PubMed:25979828). Acts as a suppressor of miRNA biogenesis by mediating the terminal uridylation of some miRNA precursors, including that of let-7 (pre-let-7). Uridylated pre-let-7 RNA is not processed by Dicer and undergo degradation. Pre-let-7 uridylation is strongly enhanced in the presence of LIN28A (PubMed:22898984). In the absence of LIN28A, TUT7 and TUT4 monouridylate group II pre-miRNAs, which includes most of pre-let7 members, that shapes an optimal 3' end overhang for efficient processing (PubMed:25979828, PubMed:28671666). Add oligo-U tails to truncated pre-miRNAs with a 5' overhang which may promote rapid degradation of non-functional pre-miRNA species (PubMed:25979828). Does not play a role in replication-dependent histone mRNA degradation (PubMed:18172165). Due to functional redundancy between TUT4 and TUT7, the identification of the specific role of each of these proteins is difficult (PubMed:25979828, PubMed:25480299, PubMed:19703396, PubMed:22898984, PubMed:18172165, PubMed:28671666). TUT4 and TUT7 restrict retrotransposition of long interspersed element-1 (LINE-1) in cooperation with MOV10 counteracting the RNA chaperone activity of L1RE1. TUT7 uridylates LINE-1 mRNAs in the cytoplasm which inhibits initiation of reverse transcription once in the nucleus, whereas uridylation by TUT4 destabilizes mRNAs in cytoplasmic ribonucleoprotein granules (PubMed:30122351).[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (2) uses a different splice site in the 5' UTR, compared to variant 1. Variants 1 and 2 encode the same protein isoform.