

Product datasheet for **SC329437**

CSRP1 (NM_001193572) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: CSRP1 (NM_001193572) Human Untagged Clone
Tag: Tag Free
Symbol: CSRP1
Synonyms: CRP; CRP1; CSRP; CYRP; D1S181E; HEL-141; HEL-S-286
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC329437 representing NM_001193572.
Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
ATGCCGAAGTGGGGAGGAGGCAAGAAATGTGGGGTGTGTCAGAAGACGGTTTACTTTGCCGAAGAGGTT
CAGTGCGAAGGCAACAGCTTCCATAAAATCCTGCTTCTGTGCATGGTCTGCAAGAAGAATCTGGACAGT
ACCACTGTGGCCGTGCATGGTGAGGAGATTTACTGCAAGTCTGCTACGGCAAGAAGTATGGGCCAAA
GGCTATGGCTACGGGCAGGGCGCAGGCACCCTCAGCACTGACAAGGGGGAGTCGCTGGGTATCAAGCAC
GAGGAAGCCCCTGGCCACAGGCCACCACCAACCCCAATGCATCAAATTTGCCAGAAGATTGGTGCC
TCCGAGCGCTGCCCCGATGCAGCCAGGCAGTCTATGCTGCGGAGAAGGTGATTGGTGTCTGGGAAGTCC
TGGCATAAGGCCTGCTTTTCGATGTGCCAAGTGTGGCAAAGGCCTTGAGTCAACCACCTGGCAGACAAG
GATGGCGAGATTTACTGCAAAGGATGTTATGCTAAAAACTTCGGGCCCAAGGGCTTTGGTTTTGGGCAA
GGAGCTGGGGCCTTGGTCCACTCTGAGTGA
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

Plasmid Map:



ACCN: NM_001193572

Insert Size: 582 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001193572.1](#)

RefSeq Size: 1932 bp

RefSeq ORF: 582 bp

Locus ID: 1465

UniProt ID: [P21291](#)

Cytogenetics: 1q32.1

MW: 20.6 kDa

Gene Summary: This gene encodes a member of the cysteine-rich protein (CSRP) family. This gene family includes a group of LIM domain proteins, which may be involved in regulatory processes important for development and cellular differentiation. The LIM/double zinc-finger motif found in this gene product occurs in proteins with critical functions in gene regulation, cell growth, and somatic differentiation. Alternatively spliced transcript variants have been described. [provided by RefSeq, Aug 2010]
Transcript Variant: This variant (5) differs in the 5' UTR compared to variant 1. Variants 1, 4 and 5 encode the same isoform (1).