

Product datasheet for **SC329457**

TRPC5OS (NM_001195578) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: TRPC5OS (NM_001195578) Human Untagged Clone
Tag: Tag Free
Symbol: TRPC5OS
Synonyms: TRPC5-AS1
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC329457 representing NM_001195578.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

ATGGATTCTGTGTTAATTCATGTAAGTACTGATGGAAGTGTGCTGTGTAGCCAGTTAATAAGAATA
GCTGATGAGCTTTACAATTCATTCTACAAGTACAAGAAGTTCCTTATGTAGAAGAAATGGTAGAGCA
GAAGAGACTGAAGCAGATGCACCTCTCCCGAGGAGCCTTCGCTACCTGATCTCCCTGATCTCTCAGAC
TTAGACTCAATACTTACACCAAGAGAGGATGAAGACCTAATTTGATATAGATCAGGCTATGTTAGAC
ATGGATAACTTATGAAGATACAGTCTCTGGTATAAATGATGACTTAACAGGTGACTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001195578
Insert Size: 336 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq: [NM_001195578.1](#)



[View online »](#)

RefSeq Size: 2770 bp
RefSeq ORF: 336 bp
Locus ID: 100329135
UniProt ID: [A6NMA1](#)
Cytogenetics: Xq23
MW: 12.3 kDa