

Product datasheet for **SC329470**

TLCD5 (NM_001198675) Human Untagged Clone


Product data:

Product Type:	Expression Plasmids
Product Name:	TLCD5 (NM_001198675) Human Untagged Clone
Tag:	Tag Free
Symbol:	TLCD5
Synonyms:	TMEM136
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC329470 representing NM_001198675. Blue=Insert sequence Red=Cloning site Green=Tag(s)

ATGCTGGCTCATCACACATTGAGTATCTTGGGCATTATCATGGCCCTTGTGCTTGGGGAGTCTGGCACA
GAGGTCAATGCAGTCCTCTTTGGAAGTGAGCTTACCAACCCCTTGCTACAGATGCGCTGGTTTCTCCGG
GAAACAGGGCACTATCACAGTTTCACTGGAGATGTAGTGGACTTCTCTTTGTGGCTCTGTTACAGGA
GTGAGGATTGGTGTGGGAGCTTGCCCTCTTTTCTGTGAAATGGTCTCCCCACGCCTAAGTGTTTGTG
AAGGCTGGGGAGTAGCGATGTATGCTGTGTCTTGGTGTTCATGTTTAGCATCTGGCGCTTTCATGG
AGGAAGAGCATCAAGAAGTACCATGCTTGAGAAGCAGCGGAGTGAGGAACGGCAGCTGAAACACAAC
GGACATCTCAAAATACACTAG

Restriction Sites: SgfI-MluI

Plasmid Map:



[View online »](#)

ACCN:	NM_001198675
Insert Size:	435 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001198675.1</u>
RefSeq Size:	3943 bp
RefSeq ORF:	435 bp
Locus ID:	219902
UniProt ID:	<u>Q6ZRR5</u>
Cytogenetics:	11q23.3
Protein Families:	Transmembrane
MW:	16.5 kDa