

Product datasheet for **SC329479**

YY1AP1 (NM_001198900) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	YY1AP1 (NM_001198900) Human Untagged Clone
Tag:	Tag Free
Symbol:	YY1AP1
Synonyms:	GRNG; HCCA1; HCCA2; YY1AP
Vector:	pCMV6-Entry (PS100001)




[View online »](#)

Fully Sequenced ORF: >SC329479 representing NM_001198900.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGGGATTCTCCAACATGGAAGATGATGGCCCAGAAGAGGAGGAGCGTGTGGCTGAGCCTCAAGCTAAC
TTTAACACCCCTCAAGCTCTACGGTTTGAGGAACACTGGCCAACCTACTAAATGAACAACATCAGATA
GCGAAGGAACATTTTGAACAGCTGAAGATGAAGAAACCTTCAGCCAAACAGCAGAAGGAGGTAGAGAAG
GTTAAACCCAGTGTAAGGAAGTTCATCAGACCCTGATTCTGGACCCAGCACAAGGAAGAGACTCCAG
CAGCAGATGCAGCAGCATGTTTCAGCTCTTGACACAAATCCACCTTCTTGCCACCTGCAACCCCAATCTC
AATCCGGAGGCCAGTAGCACCAGGATATGTCTTAAAGAGCTGGGAACCTTTGCTCAAAGCTCCATCGCC
CTTCACCATCAGTACAACCCCAAGTTTCAGACCCTGTTCCAACCTGTAACCTTGATGGGAGCTATGCAG
CTGATTGAAGACTTCAGCACACATGTCAGCATTGACTGCAGCCCTCATAAACTGTCAAGAAGACTGCC
AATGAATTTCCCTGTTTGCCAAAGCAAGTGGCTTGGATCTGGCCACAAGCAAGGTTTTTCATGTATCCA
GAGTTACTCCAGTGTGTTCCCTGAAGGCAAAGAATCCCAGGATAAGATCCTCTTACCAAGGCTGAG
GACAATTTGTTAGCTTTAGGACTGAAGCATTTTGAAGGGACTGAGTTTCTAACCCCTAATCAGCAAG
TACCTTCTAACCTGCAAGACTGCCCGCCAACTGACAGTGAGAATCAAGAACCTCAACATGAACAGAGCT
CCTGACAACATCATTAAATTTTATAAGAAGACCAAACAGCTGCCAGTCTAGGAAAATGCTGTGAAGAG
ATCCAGCCACATCAGTGGGAAGCCACCTATAGAGAGAGAAGAACACCGGCTCCCATTCTGGTTAAAGGCC
AGTCTGCCATCCATCCAGGAAGAACTGCGGCACATGGCTGATGGTGTAGAGAGGTAGGAAAATGACT
GGAACCACTGAGATCAACTCAGATCAAGGCCTAGAAAAAGACAACCTCAGAGTTGGGGAGTGAAACTCGG
TACCCACTGCTATTGCCTAAGGGTGTAGTCTGAAACTGAAGCCAGTTGCCGACCGTTTTCCCAAGAAG
GCTTGGAGACAGAAGCGTTCATCAGTCTGAAACCCCTCCTTATCCAACCCAGCCCTCTCTCCAGCCC
AGCTTCAACCCCTGGGAAAACACCAGCCCAATCAACTCATTGAGAAGCCCTCCGAGCAAAATGGTGCTC
CGGATTCTCACCAATACAGCCAGCCACTGTTTTACAGACAGTCCAGGTGTCCCTCCACTGGGGGTC
AGTGGAGGTGAGAGTTTTGAGTCTCCTGCAGCACTGCCTGCTATGCCCCCTGAGGCCAGGACAAGCTTC
CCTCTGTCTGAGTCCCAGACTTTGCTCTCTTCTGCCCTGTGCCAAGGTAATGATGCCCTCCCCTGCC
TCTTCCATGTTTCGAAAGCCATATGTGAGACGGAGACCCTCAAAAAGAAGGGGAGCCAGGGCCTTTCGC
TGTATCAAACCTGCCCTGTTATCCACCCTGCATCTGTTATCTTCACTGTTCTGCTACCACTGTGAAG
ATTGTGAGCCTTGGCGGTGGCTGTAACATGATCCAGCCTGTCAATGCGGCTGTGGCCAGAGTCCCAG
ACTATCCCATCGCCACCCTCTGGTTAACCTACTTCTTCCCCTGTCCATTGAACCAGCCCTTGTG
GCCTCCTCTGTCTCACCTTAATTGTTTCTGGCAATTCTGTGAATCTTCTATACCATCCACCCCTGAA
GATAAGGCCACATGAATGTGGACATTGCTTGTGCTGTGGCTGATGGGAAAATGCCTTTCAGGGCCTA
GAACCCAAATTAGAGCCCAGGAACTATCTCTCTCTGCTACTGTTTTCCCAAGTGGAAACATAGC
CCAGGGCCTCCACCAGTCGATAAACAGTGCCAAGAAGGATTGTCAGAGAACAGTGCCTATCGCTGGACC
GTTGTGAAAACAGAGGAGGGAAGGCAAGCTCTGGAGCCGCTCCCTCAGGGCATCCAGGAGTCTCTAAAC
AACTCTTCCCCTGGGGATTTAGAGGAAGTTGTCAAGATGGAACCTGAAGATGCTACAGAGGAAATCAGT
GGATTTCTTGA
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001198900

Insert Size: 2220 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001198900.1](#)

RefSeq Size: 2714 bp

RefSeq ORF: 2220 bp

Locus ID: 55249

UniProt ID: [Q9H869](#)

Cytogenetics: 1q22

MW: 81.7 kDa

Gene Summary:

The encoded gene product presumably interacts with YY1 protein; however, its exact function is not known. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]

Transcript Variant: This variant (7) has an alternate 5' sequence, resulting in a downstream AUG start codon, as compared to variant 10. The resulting isoform (1) has a shorter N-terminus, as compared to isoform 6. Variants 1, 6 and 7 encode the same isoform 1.