

Product datasheet for **SC329486**

AMMECR1L (NM_001199140) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	AMMECR1L (NM_001199140) Human Untagged Clone
Tag:	Tag Free
Symbol:	AMMECR1L
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC329486 representing NM_001199140. Blue=Insert sequence Red=Cloning site Green=Tag(s)

ATGGGAAAAGACGTTGTGTTCCCTCCACTCGAGCCCAAGTTGGCAGCAGGCTGTTGTGGGGTCAAGAAG
 CCCAAATTATCTGGAAGTGAACGCACAGTCACGGGAATCAGTCCACAAGTGTCCCGGCTCTAGTTCA
 GGACCTCTTCAAACACCAGCATGTGGACAGCAGCAGTGGACGGGAGAATGTGTCAGACTTAACCTG
 GGACCTGGAACTCTCCCATCACACGAATGAATCCCGCATCGGGAGCGCTGAGCCCTCTTCCCGGCCT
 AATGGAAGTGGCAACACCACTAAGAATCTGGTGGTGACTGCAGAGATGTGCTGCTACTGCTTCGACGTA
 CTCTACTGTCACCTCTATGGCTTCCCACAGCCACGACTTCCTAGATTACCAATGACCCCTATCCGCTC
 TTTGTGACGTGGAAGACAGGGCGGGACAAGCGGCTTCGTGGCTGCATTGGGACCTTCTCAGCCATGAAT
 CTTCAATCAGGACTCAGGGAATACACGTTAACCAGTGCCTTAAGGACAGCCGATTTCCCCCCTGACC
 CGAGAGGAGCTGCCTAACTTTTCTGCTCTGTCTCCCTCCTTACTAATTTGAGGATGCCAGTGATTAC
 CTGGACTGGGAGGTAGGGGTCCATGGGATTGGAATTGAATTCATTAATGAAAAGGTGTCAAACGCACA
 GCCACATATTTACCTGAGGTTGCTAAGGAACAAGACTGGGATCAGATCCAGACAATAGACTCCTTGCTC
 AGGAAAGGTGGCTTTAAGCTCCAATTACCAGTGAATTCAGAAAAACGATCAAACCTACCAGGTACCGA
 AGTGAGAAGGTGACAATCAGTTACGCAGAGTATATTGCTTCCCGACAGCACTGTTCCAGAACGGCACT
 CTTCAATGCCCCGCCCTCTACAATCATTACTCC

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001199140
Insert Size:	933 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001199140.1](#)

RefSeq Size: 4726 bp

RefSeq ORF: 933 bp

Locus ID: 83607

UniProt ID: [Q6DCA0](#)

Cytogenetics: 2q14.3

MW: 34.5 kDa